

# Package ‘sesame’

October 17, 2020

**Type** Package

**Title** Tools For Analyzing Illumina Infinium DNA Methylation Arrays

**Description** Tools For analyzing Illumina Infinium DNA methylation arrays.

**Version** 1.6.0

**Depends** R (>= 3.6), sesameData, methods

**License** MIT + file LICENSE

**RoxygenNote** 7.0.2

**Imports** BiocParallel, grDevices, utils, illuminaio, MASS,  
GenomicRanges, IRanges, grid, preprocessCore, stats, S4Vectors,  
randomForest, wheatmap, ggplot2, parallel, matrixStats,  
DNAcopy, HDF5Array, SummarizedExperiment

**Suggests** scales, knitr, rmarkdown, testthat, minfi,  
FlowSorted.CordBloodNorway.450k, FlowSorted.Blood.450k, dplyr,  
tidyverse, BiocStyle

**Encoding** UTF-8

**VignetteBuilder** knitr

**URL** <https://github.com/zwdzwd/sesame>

**BugReports** <https://github.com/zwdzwd/sesame/issues>

**biocViews** DNAMethylation, MethylationArray, Preprocessing,  
QualityControl

**Collate** 'GEO.R' 'QC.R' 'SigSetList.R' 'sesame.R' 'SigSetMethods.R'  
'age.R' 'background\_correction.R' 'cell\_composition.R' 'cnv.R'  
'controls.R' 'detection.R' 'differential\_methylation.R'  
'dye\_bias.R' 'fileSet.R' 'open.R' 'sesamize.R'  
'simulate\_data.R' 'snp.R' 'utils.R' 'visualize.R'

**git\_url** <https://git.bioconductor.org/packages/sesame>

**git\_branch** RELEASE\_3\_11

**git\_last\_commit** 92c330a

**git\_last\_commit\_date** 2020-04-27

**Date/Publication** 2020-10-16

**Author** Wanding Zhou [aut, cre],  
Hui Shen [aut],  
Timothy Triche [ctb],  
Bret Barnes [ctb]

**Maintainer** Wanding Zhou <[zhouwanding@gmail.com](mailto:zhouwanding@gmail.com)>

## R topics documented:

sesame-package . . . . .	4
as.data.frame.sesameQC . . . . .	4
BetaValueToMValue . . . . .	5
binSignals . . . . .	6
bisConversionControl . . . . .	6
buildControlMatrix450k . . . . .	7
chipAddressToSignal . . . . .	7
cnSegmentation . . . . .	8
ctl . . . . .	9
ctl<- . . . . .	9
detectionPfixedNorm . . . . .	10
detectionPnegEcdf . . . . .	11
detectionPnegNorm . . . . .	11
detectionPnegNormGS . . . . .	12
detectionPnegNormTotal . . . . .	13
detectionPoobEcdf . . . . .	13
detectionZero . . . . .	14
diffRefSet . . . . .	15
DML . . . . .	15
DMR . . . . .	16
dyeBiasCorr . . . . .	17
dyeBiasCorrMostBalanced . . . . .	18
dyeBiasCorrTypeINorm . . . . .	18
estimateCellComposition . . . . .	19
estimateLeukocyte . . . . .	19
getAFTypeIbySumAlleles . . . . .	20
getBetas . . . . .	21
getBinCoordinates . . . . .	22
getNormCtls . . . . .	22
getProbesByGene . . . . .	23
getProbesByRegion . . . . .	23
getProbesByTSS . . . . .	24
getRefSet . . . . .	25
getSegment . . . . .	26
getSexInfo . . . . .	26
IG . . . . .	27
IG<- . . . . .	27
II . . . . .	28
II<- . . . . .	28
inferEthnicity . . . . .	29
inferSex . . . . .	30
inferSexKaryotypes . . . . .	30
inferTypeIChannel . . . . .	31
initFileSet . . . . .	32
IR . . . . .	32
IR<- . . . . .	33
makeExampleSeSAMeDataSet . . . . .	34
makeExampleTinyEPICDataSet . . . . .	34
mapFileSet . . . . .	35
meanIntensity . . . . .	35

MValueToBetaValue	36
noob	36
noobsb	37
oobG	38
oobG<-	38
oobR	39
oobR<-	39
openSesame	40
openSesameToFile	41
parseGEOSignalABFile	42
predictAgeHorvath353	42
predictAgePheno	43
predictAgeSkinBlood	43
print.fileSet	44
print.sesameQC	45
probeNames	45
pval	46
pval<-	46
readFileSet	47
readIDATpair	48
reopenSesame	48
RGChannelSetToSigSets	49
searchIDATprefixes	49
segmentBins	50
sesameQC	51
sesamize	51
show,SigSet-method	52
SigSet-class	52
SigSetList	53
SigSetList-class	54
SigSetList-methods	54
SigSetListFromIDATs	55
SigSetListFromPath	55
SigSetsToRGChannelSet	56
SigSetToRatioSet	56
sliceFileSet	57
SNPcheck	58
subsetSignal	58
topLoci	59
topSegments	59
totalIntensities	60
totalIntensityZscore	60
twoCompsEst2	61
visualizeGene	61
visualizeProbes	62
visualizeRegion	63
visualizeSegments	64

**sesame-package**      *Analyze DNA methylation data*

## Description

SEnsible and step-wise analysis of DNA methylation data

## Details

This package complements array functionalities that allow processing >10,000 samples in parallel on clusters.

## Author(s)

Wanding Zhou <[Wanding.Zhou@vai.org](mailto:Wanding.Zhou@vai.org)>, Hui Shen <[Hui.Shen@vai.org](mailto:Hui.Shen@vai.org)> Timothy J Triche Jr <[Tim.Triche@vai.org](mailto:Tim.Triche@vai.org)>

## See Also

Useful links:

- <https://github.com/zwdzwd/sesame>
- Report bugs at <https://github.com/zwdzwd/sesame/issues>

## Examples

```
sset <- readIDATpair(sub('_Grn.idat',' ',system.file('extdata','4207113116_A_Grn.idat',package='sesameData')))

## The OpenSesame pipeline
betas <- openSesame(sset)
```

**as.data.frame.sesameQC**  
*Coerce a sesameQC into a dataframe*

## Description

Coerce a sesameQC into a dataframe

## Usage

```
## S3 method for class 'sesameQC'
as.data.frame(x, row.names = NULL, optional = FALSE, ...)
```

**Arguments**

x	a sesameQC object
row.names	see as.data.frame
optional	see as.data.frame
...	see as.data.frame

**Value**

a data.frame

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
qc <- sesameQC(sset)
df <- as.data.frame(qc)
```

**BetaValueToMValue**

*Convert beta-value to M-value*

**Description**

Logit transform a beta value vector to M-value vector.

**Usage**

```
BetaValueToMValue(b)
```

**Arguments**

b	vector of beta values
---	-----------------------

**Details**

Convert beta-value to M-value (aka logit transform)

**Value**

a vector of M values

**Examples**

```
BetaValueToMValue(c(0.1, 0.5, 0.9))
```

<code>binSignals</code>	<i>Bin signals from probe signals</i>
-------------------------	---------------------------------------

### Description

require GenomicRanges

### Usage

```
binSignals(probe.signals, bin.coords, probe.coords)
```

### Arguments

<code>probe.signals</code>	probe signals
<code>bin.coords</code>	bin coordinates
<code>probe.coords</code>	probe coordinates

### Value

bin signals

<code>bisConversionControl</code>	<i>Compute internal bisulfite conversion control</i>
-----------------------------------	--

### Description

Compute GCT score for internal bisulfite conversion control. The function takes a SigSet as input. The higher the GCT score, the more likely the incomplete conversion. The lower the GCT score, the more likely over-conversion.

### Usage

```
bisConversionControl(sset, use.median = FALSE)
```

### Arguments

<code>sset</code>	signal set
<code>use.median</code>	use median to compute GCT instead of mean

### Value

GCT score (the higher, the more incomplete conversion)

### Examples

```
sset <- makeExampleSeSAMeDataSet('HM450')
bisConversionControl(sset)
```

---

**buildControlMatrix450k**

*Build control summary matrix*

---

### Description

The function takes a `SigSet` as input and outputs the control matrix summary vector. This vector summarizes one single QC metric for the array control. This includes bisulfite control, stain signal extension efficiency and more.

### Usage

```
buildControlMatrix450k(sset)
```

### Arguments

sset	an object of class <code>SigSet</code>
------	--

### Value

a vector with control summaries

### Examples

```
sset <- makeExampleSeSAMeDataSet()
control.summary <- buildControlMatrix450k(sset)
```

---

**chipAddressToSignal**     *Lookup address in one sample*

---

### Description

Lookup address and transform address to probe

### Usage

```
chipAddressToSignal(dm, manifest, controls = NULL, readNBeads = FALSE)
```

### Arguments

dm	data frame in chip address, 2 columns: cy3/Grn and cy5/Red
manifest	a data frame with columns Probe_ID, M, U and col
controls	a data frame with columns Address and Name. This is optional but might be necessary for some preprocessing methods that depends on these control probes. This is left for backward compatibility. Updated version should have controls consolidated into manifest.
readNBeads	whether to read bead signal

## Details

Translate data in chip address to probe address. Type I probes can be separated into Red and Grn channels. The methylated allele and unmethylated allele are at different addresses. For type II probes methylation allele and unmethylated allele are at the same address. Grn channel is for methylated allele and Red channel is for unmethylated allele. The out-of-band signals are type I probes measured using the other channel.

## Value

a SigSet, indexed by probe ID address

cnSegmentation	<i>Perform copy number segmentation</i>
----------------	---

## Description

Perform copy number segmentation using the signals in the signal set. The function takes a SigSet for the target sample and a set of normal SigSet for the normal samples. An optional arguments specifies the version of genome build that the inference will operate on. The function outputs an object of class CNSegment with signals for the segments ( seg.signals), the bin coordinates ( bin.coords) and bin signals (bin.signals).

## Usage

```
cnSegmentation(sset, ssets.normal, refversion = c("hg19", "hg38"))
```

## Arguments

sset	SigSet
ssets.normal	SigSet for normalization
refversion	hg19 or hg38

## Value

an object of CNSegment

## Examples

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
ssets.normal <- sesameDataGet('EPIC.5.normal')
seg <- cnSegmentation(sset, ssets.normal)
```

---

ctl	<i>ctl getter generic</i>
-----	---------------------------

---

**Description**

ctl getter generic  
Get ctl slot of SigSet class

**Usage**

```
ctl(x)  
  
## S4 method for signature 'SigSet'  
ctl(x)
```

**Arguments**

x object of SigSet

**Value**

The ctl slot of SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset  
head(ctl(sset))
```

---

ctl<-	<i>ctl replacement generic</i>
-------	--------------------------------

---

**Description**

ctl replacement generic  
Replace ctl slot of SigSet class

**Usage**

```
ctl(x) <- value  
  
## S4 replacement method for signature 'SigSet'  
ctl(x) <- value
```

**Arguments**

x object of SigSet  
value new value

**Value**

a new SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
df <- ctl(sset)
df[1,1] <- 10
ctl(sset) <- df
```

**detectionPfixedNorm**    *Detection P-value based on normal fitting with gived parameters*

**Description**

The function takes a SigSet as input, computes detection p-value using negative control probes parametrized in a normal distribution and returns a new SigSet with an updated pval slot.

**Usage**

```
detectionPfixedNorm(
  sset,
  muG = 500,
  sdG = 200,
  muR = 500,
  sdR = 200,
  force = FALSE
)
```

**Arguments**

sset	a SigSet
muG	mean of background in Grn channel
sdG	SD of background in Grn channel
muR	mean of background in Red channel
sdR	SD of background in Red channel
force	force rerun even if result already exists

**Details**

Background of Grn and Red are estimated separately from a fixed normal distribution. p-value is taken from the minimum of the p-value of the two alleles (color depends on probe design).

**Value**

detection p-value

**Examples**

```
sset <- makeExampleSeSAMeDataSet()
sset <- detectionPfixedNorm(sset)
```

detectionPnegEcdf

*Detection P-value based on ECDF of negative control***Description**

The function takes a SigSet as input, computes detection p-value using negative control probes' empirical distribution and returns a new SigSet with an updated pval slot.

**Usage**

```
detectionPnegEcdf(sset, force = FALSE)
```

**Arguments**

sset	a SigSet
force	force rerun even if result already exists

**Value**

detection p-value

**Examples**

```
sset <- makeExampleSeSAMeDataSet()
sset <- detectionPnegEcdf(sset)
```

detectionPnegNorm

*Detection P-value based on normal fitting the negative controls***Description**

The function takes a SigSet as input, computes detection p-value using negative control probes parametrized in a normal distribution and returns a new SigSet with an updated pval slot.

**Usage**

```
detectionPnegNorm(sset, force = FALSE)
```

**Arguments**

sset	a SigSet
force	force rerun even if result already exists

**Details**

Background of Grn and Red are estimated separately from negative control probes-parameterized normal distribution. p-value is taken from the minimum of the p-value of the two alleles (color depends on probe design).

**Value**

detection p-value

**Examples**

```
sset <- makeExampleSeSAMeDataSet()
sset <- detectionPnegNorm(sset)
```

**detectionPnegNormGS**      *Detection P-value emulating Genome Studio*

**Description**

The function takes a `SigSet` as input, computes detection p-value using negative control probes parametrized in a normal distribution a la Genome Studio and returns a new `SigSet` with an updated `pval` slot.

**Usage**

```
detectionPnegNormGS(sset, force = FALSE)
```

**Arguments**

<code>sset</code>	a <code>SigSet</code>
<code>force</code>	force rerun even if result already exists

**Details**

P-value is calculated using negative control probes as the estimate of background where Grn channel and Red channel are merged. But when estimating p-value the Red and Grn are summed (non-ideal).

**Value**

detection p-value

**Examples**

```
sset <- makeExampleSeSAMeDataSet()
sset <- detectionPnegNormGS(sset)
```

---

**detectionPnegNormTotal**

*Detection P-value based on normal fitting the negative controls, channels are first summed*

---

### Description

The function takes a SigSet as input, computes detection p-value using negative control probes parametrized in a normal distribution with the two channels summed first and returns a new SigSet with an updated pval slot. The SD is summed to emulate the SD of the summed signal (not the most accurate treatment).

### Usage

```
detectionPnegNormTotal(sset, force = FALSE)
```

### Arguments

sset	a SigSet
force	force rerun even if result already exists

### Value

detection p-value

### Examples

```
sset <- makeExampleSeSAMeDataSet()
sset <- detectionPnegNormTotal(sset)
```

---

**detectionPoobEcdf**

*Detection P-value based on ECDF of out-of-band signal*

---

### Description

aka pOOBAH (p-vals by Out-Of-Band Array Hybridization)

### Usage

```
detectionPoobEcdf(sset, force = FALSE)

pOOBAH(sset, force = FALSE)
```

### Arguments

sset	a SigSet
force	force rerun even if result already exists

## Details

The function takes a SigSet as input, computes detection p-value using out-of-band probes empirical distribution and returns a new SigSet with an updated pval slot.

## Value

detection p-value

## Examples

```
sset <- makeExampleSeSAMeDataSet()
sset <- detectionPoobEcdf(sset)

sset <- makeExampleSeSAMeDataSet()
sset <- p00BAH(sset)
```

**detectionZero**

*Detection P-value set to all zero*

## Description

Detection P-value set to all zero

## Usage

```
detectionZero(sset, force = FALSE)
```

## Arguments

sset	a SigSet
force	force rerun even if result already exists

## Value

detection p-value set to all zero

## Examples

```
sset <- makeExampleSeSAMeDataSet()
sset <- detectionZero(sset)
```

---

diffRefSet	<i>Restrict refset to differentially methylated probes use with care, might introduce bias</i>
------------	--

---

**Description**

The function takes a matrix with probes on the rows and cell types on the columns and output a subset matrix and only probes that show discordant methylation levels among the cell types.

**Usage**

```
diffRefSet(g)
```

**Arguments**

g	a matrix with probes on the rows and cell types on the columns
---	--

**Value**

g a matrix with a subset of input probes (rows)

**Examples**

```
g <- diffRefSet(getRefSet(platform='HM450'))
```

---

DML	<i>Test differential methylation on each locus</i>
-----	--

---

**Description**

The function takes a beta value matrix with probes on the rows and samples on the columns. It also takes a sample information data frame (sample.data) and formula for testing. The function outputs a list of coefficient tables for each factor tested.

**Usage**

```
DML(
  betas,
  sample.data,
  formula,
  se.1b = 0.06,
  balanced = FALSE,
  cf.test = NULL
)
```

**Arguments**

<code>betas</code>	beta values
<code>sample.data</code>	data frame for sample information, column names are predictor variables (e.g., sex, age, treatment, tumor/normal etc) and are referenced in formula. Rows are samples.
<code>formula</code>	formula
<code>se.lb</code>	lower bound to standard error of slope, lower this to get more difference of small effect size.
<code>balanced</code>	whether design is balanced or not. default to FALSE, when unbalanced will use Welch's method to estimate standard error. balance=TRUE is faster.
<code>cf.test</code>	factors to test (default to all factors in formula except intercept). Use "all" for all factors.

**Value**

`cf` - a list of coefficient tables for each factor

**Examples**

```
data <- sesameDataGet('HM450.76.TCGA.matched')
cf <- DML(data$betas, data$sampleInfo, ~type)
```

DMR

*Find Differentially Methylated Region (DMR)***Description**

This subroutine uses Euclidean distance to group CpGs and then combine p-values for each segment. The function performs DML test first if `cf` is NULL. It groups the probe testing results into differential methylated regions in a coefficient table with additional columns designating the segment ID and statistical significance (P-value) testing the segment.

**Usage**

```
DMR(
  betas,
  sample.data = NULL,
  formula = NULL,
  cf = NULL,
  dist.cutoff = NULL,
  seg.per.locus = 0.5,
  platform = c("EPIC", "HM450"),
  refversion = c("hg38", "hg19"),
  ...
)
```

**Arguments**

betas	beta values for distance calculation
sample.data	data frame for sample information, column names are predictor variables (e.g., sex, age, treatment, tumor/normal etc) and are referenced in formula. Rows are samples.
formula	formula
cf	coefficient table from diffMeth, when NULL will be computed from beta. If cf is given, sample.data and formula are ignored.
dist.cutoff	distance cutoff (default to use dist.cutoff.quantile)
seg.per.locus	number of segments per locus higher value leads to more segments
platform	EPIC or HM450
refversion	hg38 or hg19
...	additional parameters to DML

**Value**

coefficient table with segment ID and segment P-value

**Examples**

```
data <- sesameDataGet('HM450.76.TCGA.matched')
cf <- DMR(data$betas, data$sampleInfo, ~type)
```

dyeBiasCorr

*Correct dye bias in by linear scaling.*

**Description**

The function takes a SigSet as input and scale both the Grn and Red signal to a reference (ref) level. If the reference level is not given, it is set to the mean intensity of all the in-band signals. The function returns a SigSet with dye bias corrected.

**Usage**

```
dyeBiasCorr(sset, ref = NULL)
```

**Arguments**

sset	a SigSet
ref	reference signal level

**Value**

a normalized SigSet

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
sset.db <- dyeBiasCorr(sset)
```

**dyeBiasCorrMostBalanced***Correct dye bias using most balanced sample as the reference***Description**

The function chose the reference signal level from a list of SigSet. The chosen sample has the smallest difference in Grn and Red signal intensity as measured using the normalization control probes. In practice, it doesn't matter which sample is chosen as long as the reference level does not deviate much. The function returns a list of SigSets with dye bias corrected.

**Usage**

```
dyeBiasCorrMostBalanced(ssets)
```

**Arguments**

ssets	a list of normalized SigSets
-------	------------------------------

**Value**

a list of normalized SigSets
------------------------------

**Examples**

```
ssets <- sesameDataGet('HM450.10.TCGA.BLCA.normal')
ssets.db <- dyeBiasCorrMostBalanced(ssets)
```

**dyeBiasCorrTypeINorm** *Dye bias correction by matching green and red to mid point***Description**

This function compares the Type-I Red probes and Type-I Grn probes and generates and mapping to correct signal of the two channels to the middle. The function takes one single SigSet and returns a SigSet with dye bias corrected.

**Usage**

```
dyeBiasCorrTypeINorm(sset)
```

**Arguments**

sset	a SigSet
------	----------

**Value**

a SigSet after dye bias correction.
-------------------------------------

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
sset.db <- dyeBiasCorrTypeINorm(sset)
```

---

estimateCellComposition*Estimate cell composition using reference*

---

**Description**

This is a reference-based cell composition estimation. The function takes a reference methylation status matrix (rows for probes and columns for cell types, can be obtained by getRefSet function) and a query beta value measurement. The length of the target beta values should be the same as the number of rows of the reference matrix. The method assumes one unknown component. It outputs a list containing the estimated cell fraction, the error of optimization and methylation status of the unknown component.

**Usage**

```
estimateCellComposition(g, q, refine = TRUE, dichotomize = FALSE, ...)
```

**Arguments**

g	reference methylation
q	target measurement: length(q) == nrow(g)
refine	to refine estimate, takes longer
dichotomize	to dichotomize query beta value before estimate, this relieves unclean background subtraction
...	extra parameters for optimization, this includes temp - annealing temperature (0.5) maxIter - maximum iteration to stop after converge (1000) delta - delta score to reset counter (0.0001) verbose - output debug info (FALSE)

**Value**

a list of fraction, min error and unknown component methylation state

---

estimateLeukocyte*Estimate leukocyte fraction using a two-component model*

---

**Description**

The method assumes only two components in the mixture: the leukocyte component and the target tissue component. The function takes the beta values matrix of the target tissue and the beta value matrix of the leukocyte. Both matrices have probes on the row and samples on the column. Row names should have probe IDs from the platform. The function outputs a single numeric describing the fraction of leukocyte.

**Usage**

```
estimateLeukocyte(
  betas.tissue,
  betas.leuko = NULL,
  betas.tumor = NULL,
  platform = c("EPIC", "HM450", "HM27")
)
```

**Arguments**

<code>betas.tissue</code>	tissue beta value matrix (#probes X #samples)
<code>betas.leuko</code>	leukocyte beta value matrix, if missing, use the SeSAMe default by infinium platform
<code>betas.tumor</code>	optional, tumor beta value matrix
<code>platform</code>	"HM450", "HM27" or "EPIC"

**Value**

leukocyte estimate, a numeric vector

**Examples**

```
betas.tissue <- sesameDataGet('HM450.1.TCGA.PAAD')$betas
estimateLeukocyte(betas.tissue)
```

**getAFTypeIbySumAlleles**

*Get allele frequency treating type I by summing alleles*

**Description**

Takes a `SigSet` as input and returns a numeric vector containing extra allele frequencies based on Color-Channel-Switching (CCS) probes. If no CCS probes exist in the `SigSet`, then an `numeric(0)` is returned.

**Usage**

```
getAFTypeIbySumAlleles(sset, known.ccs.only = TRUE)
```

**Arguments**

<code>sset</code>	<code>SigSet</code>
<code>known.ccs.only</code>	consider only known CCS probes

**Value**

beta values

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
betas <- getAFTypeIbySumAlleles(sset)
```

getBetas

*Get beta Values***Description**

sum.typeI is used for rescuing beta values on Color-Channel-Switching CCS probes. The function takes a SigSet and returns beta value except that Type-I in-band signal and out-of-band signal are combined. This prevents color-channel switching due to SNPs.

**Usage**

```
getBetas(
  sset,
  quality.mask = TRUE,
  nondetection.mask = TRUE,
  correct.switch = TRUE,
  mask.use.tcga = FALSE,
  pval.threshold = 0.05,
  pval.method = NULL,
  sum.TypeI = FALSE
)
```

**Arguments**

sset	SigSet
quality.mask	whether to mask low quality probes
nondetection.mask	whether to mask nondetection
correct.switch	whether to correct switch
mask.use.tcga	whether to use TCGA masking, only applies to HM450
pval.threshold	p-value threshold for nondetection mask
pval.method	method for detection threshold, like pOOBAH, PnegEcdf
sum.TypeI	whether to sum type I channels

**Value**

a numeric vector, beta values

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
betas <- getBetas(sset)
```

`getBinCoordinates`      *Get bin coordinates*

### Description

requires GenomicRanges, IRanges

### Usage

```
getBinCoordinates(seqInfo, gapInfo, probe.coords)
```

### Arguments

<code>seqInfo</code>	chromosome information object
<code>gapInfo</code>	chromosome gap information
<code>probe.coords</code>	probe coordinates

### Value

`bin.coords`

`getNormCtl`      *get normalization control signal*

### Description

get normalization control signal from SigSet. The function optionally takes mean for each channel.

### Usage

```
getNormCtl(sset, average = FALSE)
```

### Arguments

<code>sset</code>	a SigSet
<code>average</code>	whether to average

### Value

a data frame of normalization control signals

### Examples

```
sset <- readIDATpair(file.path(system.file(
  'extdata', '', package='sesameData'), '4207113116_B'))

df.ctl <- getNormCtl(sset)
```

---

getProbesByGene      *Get Probes by Gene*

---

## Description

Get probes mapped to a gene. All transcripts for the gene are considered. The function takes a gene name as appears in UCSC RefGene database. The platform and reference genome build can be changed with ‘platform’ and ‘refversion’ options. The function returns a vector of probes that falls into the given gene.

## Usage

```
getProbesByGene(  
  geneName,  
  platform = c("EPIC", "HM450"),  
  upstream = 0,  
  dwstream = 0,  
  refversion = c("hg38", "hg19"))  
)
```

## Arguments

geneName	gene name
platform	EPIC or HM450
upstream	number of bases to expand upstream of target gene
dwstream	number of bases to expand downstream of target gene
refversion	hg38 or hg19

## Value

probes that fall into the given gene

## Examples

```
probes <- getProbesByGene('CDKN2A', upstream=500, dwstream=500)
```

---

getProbesByRegion      *Get probes by genomic region*

---

## Description

The function takes a genomic coordinate and output the a vector of probes on the specified platform that falls in the given genomic region.

**Usage**

```
getProbesByRegion(
  chrm,
  beg = 1,
  end = -1,
  platform = c("EPIC", "HM450"),
  refversion = c("hg38", "hg19")
)
```

**Arguments**

chrm	chromosome
beg	begin, 1 if omitted
end	end, chromosome end if omitted
platform	EPIC or HM450
refversion	hg38 or hg19

**Value**

probes that fall into the given region

**Examples**

```
getProbesByRegion('chr5', 135413937, 135419936,
  refversion = 'hg19', platform = 'HM450')
```

---

getProbesByTSS

*Get Probes by Gene Transcription Start Site (TSS)*

---

**Description**

Get probes mapped to a TSS. All transcripts for the gene are considered. The function takes a gene name as appears in UCSC RefGene database. The platform and reference genome build can be changed with ‘platform’ and ‘refversion’ options. The function returns a vector of probes that falls into the TSS region of the gene.

**Usage**

```
getProbesByTSS(
  geneName,
  upstream = 1500,
  dwstream = 1500,
  platform = c("EPIC", "HM450"),
  refversion = c("hg38", "hg19")
)
```

**Arguments**

geneName	gene name
upstream	the number of base pairs to expand upstream the TSS
dwstream	the number of base pairs to expand dwstream the TSS
platform	EPIC or HM450
refversion	hg38 or hg19

**Value**

probes that fall into the given gene

**Examples**

```
probes <- getProbesByTSS('CDKN2A')
```

---

getRefSet

*Retrieve reference set*

---

**Description**

The function retrieves the curated reference DNA methylation status for a set of cell type names under the Infinium platform. Supported cell types include "CD4T", "CD19B", "CD56NK", "CD14Monocytes", "granulocytes", "scFat", "skin" etc. See package sesameData for more details. The function output a matrix with probes on the rows and specified cell types on the columns. 0 suggests unmethylation and 1 suggests methylation. Intermediate methylation and nonclusive calls are left with NA.

**Usage**

```
getRefSet(cells = NULL, platform = c("EPIC", "HM450"))
```

**Arguments**

cells	reference cell types
platform	EPIC or HM450

**Value**

g, a 0/1 matrix with probes on the rows and specified cell types on the columns.

**Examples**

```
betas <- getRefSet('CD4T', platform='HM450')
```

<code>getSegment</code>	<i>Select segment from coefficient table</i>
-------------------------	--

### Description

This function takes a coefficient table and returns a subset of the table targeting only the specified segment using segment ID.

### Usage

```
getSegment(cf1, seg.id)
```

### Arguments

<code>cf1</code>	coefficient table of one factor from DMR
<code>seg.id</code>	segment ID

### Value

coefficient table from given segment

### Examples

```
data <- sesameDataGet('HM450_76.TCGA.matched')
cf <- DMR(data$betas, data$sampleInfo, ~type)
getSegment(cf[[1]], cf[[1]][['Seg.ID']][1])
```

<code>getSexInfo</code>	<i>Get sex-related information</i>
-------------------------	------------------------------------

### Description

The function takes a `SigSet` and returns a vector of three numerics: the median intensity of chrY probes; the median intensity of chrX probes; and fraction of intermediate chrX probes. chrX and chrY probes excludes pseudo-autosomal probes.

### Usage

```
getSexInfo(sset)
```

### Arguments

<code>sset</code>	a <code>SigSet</code>
-------------------	-----------------------

### Value

medianY and medianX, fraction of XCI, methylated and unmethylated X probes, median intensities of auto-chromosomes.

### Examples

```
sset <- makeExampleSeSAMeDataSet()
getSexInfo(sset)
```

---

IG

*IG getter generic*

---

### Description

IG getter generic

Get IG slot of SigSet class

### Usage

IG(x)

```
## S4 method for signature 'SigSet'  
IG(x)
```

### Arguments

x object of SigSet

### Value

The IG slot of SigSet

### Examples

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset  
head(IG(sset))
```

---

IG<-

*IG replacement generic*

---

### Description

IG replacement generic

Replace IG slot of SigSet class

### Usage

IG(x) <- value

```
## S4 replacement method for signature 'SigSet'  
IG(x) <- value
```

### Arguments

x object of SigSet  
value new value

**Value**

a new SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
df <- IG(sset)
df[1,1] <- 10
IG(sset) <- df
```

II

*II getter generic***Description**

II getter generic

Get II slot of SigSet class

**Usage**

```
II(x)

## S4 method for signature 'SigSet'
II(x)
```

**Arguments**

x object of SigSet

**Value**

The II slot of SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
head(II(sset))
```

II&lt;-

*II replacement generic***Description**

II replacement generic

Replace II slot of SigSet class

**Usage**

```
II(x) <- value  
  
## S4 replacement method for signature 'SigSet'  
II(x) <- value
```

**Arguments**

x	object of SigSet
value	new value

**Value**

a new SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset  
df <- II(sset)  
df[1,1] <- 10  
II(sset) <- df
```

---

inferEthnicity	<i>Infer Ethnicity</i>
----------------	------------------------

---

**Description**

This function uses both the built-in rsprobes as well as the type I Color-Channel-Switching probes to infer ethnicity.

**Usage**

```
inferEthnicity(sset)
```

**Arguments**

sset	a SigSet
------	----------

**Details**

sset better be background subtracted and dyebias corrected for best accuracy

**Value**

string of ethnicity

**Examples**

```
sset <- makeExampleSeSAMeDataSet("HM450")  
inferEthnicity(sset)
```

**inferSex***Infer Sex***Description**

Infer Sex

**Usage**`inferSex(sset)`**Arguments**

sset            a SigSet

**Value**

'F' or 'M' We established our sex calling based on the median intensity of chromosome X, Y and the fraction of intermediately methylated probes among the identified X-linked probes. This is similar to the approach by Minfi (Aryee et al., 2014) but also different in that we used the fraction of intermediate beta value rather than median intensity for all chromosome X probes. Instead of using all probes from the sex chromosomes, we used our curated set of Y chromosome probes and X-linked probes which exclude potential cross-hybridization and pseudo-autosomal effect.

XXY male (Klinefelter's), 45,X female (Turner's) can confuse the model sometimes. Our function works on a single sample.

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
inferSex(sset)
```

**inferSexKaryotypes***Infer Sex Karyotype***Description**

The function takes a SigSet and infers the sex chromosome Karyotype and presence/absence of X-chromosome inactivation (XCI). chrX, chrY and XCI are inferred relatively independently. This function gives a more detailed look of potential sex chromosome aberrations.

**Usage**`inferSexKaryotypes(sset)`**Arguments**

sset            a SigSet

**Value**

Karyotype string, with XCI

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset  
inferSexKaryotypes(sset)
```

---

**inferTypeIChannel**

*Infer and reset color channel for Type-I probes instead of using what is specified in manifest*

---

**Description**

Infer and reset color channel for Type-I probes instead of using what is specified in manifest

**Usage**

```
inferTypeIChannel(  
  sset,  
  switch_failed = FALSE,  
  verbose = FALSE,  
  summary = FALSE  
)
```

**Arguments**

sset	a SigSet
switch_failed	whether to switch failed probes (default to FALSE)
verbose	whether to print correction summary
summary	return summarized numbers only.

**Value**

a SigSet, or numerics if summary == TRUE

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset  
inferTypeIChannel(sset)
```

<code>initFileSet</code>	<i>initialize a fileSet class by allocating appropriate storage</i>
--------------------------	---

### Description

initialize a fileSet class by allocating appropriate storage

### Usage

```
initFileSet(map_path, platform, samples, probes = NULL, inc = 4)
```

### Arguments

<code>map_path</code>	path of file to map
<code>platform</code>	EPIC, HM450 or HM27, consistent with sset@platform
<code>samples</code>	sample names
<code>probes</code>	probe names
<code>inc</code>	bytes per unit data storage

### Value

a sesame::fileSet object

### Examples

```
fset <- initFileSet('mybetas2', 'HM27', c('s1', 's2'))
```

<code>IR</code>	<i>IR getter generic</i>
-----------------	--------------------------

### Description

IR getter generic

Get IR slot of SigSet class

### Usage

```
IR(x)
```

```
## S4 method for signature 'SigSet'
IR(x)
```

### Arguments

<code>x</code>	object of SigSet
----------------	------------------

**Value**

The IR slot of SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
head(IR(sset))
```

---

IR<-

---

*IR replacement generic*

---

**Description**

IR replacement generic

Replace IR slot of SigSet class

**Usage**

```
IR(x) <- value

## S4 replacement method for signature 'SigSet'
IR(x) <- value
```

**Arguments**

x	object of SigSet
value	new value

**Value**

a new SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
df <- IR(sset)
df[1,1] <- 10
IR(sset) <- df
```

---

**makeExampleSeSAMeDataSet**

*Make a simulated SeSAMe data set*

---

## Description

Constructs a simulated SigSet dataset. For the given platform, randomly simulate methylated and unmethylated allele signals. In-band signals were simulated using a  $N(4000, 200)$  normal distribution. Out-of-band signals were simulated using a  $N(400, 200)$  normal distribution. Control signals were simulated using a  $N(400, 300)$  normal distribution.

## Usage

```
makeExampleSeSAMeDataSet(platform = c("HM450", "EPIC", "HM27"))
```

## Arguments

platform	optional, HM450, EPIC or HM27
----------	-------------------------------

## Value

Object of class `SigSet`

## Examples

```
sset <- makeExampleSeSAMeDataSet()
```

---

**makeExampleTinyEPICDataSet**

*Make a tiny toy simulated EPIC data set*

---

## Description

Construct a tiny EPIC SigSet of only 6 probes. In-band signals were simulated using a  $N(4000, 200)$  normal distribution. Out-of-band signals were simulated using a  $N(400, 200)$  normal distribution. Control signals were simulated using a  $N(400, 300)$  normal distribution.

## Usage

```
makeExampleTinyEPICDataSet()
```

## Value

Object of class `SigSet`

## Examples

```
sset <- makeExampleTinyEPICDataSet()
```

---

**mapFileSet***Deposit data of one sample to a fileSet (and hence to file)*

---

**Description**

Deposit data of one sample to a fileSet (and hence to file)

**Usage**

```
mapFileSet(fset, sample, named_values)
```

**Arguments**

fset	a sesame::fileSet, as obtained via readFileSet
sample	sample name as a string
named_values	value vector named by probes

**Value**

a sesame::fileSet

**Examples**

```
## create two samples
fset <- initFileSet('mybetas2', 'HM27', c('s1','s2'))

## a hypothetical numeric array (can be beta values, intensities etc)
hypothetical <- setNames(runif(fset$n), fset$probes)

## map the numeric to file
mapFileSet(fset, 's1', hypothetical)

## get data
sliceFileSet(fset, 's1', 'cg00000292')
```

---

**meanIntensity***Mean Intensity*

---

**Description**

The function takes one single SigSet and computes mean intensity of all the in-band measurements. This includes all Type-I in-band measurements and all Type-II probe measurements. Both methylated and unmethylated alleles are considered. This function outputs a single numeric for the mean.

**Usage**

```
meanIntensity(sset)
```

**Arguments**

sset	a SigSet
------	----------

**Value**

mean of all intensities

**Examples**

```
sset <- makeExampleSeSAMeDataSet()
meanIntensity(sset)
```

MValueToBetaValue	<i>Convert M-value to beta-value</i>
-------------------	--------------------------------------

**Description**

Convert M-value to beta-value (aka inverse logit transform)

**Usage**

```
MValueToBetaValue(m)
```

**Arguments**

m	a vector of M values
---	----------------------

**Value**

a vector of beta values

**Examples**

```
MValueToBetaValue(c(-3, 0, 3))
```

noob	<i>Noob background correction</i>
------	-----------------------------------

**Description**

The function takes a SigSet and returns a modified SigSet with background subtracted. Background was modelled in a normal distribution and true signal in an exponential distribution. The Norm-Exp deconvolution is parameterized using Out-Of-Band (oob) probes

**Usage**

```
noob(sset, offset = 15)
```

**Arguments**

sset	a SigSet
offset	offset

**Value**

a new SigSet with noob background correction

**Examples**

```
sset <- makeExampleTinyEPICDataSet()
sset.nb <- noob(sset)
```

noobsb

*Background subtraction with bleeding-through subtraction*

**Description**

The function takes a SigSet and returns a modified SigSet with background subtracted. Signal bleed-through was modelled using a linear model with error estimated from cross-channel regression. Norm-Exp deconvolution using Out-Of-Band (oob) probes.

**Usage**

```
noobsb(sset, offset = 15, detailed = FALSE)
```

**Arguments**

sset	a SigSet
offset	offset
detailed	if TRUE, return a list of SigSet and regression function

**Value**

a modified SigSet with background correction

**Examples**

```
sset <- makeExampleSeSAMeDataSet('HM450')
sset.nb <- noobsb(sset)
```

oobG	<i>oobG getter generic</i>
------	----------------------------

### Description

*oobG* getter generic  
Get *oobG* slot of *SigSet* class

### Usage

```
oobG(x)

## S4 method for signature 'SigSet'
oobG(x)
```

### Arguments

x object of *SigSet*

### Value

The *oobG* slot of *SigSet*

### Examples

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
head(oobG(sset))
```

oobG<-	<i>oobG replacement generic</i>
--------	---------------------------------

### Description

*oobG* replacement generic  
Replace *oobG* slot of *SigSet* class

### Usage

```
oobG(x) <- value

## S4 replacement method for signature 'SigSet'
oobG(x) <- value
```

### Arguments

x object of *SigSet*  
value new value

**Value**

a new SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
df <- oobG(sset)
df[1,1] <- 10
oobG(sset) <- df
```

---

oobR

*oobR getter generic*

---

**Description**

oobR getter generic

Get oobR slot of SigSet class

**Usage**

```
oobR(x)

## S4 method for signature 'SigSet'
oobR(x)
```

**Arguments**

x object of SigSet

**Value**

The oobR slot of SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
head(oobR(sset))
```

---

oobR<-

*oobR replacement generic*

---

**Description**

oobR replacement generic

Replace oobR slot of SigSet class

**Usage**

```
oobR(x) <- value

## S4 replacement method for signature 'SigSet'
oobR(x) <- value
```

**Arguments**

x	object of <code>SigSet</code>
value	new value

**Value**

a new `SigSet`

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
df <- oobR(sset)
df[1,1] <- 10
oobR(sset) <- df
```

---

openSesame

*The openSesame pipeline*

---

**Description**

This function is a simple wrapper of noob + nonlinear dye bias correction + pOOBAH masking.

**Usage**

```
openSesame(
  x,
  platform = "",
  manifest = NULL,
  what = "beta",
  BPPARAM = SerialParam(),
  ...
)
```

**Arguments**

x	<code>SigSet(s)</code> , <code>IDAT</code> prefix(es), <code>minfi GenomicRatioSet(s)</code> , or <code>RGChannelSet(s)</code>
platform	optional platform string
manifest	optional dynamic manifest
what	either <code>'sigset'</code> or <code>'beta'</code>
BPPARAM	get parallel with <code>MulticoreParam(n)</code>
...	parameters to <code>getBetas</code>

## Details

If the input is an IDAT prefix or a SigSet, the output is the beta value numerics. If the input is a minfi GenomicRatioSet or RGChannelSet, the output is the sesamized GenomicRatioSet.

## Value

a numeric vector for processed beta values

## Examples

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
IDATprefixes <- searchIDATprefixes(
  system.file("extdata", "", package = "sesameData"))
betas <- openSesame(IDATprefixes)
```

`openSesameToFile`

*openSesame pipeline with file-backed storage*

## Description

openSesame pipeline with file-backed storage

## Usage

```
openSesameToFile(map_path, idat_dir, BPPARAM = SerialParam(), inc = 4)
```

## Arguments

<code>map_path</code>	path of file to be mapped (beta values file)
<code>idat_dir</code>	source IDAT directory
<code>BPPARAM</code>	get parallel with MulticoreParam(2)
<code>inc</code>	bytes per item data storage. increase to 8 if precision is important. Most cases 32-bit representation is enough.

## Value

a sesame::fileSet

## Examples

```
openSesameToFile('mybetas',
  system.file('extdata', package='sesameData'))
```

`parseGEOSignalABFile`    *Parse GEO signal-A/B File into a SigSet List*

## Description

This function is meant to be a convenience function for parsing data from Signal\_A and Signal\_B file provided by GEO. In many cases, this function generates a "partial" SigSet due to lack of out-of-band signal and control probe measurement in those Signal\_A/B files. The detection p-value is based on a fixed normal distribution rather than from negative control or OOB probes.

## Usage

```
parseGEOSignalABFile(path, platform = "HM450", drop = TRUE, parallel = TRUE)
```

## Arguments

path	path to Signal-A/B file downlaoded from GEO. The file can remain gzipped.
platform	HM450, EPIC or HM27
drop	whether to reduce to SigSet when there is only one sample.
parallel	whether to use multiple cores.

## Value

a SigSetList or a SigSet

## Examples

```
path = system.file(
  'extdata',
  'GSE36369_NonEBV_SignalA_SignalB_3samples_1k.txt.gz',
  package='sesame')
ssets <- parseGEOSignalABFile(path)
```

`predictAgeHorvath353`    *Horvath 353 age predictor*

## Description

The function takes a named numeric vector of beta values. The name attribute contains the probe ID (cg, ch or rs IDs). The function looks for overlapping probes and estimate age using Horvath aging model (Horvath 2013 Genome Biology). The function outputs a single numeric of age in years.

## Usage

```
predictAgeHorvath353(betas)
```

## Arguments

betas	a probeID-named vector of beta values
-------	---------------------------------------

**Value**

age in years

**Examples**

```
betas <- sesameDataGet('HM450.1.TCGA.PAAD')$betas
predictAgeHorvath353(betas)
```

---

predictAgePheno

*Phenotypic age predictor*

---

**Description**

The function takes a named numeric vector of beta values. The name attribute contains the probe ID (cg, ch or rs IDs). The function looks for overlapping probes and estimate age using Horvath aging model (Levine et al. 2018 Aging, 513 probes). The function outputs a single numeric of age in years.

**Usage**

```
predictAgePheno(betas)
```

**Arguments**

betas            a probeID-named vector of beta values

**Value**

age in years

**Examples**

```
betas <- sesameDataGet('HM450.1.TCGA.PAAD')$betas
predictAgePheno(betas)
```

---

predictAgeSkinBlood

*Horvath Skin and Blood age predictor*

---

**Description**

The function takes a named numeric vector of beta values. The name attribute contains the probe ID (cg, ch or rs IDs). The function looks for overlapping probes and estimate age using Horvath aging model (Horvath et al. 2018 Aging, 391 probes). The function outputs a single numeric of age in years.

**Usage**

```
predictAgeSkinBlood(betas)
```

**Arguments**

**betas** a probeID-named vector of beta values

**Value**

age in years

**Examples**

```
betas <- sesameDataGet('HM450.1.TCGA.PAAD')$betas
predictAgeSkinBlood(betas)
```

**print.fileSet** *Print a fileSet*

**Description**

Print a fileSet

**Usage**

```
## S3 method for class 'fileSet'
print(x, ...)
```

**Arguments**

**x** a sesame::fileSet  
**...** stuff for print

**Value**

string representation

**Examples**

```
fset <- initFileSet('mybetas2', 'HM27', c('s1','s2'))
fset
```

---

print.sesameQC	<i>Print sesameQC object</i>
----------------	------------------------------

---

**Description**

Print sesameQC object

**Usage**

```
## S3 method for class 'sesameQC'  
print(x, ...)
```

**Arguments**

x	a sesameQC object
...	extra parameter for print

**Value**

print sesameQC result on screen

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset  
sesameQC(sset)
```

---

---

probeNames	<i>Get Probe Names of SigSet class</i>
------------	--

---

**Description**

Get Probe Names of SigSet class

**Usage**

```
probeNames(x)  
  
## S4 method for signature 'SigSet'  
probeNames(x)
```

**Arguments**

x	object of Sigset
---	------------------

**Value**

A char vector

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset  
head(probeNames(sset))
```

**pval***pval getter generic***Description****pval** *getter generic*

Get pval slot of SigSet class

**Usage**

```
pval(x)

## S4 method for signature 'SigSet'
pval(x)
```

**Arguments**

<b>x</b>	object of SigSet
----------	------------------

**Value**

The pval slot of SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
head(pval(sset))
```

**pval<-***pval replacement generic***Description****pval** *replacement generic*

Replace pval slot of SigSet class

**Usage**

```
pval(x) <- value

## S4 replacement method for signature 'SigSet'
pval(x) <- value
```

**Arguments**

<b>x</b>	object of SigSet
<b>value</b>	new value

**Value**

a new SigSet

**Examples**

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
df <- pval(sset)
df[1] <- 0.01
pval(sset) <- list(p00BAH=df)
```

---

readFileSet

*Read an existing fileSet from storage*

---

**Description**

This function only reads the meta-data.

**Usage**

```
readFileSet(map_path)
```

**Arguments**

map_path	path of file to map (should contain valid _idx.rds index)
----------	---

**Value**

a sesame::fileSet object

**Examples**

```
## create two samples
fset <- initFileSet('mybetas2', 'HM27', c('s1','s2'))

## a hypothetical numeric array (can be beta values, intensities etc)
hypothetical <- setNames(runif(fset$n), fset$probes)

## map the numeric to file
mapFileSet(fset, 's1', hypothetical)

## read it from file
fset <- readFileSet('mybetas2')

## get data
sliceFileSet(fset, 's1', 'cg00000292')
```

<code>readIDATpair</code>	<i>Import a pair of IDATs from one sample</i>
---------------------------	---

## Description

The function takes a prefix string that are shared with \_Grn.idat and \_Red.idat. The function returns a SigSet.

## Usage

```
readIDATpair(
  prefix.path,
  platform = "",
  manifest = NULL,
  controls = NULL,
  readNBeads = FALSE,
  verbose = FALSE
)
```

## Arguments

prefix.path	sample prefix without _Grn.idat and _Red.idat
platform	EPIIC, HM450 and HM27 etc.
manifest	optional design manifest file
controls	optional control probe manifest file
readNBeads	whether to read number of beads
verbose	be verbose? (FALSE)

## Value

a SigSet

## Examples

```
sset <- readIDATpair(sub('_Grn.idat',' ',system.file(
  "extdata", "4207113116_A_Grn.idat", package = "sesameData")))
```

<code>reopenSesame</code>	<i>re-compute beta value for GenomicRatioSet</i>
---------------------------	--

## Description

re-compute beta value for GenomicRatioSet

## Usage

```
reopenSesame(x, naFrac = 0.2)
```

**Arguments**

x	GenomicRatioSet
naFrac	maximum NA fraction for a probe before it gets dropped (1)

**Value**

a GenomicRatioSet

RGChannelSetToSigSets *Convert RGChannelSet (minfi) to a list of SigSet (SeSAmE)*

**Description**

Notice the colData() and rowData() is lost. Most cases, rowData is empty anyway.

**Usage**

```
RGChannelSetToSigSets(rgSet, BPPARAM = SerialParam())
```

**Arguments**

rgSet	a minfi::RGChannelSet
BPPARAM	get parallel with MulticoreParam(n)

**Value**

a list of sesame::SigSet

**Examples**

```
if (require(FlowSorted.Blood.450k)) {
  rgSet <- FlowSorted.Blood.450k[,1:2]
  ssets <- RGChannelSetToSigSets(rgSet)
}
```

searchIDATprefixes *Identify IDATs from a directory*

**Description**

The input is the directory name as a string. The function identifies all the IDAT files under the directory. The function returns a vector of such IDAT prefixes under the directory.

**Usage**

```
searchIDATprefixes(dir.name, recursive = TRUE, use.basename = TRUE)
```

**Arguments**

- `dir.name` the directory containing the IDAT files.  
`recursive` search IDAT files recursively  
`use.basename` basename of each IDAT path is used as sample name This won't work in rare situation where there are duplicate IDAT files.

**Value**

the IDAT prefixes (a vector of character strings).

**Examples**

```
## only search what are directly under
IDATprefixes <- searchIDATprefixes(
  system.file("extdata", "", package = "sesameData"))

## search files recursively is by default
IDATprefixes <- searchIDATprefixes(
  system.file(package = "sesameData"), recursive=TRUE)
```

---

`segmentBins`

*Segment bins using DNAcopy*

---

**Description**

Segment bins using DNAcopy

**Usage**

```
segmentBins(bin.signals, bin.coords)
```

**Arguments**

- `bin.signals` bin signals (input)  
`bin.coords` bin coordinates

**Value**

segment signal data frame

sesameQC

*Generate summary numbers that indicative of experiment quality***Description**

Generate summary numbers that indicative of experiment quality

**Usage**

```
sesameQC(sset, betas = NULL)
```

**Arguments**

sset	a SigSet object
betas	processed beta values

**Value**

a sesameQC class object

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
sesameQC(sset)
```

sesamize

*"fix" an RGChannelSet (for which IDATs may be unavailable) with Sesame The input is an RGSet and the output is a sesamized minfi::GenomicRatioSet*

**Description**

"fix" an RGChannelSet (for which IDATs may be unavailable) with Sesame The input is an RGSet and the output is a sesamized minfi::GenomicRatioSet

**Usage**

```
sesamize(rgSet, naFrac = 1, BPPARAM = SerialParam(), HDF5 = NULL)
```

**Arguments**

rgSet	an RGChannelSet, perhaps with colData of various flavors
naFrac	maximum NA fraction for a probe before it gets dropped (1)
BPPARAM	get parallel with MulticoreParam(n)
HDF5	is the rgSet HDF5-backed? if so, avoid eating RAM (perhaps)

**Value**

a sesamized GenomicRatioSet

`show, SigSet-method`      *The display method for SigSet*

### Description

The function outputs the number of probes in each category and the first few signal measurements. NBeads slots are not shown here.

### Usage

```
## S4 method for signature 'SigSet'
show(object)
```

### Arguments

object	displayed object
--------	------------------

### Value

message of number of probes in each category.

### Examples

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
print(sset)
```

*SigSet-class*      *SigSet class*

### Description

This is the main data class for SeSAMe. The class holds different classes of signal intensities.

The function takes a string describing the platform of the data. It can be one of "HM27", "HM450" or "EPIC".

The function takes a string describing the platform of the data. It can be one of "HM27", "HM450" or "EPIC".

### Usage

```
## S4 method for signature 'SigSet'
initialize(.Object, platform, ...)
SigSet(...)
```

### Arguments

.Object	target object
platform	"EPIC", "HM450", "HM27" or other strings for custom arrays
...	additional arguments

**Details**

The NBeads\* slots are normally left empty but can be optionally turned on.

**Value**

a SigSet object  
a SigSet object

**Slots**

IG intensity table for type I probes in green channel  
IR intensity table for type I probes in red channel  
IGG Type-I green that is inferred to be green  
IRR Type-I red that is inferred to be red  
II intensity table for type II probes  
oobG out-of-band probes in green channel  
oobR out-of-band probes in red channel  
NBeadsIG Number of Beads for Infinium I green channel  
NBeadsIR Number of Beads for Infinium I red channel  
NBeadsII Number of Beads for Infinium II  
ctl all the control probe intensities  
pval named numeric vector of p-values  
platform "EPIC", "HM450" or "HM27"

**Examples**

```
## Create an empty EPIC object.  
SigSet("EPIC")  
SigSet('EPIC')
```

---

SigSetList                   *constructor*

---

**Description**

constructor

**Usage**

```
SigSetList(...)
```

**Arguments**

...                           the SigSet objects that will be the List elements

**Value**

a SigSetList

**Examples**

```
sset1 <- readIDATpair(file.path(system.file(
  'extdata', '', package='sesameData'), '4207113116_A'))

sset2 <- readIDATpair(file.path(system.file(
  'extdata', '', package='sesameData'), '4207113116_B'))

SigSetList(sset1, sset2)
```

**SigSetList-class***a List of SigSets with some methods of its own***Description***a List of SigSets with some methods of its own***SigSetList-methods***SigSetList methods (centralized). Currently scarce... ‘show’ print a summary of the SigSetList.***Description**

*SigSetList methods (centralized). Currently scarce...  
 ‘show’ print a summary of the SigSetList.*

**Usage**

```
## S4 method for signature 'SigSetList'
show(object)
```

**Arguments**

object	a SigSetList
--------	--------------

**Value**

*Description of SigSetList*

**Examples**

```
SigSetListFromPath(system.file("extdata", "", package = "sesameData"))
```

---

`SigSetListFromIDATs`    *read IDATs into a SigSetList*

---

**Description**

FIXME: switch from ‘parallel’ to BiocParallel

**Usage**

```
SigSetListFromIDATs(stubs, parallel = FALSE)
```

**Arguments**

<code>stubs</code>	the IDAT filename stubs
<code>parallel</code>	run in parallel? (default FALSE)

**Value**

a `SigSetList`

**Examples**

```
## a SigSetList of length 1
ssets <- SigSetListFromIDATs(file.path(
  system.file("extdata", "", package = "sesameData"), "4207113116_A"))
```

---

`SigSetListFromPath`    *read an entire directory’s worth of IDATs into a SigSetList*

---

**Description**

read an entire directory’s worth of IDATs into a `SigSetList`

**Usage**

```
SigSetListFromPath(path = ".", parallel = FALSE, recursive = TRUE)
```

**Arguments**

<code>path</code>	the path from which to read IDATs (default “.”)
<code>parallel</code>	run in parallel? (default FALSE)
<code>recursive</code>	whether to search recursively

**Value**

a `SigSetList`

**Examples**

```
## Load all IDATs from directory
ssets <- SigSetListFromPath(
  system.file("extdata", "", package = "sesameData"))
```

**SigSetsToRGChannelSet** *Convert sesame::SigSet to minfi::RGChannelSet*

### Description

Convert sesame::SigSet to minfi::RGChannelSet

### Usage

```
SigSetsToRGChannelSet(ssets, BPPARAM = SerialParam(), annotation = NA)
```

### Arguments

ssets	a list of sesame::SigSet
BPPARAM	get parallel with MulticoreParam(n)
annotation	the minfi annotation string, guessed if not given

### Value

a minfi::RGChannelSet

### Examples

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
rgSet <- SigSetsToRGChannelSet(sset)
```

**SigSetToRatioSet** *Convert one sesame::SigSet to minfi::RatioSet*

### Description

Convert one sesame::SigSet to minfi::RatioSet

### Usage

```
SigSetToRatioSet(sset, annotation = NA)
```

### Arguments

sset	a sesame::SigSet
annotation	minfi annotation string

### Value

a minfi::RatioSet

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
ratioSet <- SigSetToRatioSet(sset)
```

**sliceFileSet***Slice a fileSet with samples and probes***Description**

Slice a fileSet with samples and probes

**Usage**

```
sliceFileSet(fset, samples = fset$samples, probes = fset$probes, memmax = 10^5)
```

**Arguments**

fset	a sesame::fileSet, as obtained via readFileSet
samples	samples to query (default to all samples)
probes	probes to query (default to all probes)
memmax	maximum items to read from file to memory, to protect from accidental memory congestion.

**Value**

a numeric matrix of length(samples) columns and length(probes) rows

**Examples**

```
## create two samples
fset <- initFileSet('mybetas2', 'HM27', c('s1', 's2'))

## a hypothetical numeric array (can be beta values, intensities etc)
hypothetical <- setNames(runif(fset$n), fset$probes)

## map the numeric to file
mapFileSet(fset, 's1', hypothetical)

## get data
sliceFileSet(fset, 's1', 'cg00000292')
```

SNPcheck

*Check sample identity using SNP probes***Description**

Check sample identity using SNP probes

**Usage**

```
SNPcheck(betas)
```

**Arguments**

betas	numeric matrix (row: probes, column: samples)
-------	---

**Value**

grid object plotting SNP clustering

**Examples**

```
betas <- sesameDataGet('HM450.10.TCGA.PAAD.normal')
SNPcheck(betas)
```

subsetSignal

*Select a subset of probes***Description**

The function takes a SigSet as input and output another SigSet with probes from the given probe selection.

**Usage**

```
subsetSignal(sset, probes)
```

**Arguments**

sset	a SigSet
probes	target probes

**Value**

another sset with probes specified

**Examples**

```
sset <- sesameDataGet('EPIC.1.LNCaP')$sset
subsetSignal(sset, rownames(slot(sset, 'IR')))
```

---

topLoci	<i>Top loci in differential methylation</i>
---------	---

---

### Description

This is a convenience function to show top differential methylated segments. The function takes a coefficient table as input and output the same table ordered by the significance of the locus.

### Usage

```
topLoci(cf1)
```

### Arguments

cf1 coefficient table of one factor from diffMeth

### Value

coefficient table ordered by p-value of each locus

### Examples

```
data <- sesameDataGet('HM450.76.TCGA.matched')
cf <- DMR(data$betas, data$sampleInfo, ~type)
topLoci(cf[[1]])
```

---

topSegments	<i>Top segments in differential methylation</i>
-------------	---

---

### Description

This is a utility function to show top differential methylated segments. The function takes a coefficient table as input and output the same table ordered by the significance of the segments.

### Usage

```
topSegments(cf1)
```

### Arguments

cf1 coefficient table of one factor from DMR

### Value

coefficient table ordered by adjusted p-value of segments

### Examples

```
data <- sesameDataGet('HM450.76.TCGA.matched')
cf <- DMR(data$betas, data$sampleInfo, ~type)
topSegments(cf[[1]])
```

**totalIntensities**      *M+U Intensities for All Probes*

### Description

The function takes one single SigSet and computes total intensity of all the in-band measurements by summing methylated and unmethylated alleles. This function outputs a single numeric for the mean.

### Usage

```
totalIntensities(sset)
```

### Arguments

sset            a SigSet

### Value

a vector of M+U signal for each probe

### Examples

```
sset <- makeExampleSeSAMeDataSet()
intensities <- totalIntensities(sset)
```

**totalIntensityZscore**    *Calculate intensity Z-score*

### Description

This function compute intensity Z-score with respect to the mean. Log10 transformation is done first. Probes of each design type are grouped before Z-scores are computed.

### Usage

```
totalIntensityZscore(sset)
```

### Arguments

sset            a SigSet

### Value

a vector of Z-score for each probe

### Examples

```
sset <- sesameDataGet('HM450.1.TCGA.PAAD')$sset
head(totalIntensityZscore(sset))
```

---

<code>twoCompsEst2</code>	<i>Estimate the fraction of the 2nd component in a 2-component mixture</i>
---------------------------	--

---

## Description

Estimate the fraction of the 2nd component in a 2-component mixture

## Usage

```
twoCompsEst2(
  pop1,
  pop2,
  target,
  use.ave = TRUE,
  diff_1m2u = NULL,
  diff_1u2m = NULL
)
```

## Arguments

<code>pop1</code>	Reference methylation level matrix for population 1
<code>pop2</code>	Reference methylation level matrix for population 2
<code>target</code>	Target methylation level matrix to be analyzed
<code>use.ave</code>	use population average in selecting differentially methylated probes
<code>diff_1m2u</code>	A vector of differentially methylated probes (methylated in population 1 but unmethylated in population 2)
<code>diff_1u2m</code>	A vector of differentially methylated probes (unmethylated in population 1 but methylated in population 2)

## Value

Estimate of the 2nd component in the 2-component mixture

---

<code>visualizeGene</code>	<i>Visualize Gene</i>
----------------------------	-----------------------

---

## Description

Visualize the beta value in heatmaps for a given gene. The function takes a gene name which is taken from the UCSC refGene. It searches all the transcripts for the given gene and optionally extend the span by certain number of base pairs. The function also takes a beta value matrix with sample names on the columns and probe names on the rows. The function can also work on different genome builds (default to hg38, can be hg19).

**Usage**

```
visualizeGene(
  geneName,
  betas,
  platform = c("EPIC", "HM450"),
  upstream = 2000,
  dwstream = 2000,
  refversion = c("hg38", "hg19"),
  ...
)
```

**Arguments**

geneName	gene name
betas	beta value matrix (row: probes, column: samples)
platform	HM450 or EPIC (default)
upstream	distance to extend upstream
dwstream	distance to extend downstream
refversion	hg19 or hg38 (default)
...	additional options, see visualizeRegion

**Value**

None

**Examples**

```
betas <- sesameDataGet('HM450.76.TCGA.matched')$betas
visualizeGene('ADA', betas, 'HM450')
```

**visualizeProbes**

*Visualize Region that Contains the Specified Probes*

**Description**

Visualize the beta value in heatmaps for the genomic region containing specified probes. The function works only if specified probes can be spanned by a single genomic region. The region can cover more probes than specified. Hence the plotting heatmap may encompass more probes. The function takes as input a string vector of probe IDs (cg/ch/rs-numbers). If draw is FALSE, the function returns the subset beta value matrix otherwise it returns the grid graphics object.

**Usage**

```
visualizeProbes(
  probeNames,
  betas,
  platform = c("EPIC", "HM450"),
  refversion = c("hg38", "hg19"),
  upstream = 1000,
  dwstream = 1000,
  ...
)
```

**Arguments**

probeNames	probe names
betas	beta value matrix (row: probes, column: samples)
platform	HM450 or EPIC (default)
refversion	hg19 or hg38 (default)
upstream	distance to extend upstream
dwstream	distance to extend downstream
...	additional options, see visualizeRegion

**Value**

None

**Examples**

```
betas <- sesameDataGet('HM450.76.TCGA.matched')$betas
visualizeProbes(c('cg22316575', 'cg16084772', 'cg20622019'), betas, 'HM450')
```

---

visualizeRegion	<i>Visualize Region</i>
-----------------	-------------------------

---

**Description**

The function takes a genomic coordinate (chromosome, start and end) and a beta value matrix (probes on the row and samples on the column). It plots the beta values as a heatmap for all probes falling into the genomic region. If ‘draw=TRUE’ the function returns the plotted grid graphics object. Otherwise, the selected beta value matrix is returned. ‘cluster.samples=TRUE/FALSE’ controls whether hierarchical clustering is applied to the subset beta value matrix.

**Usage**

```
visualizeRegion(
  chrm,
  plt.beg,
  plt.end,
  betas,
  platform = c("EPIC", "HM450"),
  refversion = c("hg38", "hg19"),
  sample.name.fontsize = 10,
  heat.height = NULL,
  draw = TRUE,
  show.sampleNames = TRUE,
  show.samples.n = NULL,
  show.probeNames = TRUE,
  cluster.samples = FALSE,
  nprobes.max = 1000,
  na.rm = FALSE,
  dmin = 0,
  dmax = 1
)
```

**Arguments**

chr	chromosome
plt.beg	begin of the region
plt.end	end of the region
betas	beta value matrix (row: probes, column: samples)
platform	EPIC or HM450
refversion	hg38 or hg19
sample.name.fontsize	sample name font size
heat.height	heatmap height (auto inferred based on rows)
draw	draw figure or return betas
show.sampleNames	whether to show sample names
show.samples.n	number of samples to show (default: all)
show.probeNames	whether to show probe names
cluster.samples	whether to cluster samples
nprobes.max	maximum number of probes to plot
na.rm	remove probes with all NA.
dmin	data min
dmax	data max

**Value**

graphics or a matrix containing the captured beta values

**Examples**

```
betas <- sesameDataGet('HM450.76.TCGA.matched')$betas
visualizeRegion('chr20', 44648623, 44652152, betas, 'HM450')
```

**visualizeSegments**      *Visualize segments*

**Description**

The function takes a CNSegment object obtained from cnSegmentation and plot the bin signals and segments (as horizontal lines).

**Usage**

```
visualizeSegments(seg, to.plot = NULL)
```

**Arguments**

seg	a CNSegment object
to.plot	chromosome to plot (by default plot all chromosomes)

**Details**

```
require ggplot2, scales
```

**Value**

```
plot graphics
```

**Examples**

```
seg <- sesameDataGet('EPIC.1.LNCaP')$seg  
visualizeSegments(seg)
```

# Index

\* **DNA**Methylation  
    sesame-package, 4

\* **Microarray**  
    sesame-package, 4

\* **QualityControl**  
    sesame-package, 4

\_PACKAGE (sesame-package), 4

as.data.frame.sesameQC, 4

BetaValueToMValue, 5

binSignals, 6

bisConversionControl, 6

buildControlMatrix450k, 7

chipAddressToSignal, 7

cnSegmentation, 8

ctl, 9

ctl, SigSet-method (ctl), 9

ctl<-, 9

ctl<-, SigSet-method (ctl<-), 9

detectionPfixedNorm, 10

detectionPnegEcdf, 11

detectionPnegNorm, 11

detectionPnegNormGS, 12

detectionPnegNormTotal, 13

detectionPoobEcdf, 13

detectionZero, 14

diffRefSet, 15

DML, 15

DMR, 16

dyeBiasCorr, 17

dyeBiasCorrMostBalanced, 18

dyeBiasCorrTypeINorm, 18

estimateCellComposition, 19

estimateLeukocyte, 19

getAFTypeIbySumAlleles, 20

getBetas, 21

getBinCoordinates, 22

getNormCtls, 22

getProbesByGene, 23

getProbesByRegion, 23

getProbesByTSS, 24

getRefSet, 25

getSegment, 26

getSexInfo, 26

IG, 27

IG, SigSet-method (IG), 27

IG<-, 27

IG<-, SigSet-method (IG<-), 27

II, 28

II, SigSet-method (II), 28

II<-, 28

II<-, SigSet-method (II<-), 28

inferEthnicity, 29

inferSex, 30

inferSexKaryotypes, 30

inferTypeIChannel, 31

initFileSet, 32

initialize, SigSet-method  
    (SigSet-class), 52

IR, 32

IR, SigSet-method (IR), 32

IR<-, 33

IR<-, SigSet-method (IR<-), 33

makeExampleSeSAMeDataSet, 34

makeExampleTinyEPICDataSet, 34

mapFileSet, 35

meanIntensity, 35

MValueToBetaValue, 36

noob, 36

noobsb, 37

oobG, 38

oobG, SigSet-method (oobG), 38

oobG<-, 38

oobG<-, SigSet-method (oobG<-), 38

oobR, 39

oobR, SigSet-method (oobR), 39

oobR<-, 39

oobR<-, SigSet-method (oobR<-), 39

openSesame, 40

openSesameToFile, 41

parseGEOSignalABFile, 42  
p00BAH (detectionPoobEcdf), 13  
predictAgeHorvath353, 42  
predictAgePheno, 43  
predictAgeSkinBlood, 43  
print.fileSet, 44  
print.sesameQC, 45  
probeNames, 45  
probeNames, SigSet-method (probeNames),  
    45  
pval, 46  
pval, SigSet-method (pval), 46  
pval<-, 46  
pval<-, SigSet-method (pval<-), 46  
  
readFileSet, 47  
readIDATpair, 48  
reopenSesame, 48  
RGChannelSetToSigSets, 49  
  
searchIDATprefixes, 49  
segmentBins, 50  
sesame (sesame-package), 4  
sesame-package, 4  
sesameQC, 51  
sesamize, 51  
show, SigSet-method, 52  
show, SigSetList-method  
    (SigSetList-methods), 54  
SigSet (SigSet-class), 52  
SigSet-class, 52  
SigSetList, 53  
SigSetList-class, 54  
SigSetList-methods, 54  
SigSetListFromIDATs, 55  
SigSetListFromPath, 55  
SigSetsToRGChannelSet, 56  
SigSetToRatioSet, 56  
sliceFileSet, 57  
SNPcheck, 58  
subsetSignal, 58  
  
topLoci, 59  
topSegments, 59  
totalIntensities, 60  
totalIntensityZscore, 60  
twoCompsEst2, 61  
  
visualizeGene, 61  
visualizeProbes, 62  
visualizeRegion, 63  
visualizeSegments, 64