oneChannelGUI

November 11, 2009

R topics documented:

1	3
biomartFilter	3
buildingLocalAnnotation	4
	4
combineGeoMSF	5
consistentFilters	5
cosieHscore	6
cosieHsfull	6
cosieMmcore	7
cosieMmfull	7
cosieWrapper	8
createGeoTarget	8
	9
crosshybhuex.annotation	9
crosshybmoex.annotation	0
crosshybraex.annotation	1
dfMAplot	1
EG2probeset	2
erankProdAltSplFilter	2
erankProdAltSpl	3
exonContrasts	4
exonsSpecific2as	4
exonTopTableExtract	5
extractAffyids	5
filteringTable	6
geneExonLibs	6
geoVSbioc	7
GOenrichment	7
huex.annotation	8
huex.variantexons	8
hugene.annotation	9
inspecting.one.splice.index	9
inspecting.splice.index	0
intensityFilter	0
IPAlistFilter	1
iqrFilter	1
limma Evons	2

49

Index

listFilter	22
	23
11 6	23
11 6	24
	24
Ci C	25
CI	25
Ci	26
•	26
	27
ϵ	27
moex.annotation	28
moex.variantexons	28
mogene.annotation	29
myExpresso	29
normBoxplot	30
ocPlotHist	30
ocPlotPCA	30
oneChannelGUI-package	31
oneChannelGUI	31
	33
OpenLargefiles	34
plierToZero	34
	35
PlotOptionsv1	35
plotVariantSI	36
raex.annotation	36
raex.variantexons	37
ragene.annotation	37
	38
	38
	39
1	39
•	40
•	40
•	40
<u>•</u>	41
	42
sample.size.evaluation	42
1	43
	43
	44
••	44
	45
1	45
	46
	46
updateLibs	47
variantExons	47
variantSI	47
VennDiagram	48
	10

AptMidas 3

AptMidas

Graphical interface to APT midas

Description

This is a graphical interface to the midas program for detection of alternative splicing detection implemented in the Affymetrix APT tools

Usage

```
AptMidas()
```

Note

For more information see Affymetrix Alternative Transcript Analysis Methods for Exon Arrays whitepaper. Before using MiDAS is strongly recomanded to filter out gene level probe sets with low intensity values to avoid searching alternative splicing for probe sets which are not expressed. This can be done using filtering method implemented in oneChannelGUI which define a background intensity threshold on the basis of the intron exon signals of a set of housekeeping genes present in the exon arrays. However it is also possible to use a filter based on the dabg p-value calculated using Affymetrix APT tools. This function will also calculate Splice Index

Author(s)

Raffaele A Calogero

See Also

erankProdAltSpl

biomartFilter

Filtering only gene-level probe sets with multiple ensembl transcripts

Description

This function allows to filter exon array data to selected only those gene-level probe sets associated to multiple mRNAs annotated in ensembl data base

Usage

```
biomartFilter()
```

Author(s)

4 colExtract

buildingLocalAnnotation

Updates local gene-level annotation data for gene and exon arrays using the netaffx database

Description

Internal oneChannelGUI Gene/Exon gene-level annotation data can be upgraded using this function, which queries netaffx database. annotation RDA files are saved in the data subdir of oneChannel-GUI dir. Windows users need to exchange the older copies present in Rdata.zip, simply dragging them in the zip file.

Usage

```
buildingLocalAnnotation()
```

Author(s)

Raffaele A Calogero

colExtract

Extract a column from a tab delimited file with header

Description

This function allows to extract a specific column from a tab delimited file generated by oneChannelGUI. The file should contain an header. This function is useful to extract probe set ids to be used fro ven diagram representations

Usage

```
colExtract()
```

Author(s)

combineGeoMSF 5

combineGeoMSF	This function allows to combine GEO Matrix Series Files belonging to
	the same experiment.

Description

The function combines in a unique ExpressionSet the data derived from multiple Matrix Series Files belongig to a GEO experiment containing more than 255 arrays.

Usage

```
combineGeoMSF()
```

Note

see oneChannelGUI vignette for more info

Author(s)

Raffaele A Calogero

consistentFilters This function allows filtering using the combination of multiple paramenters, e.g. MiDAS p-values and Rank Product p-values

Description

This filter can be used to moderate multiple tests errors. E.g. finding the intersection between MiDAS p-values and Rank Product p-values user will remove some of the false positive produced by the two methods. A filter on the size of delta Splice Index associated to MiDAS p-values filter will will allow to remove statistical significant splicing events which are characterized by a very limited variation.

Usage

```
consistentFilters()
```

Note

This fuction needs the presence of Splice Index data, MiDAS p-values and RP p-values. It works for two groups only

Author(s)

Raffaele A Calogero

See Also

erankProdAltSpl, AptMidas

6 cosieHsfull

cosieHscore Parameters precalculate

Parameters precalculated for the Affymetrix core human annotation subgroup from FMI Bioinformatics group and needed to correct SI

Description

The description of how these data were generated is described in Gaidatzis et al. Nucleic Acids Research, 2009, 1

Usage

cosieHscore

Format

A matrix with 7 observations: status, geneTrainFrom, geneTrainTo, a, b, c, d

References

www.fmi.ch/groups/gbioinfo/cosie/cosie.html

cosieHsfull

Parameters precalculated for the Affymetrix full human annotation subgroup from FMI Bioinformatics group and needed to correct SI

Description

The description of how these data were generated is described in Gaidatzis et al. Nucleic Acids Research, 2009, 1. Authors recommend the use of the core set since the full set is not as well characterized.

Usage

cosieHsfull

Format

A matrix with 7 observations: status, geneTrainFrom, geneTrainTo, a, b, c, d

References

www.fmi.ch/groups/gbioinfo/cosie/cosie.html

cosieMmcore 7

cosieMmcore	Parameters precalculated for the Affymetrix core mouse annotation
	subgroup from FMI Bioinformatics group and needed to correct SI

Description

The description of how these data were generated is described in Gaidatzis et al. Nucleic Acids Research, 2009, 1

Usage

cosieMmcore

Format

A matrix with 7 observations: status, geneTrainFrom, geneTrainTo, a, b, c, d

References

www.fmi.ch/groups/gbioinfo/cosie/cosie.html

cosieMmfull Parameters precalculated for the Affymetrix full mouse annotation subgroup from FMI Bioinformatics group and needed to correct SI

Description

The description of how these data were generated is described in Gaidatzis et al. Nucleic Acids Research, 2009, 1. Authors recommend the use of the core set since the full set is not as well characterized.

Usage

cosieMmfull

Format

A matrix with 7 observations: status, geneTrainFrom, geneTrainTo, a, b, c, d

References

www.fmi.ch/groups/gbioinfo/cosie/cosie.html

8 createGeoTarget

cosieWrapper	A wrapper to cosie, Corrected Splicing Indices for Exon arrays, that
	for any given set of new exon array experiments corrects for the ob-
	served bias and improves the detection of alternative splicing

Description

COSIE is a function that for a given set of exon arrays corrects for the observed bias and improves the detection of alternative splicing. It adjusts splicing indices for exons, especially for those that belong to differentially expressed genes. For this adjustment, COSIE uses parameters that are specific for each probeset which were trained from a large number of published exon arrays. The downside of this approach is that such parameters cannot be estimated for all probesets on the microarray. Based on our training set, COSIE corrects 95.1 percent of the probesets. Separate parameter files are provided for both the full and core sets, including all probesets that are linked to transcripts. We recommend the use of the core set that was also used in the cited study below. The full set is not as well characterized.

Usage

```
cosieWrapper()
```

Author(s)

Gaidatzis et al. Nucleic Acids Research, 2009, 1

See Also

cosieHscore, cosieMmcore, cosieHsfull, cosieMmfull

createGeoTarget

Creating a affylmGUI Target starting from a GEO Matrix series file

Description

The function extracts from GEO series matrix file all the information to create a Target file, that can be used to load the GEO series matrix file into oneChannelGUI.

Usage

```
createGeoTarget()
```

Note

see oneChannelGUI vignette for Target file description

Author(s)

crosshybFilter 9

crosshybFilter	Removing from exon array gene/exon level probe sets characterized by cross hybridization with other transcripts

Description

XHYB field is mainly an indicator of weak assignment between a transcript cluster and the assigned mRNA, suggesting a potential crosshyb, CRSSHYB is a measure of the promiscuity of the probes within a probe set among transcribed sequences.

- 1. 1 unique. All probes in the probe set perfectly match only one sequence in the putatively transcribed array design content. The vast majority of probe sets are unique.
- 2. 2 similar. All the probes in the probe set perfectly match more than one sequence in the putatively transcribed array design content.
- 3. 3 mixed. The probes in the probe set either perfectly match or partially match more than one sequence in the putatively transcribed array design content.

XHYB and CRSSHYB are used to remove probe sets characterized by multiple hybridization of exon probes

Usage

```
crosshybFilter()
```

Author(s)

Raffaele A Calogero

crosshybhuex.annotation

Cross hybriduzation data for exon CORE subset of human exon array 1.0 ST

Description

These data are derived from Affymetrix annotation file huex10stv2na23hg18. XHYB field is mainly an indicator of weak assignment between a transcript cluster and the assigned mRNA, suggesting a potential crosshyb, CRSSHYB is a measure of the promiscuity of the probes within a probe set among transcribed sequences.

- 1. 1 unique. All probes in the probe set perfectly match only one sequence in the putatively transcribed array design content. The vast majority of probe sets are unique.
- 2. 2 similar. All the probes in the probe set perfectly match more than one sequence in the putatively transcribed array design content.
- 3. 3 mixed. The probes in the probe set either perfectly match or partially match more than one sequence in the putatively transcribed array design content.

XHYB and CRSSHYB are used to remove probe sets characterized by multiple hybridization of exon probes Cross-hybridization potential of the probe set determined at the time of array design. This field is based on computational sequence alignment against all known and putatively transcribed array design content, which includes all potentially transcribed regions of the genome and other transcribed sequences that could not be mapped to the genome.

Usage

crosshybhuex.annotation

Format

A data frame with 9 observations: EPROBESETID, GPROBESETID, ACC, XHYB, CHR, STRAND, START, STOP, CROSSHYBTYPE

References

Affymetrix web site

crosshybmoex.annotation

Cross hybriduzation data for exon CORE subset of mouse exon array 1.0 ST

Description

These data are derived from Affymetrix annotation file moex10stv1na24mm8. XHYB field is mainly an indictor of weak assignment between a transcript cluster and the assigned mRNA, suggesting a potential crosshyb, CRSSHYB is a measure of the promiscuity of the probes within a probe set among transcribed sequences.

- 1. 1 unique. All probes in the probe set perfectly match only one sequence in the putatively transcribed array design content. The vast majority of probe sets are unique.
- 2. 2 similar. All the probes in the probe set perfectly match more than one sequence in the putatively transcribed array design content.
- 3. 3 mixed. The probes in the probe set either perfectly match or partially match more than one sequence in the putatively transcribed array design content.

XHYB and CRSSHYB are used to remove probe sets characterized by multiple hybridization of exon probes Cross-hybridization potential of the probe set determined at the time of array design. This field is based on computational sequence alignment against all known and putatively transcribed array design content, which includes all potentially transcribed regions of the genome and other transcribed sequences that could not be mapped to the genome.

Usage

crosshybmoex.annotation

Format

A data frame with 9 observations: EPROBESETID, GPROBESETID, ACC, XHYB, CHR, STRAND, START, STOP, CROSSHYBTYPE

References

Affymetrix web site

crosshybraex.annotation 11

crosshybraex.annotation

Cross hybridization data for exon CORE subset of rat exon array 1.0 ST

Description

These data are derived from Affymetrix annotation file raex10stv1na24rn4. XHYB field is mainly an indictor of weak assignment between a transcript cluster and the assigned mRNA, suggesting a potential crosshyb, CRSSHYB is a measure of the promiscuity of the probes within a probe set among transcribed sequences.

- 1. 1 unique. All probes in the probe set perfectly match only one sequence in the putatively transcribed array design content. The vast majority of probe sets are unique.
- 2. 2 similar. All the probes in the probe set perfectly match more than one sequence in the putatively transcribed array design content.
- 3. 3 mixed. The probes in the probe set either perfectly match or partially match more than one sequence in the putatively transcribed array design content.

XHYB and CRSSHYB are used to remove probe sets characterized by multiple hybridization of exon probes Cross-hybridization potential of the probe set determined at the time of array design. This field is based on computational sequence alignment against all known and putatively transcribed array design content, which includes all potentially transcribed regions of the genome and other transcribed sequences that could not be mapped to the genome.

Usage

crosshybraex.annotation

Format

A data frame with 9 observations: EPROBESETID, GPROBESETID, ACC, XHYB, CHR, STRAND, START, STOP, CROSSHYBTYPE

References

Affymetrix web site

dfMAplot

MA and Volcano plots from data present in a limma derived topTable

Description

MA and Volcano plots can be generated starting from limma results summarized in a topTable. Specific subsets of the topTable defined by p-value below an user-defined threshold and/or log2 fold changes over an user-defined threshold can be saved. The subset of data can be saved as a tab delimited file

Usage

12

```
dfMAplot(table1)
```

Arguments

table1

topTable data.frame generate by affylmGUI

Note

To know more about topTable see limma help

Author(s)

Raffaele A Calogero

EG2probeset

This function allows to link GeneBank and Entrez Gene ids to genelevel probe set ids

Description

This function allows to link oneChannelGUI embedded Affymetrix annotated accession numbers to gene-level probe set ids. Usnig the ACC EG are linked using the Bioconductor human, mouse or rat LLMappings annotation library

Usage

```
EG2probeset()
```

Author(s)

Raffaele A Calogero

```
erankProdAltSplFilter
```

Filtering Rank Product results for the detection of alternative splicing events

Description

This is a graphical interface to filter data on the basis of p-value generated by rank product analysis applied for the detection of alternative splicing

Usage

```
erankProdAltSpl()
```

Author(s)

erankProdAltSpl 13

See Also

erankProdAltSpl, AptMidas

erankProdAltSpl

Implementation of the Rank Product method for the detection of alternative splicing events

Description

This is a graphical interface to the RP function from RankProd package applied to detection of alternative splicing

Usage

```
erankProdAltSpl()
```

Details

Before using this method it is strongly suggested to perform a filter on the basis of DABG p-values using the filtering function available in the filtering menu. DABG values can be calculated if exon array probe set data are generated using the oneChannelGUI graphical implementation to APT tools. Affymetrix suggests to calculate probe set intensity at gene level using iterPlier and at exon level using plier. Subsequently SpliceIndex need to be calculated using the function available in the exon menu. Finally the Rank Product method could be applied exon by exon. For more details on the method see RankProd package. Selection of putative alternative splicing could be done using the filtering function available in the filtering menu of oneChannelGUI

Note

IMPORTANT we are still evaluating the efficacy of this method for detection of alternative splicing events. Use it being concious of this!

Author(s)

Raffaele A Calogero

See Also

inspecting.splice, spliceIndex

14 exonsSpecific2as

exonContrasts

Defining t-test regularized p-values

Description

This function constructs the contrasts as affylmGUI but applied to exon-level. It also perfoms Bayes regularization. Raw p-values are plotted to see if BH or BY type I error correction can be applied. Corrected p-values are saved and used for extraction of alternative spliced exons.

Usage

```
exonContrasts()
```

Details

It is important to note that if multiple contrasts should be considered after the calculation of each of them it is essential to extract the alternative spliced exons with exonTopTableExtract, because everytime exonContrasts is run it overwrites the previous results.

Author(s)

Raffaele A Calogero

exonsSpecific2as

Defining the exons associated to the various alternative isoforms

Description

This function uses the output derived from the function mapping2ensembl and produces a list of 1 and 0 for each of the alternative trasncripts associated to a specific Entrez Gene. This funciton is useful to define which splicing events are not associated to exons conserved over all the possible isoforms

Usage

```
exonsSpecific2as()
```

Author(s)

exonTopTableExtract 15

exonTopTableExtract

Extracts data on the basis of a defined t-test regularized p-value

Description

This function filters the data produced by exonContrasts to extract a list of alternative spliced exons that are saved in a file and they can be used for further analysis, i.e. extracting only variant exons. The function also filter the data present in the onechannelGUI project bot at gene and exon-level.

Usage

```
exonTopTableExtract()
```

Details

It is important to note that if multiple contrasts should be considered after the calculation of each of them it is essential to extract the alternative spliced exons with exonTopTableExtract, because everytime exonContrasts is run it overwrites the previous results.

Author(s)

Raffaele A Calogero

extractAffyids

Extracting probe ids associated to a specific Gene Ontology term

Description

It is possible to identify the affy ids associated to a specific GO term using the extractAffyids function.

Usage

```
extractAffyids()
```

Details

The function asks to the user to select a file containing probe set ids separated by carriage return. The file should contain only one column and no header. The user is also asked to select a specific GO term. The probe sets associated to the specific GO term will be annotated ans saved in a HTML file.

Note

For the annotation the annotation library associated to the raw data loaded in the affylmGUI environment is used.

Author(s)

16 geneExonLibs

See Also

GOenrichment, plotGO

filteringTable

Filtering a tab delimited file

Description

This function allows to filter a tab delimited file using a vector of data present in an other file. The two files should have an header and the column name to be used for the filtering should be equal in both files

Usage

```
filteringTable()
```

Author(s)

Raffaele A Calogero

geneExonLibs

Download the Library files for gene and exon analysis

Description

Affymetrix Gene/Exon library files are necessary to APT tools to calculate probe set summaries. The versions downloaded from www.bioinformatica.unito.it, with htis function, contain all informations needed to analyse gene exon arrays.

Usage

```
geneExonLibs()
```

Author(s)

geoVSbioc 17

geoVSbioc

linking GEO platforms to available BioC annotations libraries

Description

This data file gives the linke between GEO platforms and BioC annotation libraries. If the GEO BioC link exists the Bioconductor annotation lib is directly loaded in the annotation fild of the SespressionSet

Usage

geoVSbioc

Format

A data frame with 4 observations: GEOAcc, Organisms, Title, BiocAnLIb

References

GEO and Bioconductor

GOenrichment

Searching for Gene Ontology enriched terms within a set of differentially expressed genes

Description

In Bioconductor is available a library called GOstats, which allows the calculation of enriched GO terms within a set of differentially expressed probe sets. This is a graphical implementation of a function allowing the extraction of GO enriched term in a sub set of differentially expressed probe sets. To know more about it see GOstat library

Usage

```
GOenrichment()
```

Details

The function asks to the user to select a file containing probe set ids separated by carriage return. The file should contain only one column and no header. The set of enriched terms are plotted in red over the graph of all GO term associated to the differentially expressed genes. GO enriched terms can be also saved in a tab delimited file.

Author(s)

Raffaele A Calogero

References

Robert Gentleman GOstat package

18 huex.variantexons

See Also

extractAffyids, plotGO

huex.annotation

Annotation data for CORE subset of human exon array 1.0 ST

Description

These annotation data are derived from Affymetrix annotation file huex.annotation

Usage

huex.annotation

Format

A data frame with 5 observations: PROBESETID, SYMBOL, DESCRIPTION, CYTOBAND, ACC

References

Affymetrix web site

huex.variantexons table linking exon-level probe set ids to variant exons

Description

The table is made integrating exon-level netaffx annotaiton with UCSC derived annotation. Specifically exon-level rpobesets are mapped on variant exons, i.e. those that are specific of only a subset of all isoforms associated to a gene. This table is used to filter differentially expressed exons to select only those associated to variant isoforms

Usage

huex.variantexons

Format

A data frame with 12 observations: affyname, affystart, affyend, affywidth, affystrand, vspname, vspstart, vspend, vspwidth, vspstrand, chr, genome

References

Affymetrix web site and UCSC database

hugene.annotation 19

hugene.annotation Annotation data for human gene array 1.0 ST

Description

These annotation data are derived from Affymetrix annotation file hugene 10stv1na24hg18.

Usage

hugene.annotation

Format

A data frame with 5 observations: PROBESETID, SYMBOL, DESCRIPTION, CYTOBAND, ACC

References

Affymetrix web site

inspecting.one.splice.index

Plotting on the profiles of splice indexes for a transcript cluster ID

Description

This function plots the splice index profiles for one trnascript cluster ID

Usage

```
inspecting.one.splice.index( transcriptID, SpliceIndexDiff, PvalExon)
```

Arguments

transcriptID A character string indicating the geatureName of the transcript cluster to be evaluate

SpliceIndexDiff

a numerical value indicating the treshold to detect an alternative splicing. It represent the minimal absolute difference between the splice indexes measured for the same exon under two different experimental conditions

PvalExon

The max p-value obtainable by a t-test done on the splice indexes measured for the same exon under two different experimental conditions

Author(s)

Raffaele A Calogero

See Also

spliceIndex

20 intensityFilter

```
inspecting.splice.index
```

Plotting on a pdf file the profiles of splice indexes

Description

This function prints in a pdf file the splice index profiles of the available genes

Usage

```
inspecting.splice.index()
```

Author(s)

Raffaele A Calogero

See Also

spliceIndex

intensityFilter

intensity filtering with a mouse click

Description

This function removes all probe sets in which a certain percentage of experiments is below a user defined intensity threshold.

Usage

```
intensityFilter()
```

Details

The aim of non specific filtering is to remove the genes that are unlikely to carry information about the phenotypes under investigation. This filtering remove genes that do not have a centain level of, user defined, intensities in a set of, user defined, experiments.

Note

Factor analysis will be limited by the problem of having fewer samples than genes. Therefore, preselecting a smaller set of genes is definetively helpful.

Author(s)

Raffaele A Calogero

See Also

iqrFilter, listFilter, IPAlistFilter

IPAlistFilter 21

IPAlistFilter	Filtering an expression set using a set of Entrez genes extracted from Ingenuity Pathways analysis (IPA)
---------------	--

Description

It is possible to sub set an expression set loaded in the affylmGUI environment starting form a list of Entrez genes derived by IPA search tool.

Usage

```
IPAlistFilter()
```

Details

The function asks to the user to select a file containing Entrez genes separated by carriage return. The file should contain only one column and no header.

Author(s)

Raffaele A Calogero

See Also

iqrFilter, listFilter, intensityFilter

igrFilter

Interquantile filtering with a mouse click

Description

This function implements the interquantile filtering proposed by Heydebreck in 2004

Usage

```
iqrFilter()
```

Details

The aim of non specific filtering is to remove the genes that are unlikely to carry information about the phenotypes under investigation. This filtering remove genes that show little changes within the experimental points.

Note

Factor analysis will be limited by the problem of having fewer samples than genes. Therefore, preselecting a smaller set of genes is definetively helpful.

Author(s)

22 listFilter

References

Heydebreck et al. Bioconductor project Papers 2004

See Also

IPAlistFilter, listFilter, intensityFilter

limmaExons

graphical interface to limma for alternative splicing detection

Description

Applying the limma model fitting to eoxn-level data. Same implementation of AffylmGUI but applied to exon-level data. The first indicatin of alternative splicing detection using limma was proposed by Shah and Pallas in BMC Bioinformatics. 2009 Jan 20;10:26

Usage

```
limmaExons()
```

Details

The funciton fit the limma linear model to exon-level data

Author(s)

Raffaele A Calogero

listFilter

Subsetting an expression set using a list of Affymetrix ids

Description

This function subsets the normalized expression set present in the affylmGUI environment on the basis of a list of probe set ids passed via flat file.

Usage

```
listFilter()
```

Details

The function asks to the user to select a file containing probe set ids separated by carriage return. The file should contain only one column and no header.

Note

In transcriptional studies focusing on genes characterized by specific feature (i.e. transcription factor elements in promoters) the best filtering approach is selecting only those genes linked to the interesting biological feature.

makeBED15 23

Author(s)

Raffaele A Calogero

See Also

IPAlistFilter, iqrFilter, intensityFilter

makeBED15

This function creates files in BED15 format to be loaded on the UCSC browser

Description

This function creates files in BED15 format to be loaded on UCSC genome browser. The functin uses the data derived by variantSI filter on the basis of chromosome annotation

Usage

makeBED15()

Author(s)

Raffaele A Calogero

See Also

variantSI, variantExons, plotVariantSI

mapping2ensembl

Associating e-level probe sets to entrez gene exonic structure

Description

This function associates the statistical and expression data produced by a oneChannleGUI exonlevel analysis to the exonic structure of Entrez Gene ID. This function uses biomaRt to retrieve the sequence of EG exons. RRE database is instead used to retrieve the exon-level target sequences. Any exon-level probe set id to be associated to the EG exonic sequence need to be a perfct matching substring of the exon. In the otehr case no exon is associated to the probe set

Usage

```
mapping2ensembl()
```

Author(s)

24 mapping2RefSeq

rt- on

Description

This function retrieve from RRE the PSR sequences associated to the exon-level probe sets and all exons associated to the gene associated to PRS. Subsequently identify the exon where PSR maps and procuces a fasta file were are located exon-level PSR and target exon. The mapping is done using the countPattern function of the Biostrings package. Up to three mismathces are allowed in PSR mapping on exonic sequence.

Usage

```
mapping2exon()
```

Author(s)

Raffaele A Calogero

mapping2RefSeq

This function maps on NCBI Reference sequences spliced exons detected by the function one Channel GUI: Inspecting splice indexes

Description

This function retrieve from RRE the PSR sequences associated to the exon-level probe sets using blastn detects the best refseq associated to any of the exon-level probe sets retrieve from org.XX.eg.db the EG associated to any of the detected refseq and retrieves all the refseqs associated to the EG. Subsequently check if PSR maps on all the refseqs associated to the eg (conserved exon) or only some of them (isoform specific exon)

Usage

```
mapping2RefSeq()
```

Author(s)

masigpro 25

masigpro.edesign

The function creates an edesign object needed to run maSigPro

Description

The function creates an edesign object needed to run maSigPro. To know more about edesign object see maSigPro help. This function uses a specific configuration of Target column of the affylmGUI target file. To know more about target file see affylmGUI help. Each row of the column named Target, in the affylmGUI target file, describes the array on the basis of the experimental design. Each element needed for the construction of edesign is separated from the others by an underscore. The first three elements of the row are fixed and represent Time Replicate Control all separated by an underscore: Time_Replicate_Control. All the other elements refer to various experimental conditions. Considering two different conditions to be evaluated each row is made of 5 elements: Time_Replicate_Control_cond2 all separated by an underscore. Having an experiment made of 9 arrays, with two time points, 0h and 24h, in triplicate, and two different experimental conditions to be evaluated, the affylmGUI target file will look like:

Name	FileName	Target
mC1	M1.CEL	0_1_1_0_0
mC2	M4.CEL	0_1_1_0_0
mC3	M7.CEL	0_1_1_0_0
mE1	M3.CEL	24_2_0_1_0
mE2	M6.CEL	24_2_0_1_0
mE3	M9.CEL	24_2_0_1_0
mI1	M2.CEL	24_3_0_0_1
mI2	M5.CEL	24_3_0_0_1
mI3	M8.CEL	24_3_0_0_1

Usage

masigpro.edesign()

Author(s)

Raffaele A Calogero

See Also

masigpro, masigpro.view

masigpro

The function executes maSigPro analysis

Description

The function creates: 1. Create a regression matrix for the full regression model (make.design.matrix function). 2. Computes the p-value associated to the F-Statistic of the model, which is used to select significant genes (p.vector function). 3. Applies a variable selection procedure to find significant variables for each gene (T.fit function). This will ultimatelly be used to find which are the profile

26 metaArrayIC

differences between experimental groups. 4. Finally, it generates lists of significant genes according to R-squared of the models (get.siggenes function). To know more about the various steps see maSigPro help.

Usage

```
masigpro()
```

Author(s)

Raffaele A Calogero

See Also

masigpro.edesign, masigpro.view

masigpro.view

The function allows the visualization of maSigPro results

Description

The function is a graphical implementation of the maSigPro PlotGroups function. To know more about it see maSigPro help.

Usage

```
masigpro.view()
```

Author(s)

Raffaele A Calogero

See Also

masigpro.edesign, masigpro

metaArrayIC

Graphical interface to metaArray Integrative Correlation function

Description

The integrative correlation analysis (Parmigiani et al., 2004) is a convenient tool to monitor the interstudy concordance of within-study correlations of gene expression. The gene-specific reproducibility score takes the correlation between each gene and all other genes within individual study and calculate the average correlation of these correlations across all pairs of studies.

Usage

```
metaArrayIC()
```

metaArrayMerge 27

Author(s)

Raffaele A Calogero

References

MergeMaid package and metaArray Package

metaArrayMerge

Tool to create a merge object for metaArray package

Description

This function will create an ExpressionSet from a study starting from a tab delimite file and a target file this ExpressionSet will be merged with the NormalizedAffyData if they contain the same number of row and rownames in the same order. Data generated with this function could be analyzed using metaArrayIC function.

Usage

```
metaArrayIC()
```

Author(s)

Raffaele A Calogero

See Also

mataArrayIC

ML.edesign

The function creates an data frame containing the parameters useful for class prediction

Description

This function uses a specific configuration of Target column of the affylmGUI target file. To know more about target file see affylmGUI help. Each row of the column named Target, in the affylmGUI target file, describes the array on the basis of the experimental design. Each element needed for the construction of the data frame is separated from the others by an underscore. All the other elements refer to experimental conditions or clinical parameters. The absence of a parameter NEEDS to be described in the Target file by NA Considering two different conditions to be evaluated each row is made of 5 elements: Time_Replicate_Control_cond1_cond2 all separated by an underscore. Having an experiment made of 9 arrays with 4 different experimental parameters the affylmGUI target file will look like:

Name	FileName	Target
mC1	M1.CEL	$0_1_pos_0_NA$
mC2	M4.CEL	0_1_pos_0_yes
mC3	M7.CEL	0_1_neg_0_no

28 moex.variantexons

mE1	M3.CEL	24_2_neg_1_NA
mE2	M6.CEL	24_2_NA_1_yes
mE3	M9.CEL	24_2_neg_1_yes
mI1	M2.CEL	12_3_0_pos_yes
mI2	M5.CEL	12_3_0_pos_no
mI3	M8.CEL	12 3 0 pos no

Usage

```
ML.edesign()
```

Author(s)

Raffaele A Calogero

moex.annotation

Annotation data for CORE subset of mouse exon array 1.0 ST

Description

These annotation data are derived from Affymetrix annotation file moex10stv1na23mm8.

Usage

```
moex.annotation
```

Format

A data frame with 5 observations: PROBESETID, SYMBOL, DESCRIPTION, CYTOBAND, ACC

References

Affymetrix web site

moex.variantexons table linking exon-level probe set ids to variant exons

Description

The table is made integrating exon-level netaffx annotaiton with UCSC derived annotation. Specifically exon-level rpobesets are mapped on variant exons, i.e. those that are specific of only a subset of all isoforms associated to a gene. This table is used to filter differentially expressed exons to select only those associated to variant isoforms

Usage

```
moex.variantexons
```

mogene.annotation 29

Format

A data frame with 12 observations: affyname, affystart, affyend, affywidth, affystrand, vspname, vspstart, vspend, vspwidth, vspstrand, chr, genome

References

Affymetrix web site and UCSC database

mogene.annotation Annotation data for mouse gene array 1.0 ST

Description

These annotation data are derived from Affymetrix annotation file mogene 10stv1na24hg18.

Usage

```
mogene.annotation
```

Format

A data frame with 5 observations: PROBESETID, SYMBOL, DESCRIPTION, CYTOBAND, ACC

References

Affymetrix web site

myExpresso

Running the affy expresso function with the widget

Description

Various probe set intensity summary and normalization can be customized using the expresso function.

Usage

```
myExpresso()
```

Details

This function run expresso with the graphical interface for parameters selection. It is important to note that expresso is more slow than the C coded rma)

Author(s)

30 ocPlotPCA

normBoxplot

Box plot of the arrays data available in NormalizeAffyData slot

Description

Box plot visualization of normalized array data

Usage

```
normBoxplot()
```

Author(s)

Raffaele A Calogero

ocPlotHist

Gene/Exon level density plots

Description

This function runs a modified version of the plotHist of the affycoretools to be used to check density distribution plots for gene and exon expression data generated by expression console.

Usage

```
ocPlotHist()
```

Author(s)

Raffaele A Calogero

See Also

ocPlotPCA

ocPlotPCA

Gene/Exon level density plots

Description

This function runs a modified version of the plotPCA of the affycoretools to be used to check density distribution plots for gene and exon expression data.

Usage

```
ocPlotPCA()
```

oneChannelGUI 31

Author(s)

Raffaele A Calogero

See Also

ocPlotHist

oneChannelGUI-package

Set of functions extending the capability of affylmGUI package

Description

This package is directed to Bioconductor beginners that have little or no experience in writing R code. The package implements, as simple functions accessible over the affylmGUI graphical interface, some code useful for QC, data filtering, data output manipulation and identification of GO enriched classes.

Details

Package: oneChannelGUI

Type: Package Version: 1.0

Date: 2006-12-16

License: GPL version 2 or newer

Author(s)

Author: Raffaele A Calogero Maintainer: Raffaele A Calogero <raffaele.calogeor@unito.it>

Examples

```
# library(oneChannelGUI)
## To start the oneChannelGUI with the modifications
#oneChannelGUI()
```

oneChannelGUI

Starting oneChannelGUI package

Description

Starting oneChannelGUI package. oneChannelGUI contains a set of functions extending the capabilities of affylmGUI package

32 oneChannelGUI

Usage

```
AboutextendedaffylmGUI()
AboutaffylmGUI()
oneChannelGUI()
oneChannelGUIHelp()
maSigProHelp()
siggenesHelp()
oneChannelGUIHelp()
initialize.extAffylmGUI()
OpenExonandTargetsfiles()
GOstatsHelp()
NewLimmaFile()
OpenLimmaFile()
OpenALimmaFile(FileName)
OpenAFile(FileName)
OpenExonFile()
OpenLargeFile()
changeMenu()
oneChannelGUI.start()
libraryFilesDir()
whichKindOfArray()
intronicBq()
ExportNormalizedExpressionValues1()
ExportNormalizedExpressionValues()
ExportfeatureNames()
SaveAsLimmaFile()
addAnnLib()
OpenCDFandTargetsfiles()
ComputeLinearModelFit()
GetNormalizationMethod()
NormalizeNow()
ComputeContrasts()
midasFilter()
dabgFilter()
largedatasetNorm()
RankProdHelp()
affyPLMHelp()
genefilterHelp()
pamrHelp()
pdmclassHelp()
sizepowerHelp()
ssizeHelp()
OpenAGeoFile(FileName)
OpenGeoFile()
OpenGeoFiles()
recoverUnfiltered()
delete.ML()
affyPlotMA()
changes()
aptFolder()
deleteLocalData()
SetED()
```

OpenBeadStudioFiles 33

```
chooseEDir()
.annotation(eset)
.myrk(x,df)
```

Arguments

FileName Internal argument not to be set by the user eset Internal argument not to be set by the user x Internal argument not to be set by the user df Internal argument not to be set by the user

Details

This function launches a modify version of the Graphical User Interface by James Wettenhall for the limma package by Gordon Smyth. The GUI uses Tk widgets (via the R TclTk interface by Peter Dalgaard) in order to provide a simple interface to various tools for quality control and statistical analysis of Affymetrix gene chips.

Author(s)

Raffaele A Calogero

Examples

```
# library(oneChannelGUI)
## To start the affylmGUI with the modifications
#oneChannelGUI()
```

OpenBeadStudioFiles

Read BeadStudio expression data file

Description

Read BeadStudio expression data file

Usage

```
OpenBeadStudioFiles()
```

Details

Reads an Illumina intnesity data file produced by BeadStudio. Using BeadStudio version 'One" the file will have a 'gene profile.csv' extension and using version "Two" the file will have a .txt extension. See package vignette for more information. Multiple filenames can be specified as a vector and the data are then combined into one output file. This function should only really be used for custom analysis as the beadAnalysis() function provides easier, flexible use.

Author(s)

Derived from readBead by Gareth Elvidge \(\)gareth.elvidge@well.ox.ac.uk \(\)

34 plierToZero

OpenLargefiles

This function loads large data set made from tab delimited files

Description

The function creates and expressionSet starting from al file containing the expression data in a tab delimited format. This file is loaded together with the description of the clinical parameters present in Target This function uses a specific configuration of Target column of the affylmGUI target file. To know more about target file see affylmGUI help. Each row of the column named Target, in the affylmGUI target file, describes the clinical parameters. Each clinical parameter is separated from the others by an underscore. The affylmGUI target file will look like:

Name	FileName	Target
mC1	M1.CEL	pos_yes_1_NA_0
mC2	M4.CEL	pos_no_2_NA_0
mC3	M7.CEL	neg_no_3_pos_0
mE1	M3.CEL	neg_yes_3_neg_0
mE2	M6.CEL	neg_no_NA_1_0
mE3	M9.CEL	neg_yes_3_pos_0
mI1	M2.CEL	pos_no_2_neg_1
mI2	M5.CEL	pos_yes_2_pos_1
mI3	M8.CEL	pos_no_2_pos_1

Usage

OpenLargefiles()

Author(s)

Raffaele A Calogero

plierToZero

Setting to 0 low log2 intensity values generated with plier

Description

The calculation of log2 of probe set intensity by mean of plier generates a set of intensities very low this function will set to 0 all the log2 intensities below 1 produced by iter-plier or plier alghoritm

Usage

```
plierToZero()
```

Author(s)

plotGO 35

plotGO

Plotting parents of a GO term with few mouse click

Description

To know more on the parents of a specific GO term you can use the plotGO function

Usage

```
plotGO()
```

Details

A GO term to be investigated for its parents has to be placed in the graphical window.

Author(s)

Raffaele A Calogero

See Also

GOenrichment, extractAffyids

PlotOptionsv1

A modified version of the function used in affyPLM library

Description

As default the plots are generated on the R GUI to reduce RAM consumption.

Usage

```
PlotOptionsv1()
```

Author(s)

36 raex.annotation

plotVariantSI

This function plots on UCSC genome browser data derived by variantSI filter

Description

This function plots on UCSC genome browser data derived by variantSI filter on the basis of chromosome annotation

Usage

```
plotVariantSI()
```

Author(s)

Raffaele A Calogero

See Also

variantSI, variantExons, makeBED15

raex.annotation

Annotation data for CORE subset of rat exon array 1.0 ST

Description

These annotation data are derived from Affymetrix annotation file raex.annotation

Usage

```
raex.annotation
```

Format

A data frame with 5 observations: PROBESETID, SYMBOL, DESCRIPTION, CYTOBAND, ACC

References

Affymetrix web site

raex.variantexons 37

raex.variantexons table linking exon-level probe set ids to variant exons

Description

The table is made integrating exon-level netaffx annotaiton with UCSC derived annotation. Specifically exon-level rpobesets are mapped on variant exons, i.e. those that are specific of only a subset of all isoforms associated to a gene. This table is used to filter differentially expressed exons to select only those associated to variant isoforms

Usage

raex.variantexons

Format

A data frame with 12 observations: affyname, affystart, affyend, affywidth, affystrand, vspname, vspstart, vspend, vspwidth, vspstrand, chr, genome

References

Affymetrix web site and UCSC database

ragene.annotation Annotation data for rat gene array 1.0 ST

Description

These annotation data are derived from Affymetrix annotation file ragene 10stv1na24hg18.

Usage

ragene.annotation

Format

A data frame with 5 observations: PROBESETID, SYMBOL, DESCRIPTION, CYTOBAND, ACC

References

Affymetrix web site

38 rankProd

rankingConversion This function transforms intensity data in normalized ranks

Description

This function transforms intensity data in normalized ranks, i.e. high intensity genes will have a value near to 0 as instead low intensity genes a normalized rank near to 1.

Usage

```
rankingConversion()
```

Author(s)

Raffaele A Calogero

rankProd

graphical interface to rank product method implemented in RankProd Bioconductor library.

Description

To know more about rank product method see RankProd help.

Usage

```
rankProd()
```

Details

The target file for the RankProd implementation contain the origin of the data as a number separated by an under score from the corresponding covariate. If all data are from the same origin the origin definition is not needed. Therefore target will contain only the covariates.

Name	FileName	Target
mC1	M1.CEL	CTRL_1
mC2	M4.CEL	CTRL_1
mC3	M7.CEL	CTRL_2
mE1	M3.CEL	CTRL_2
mE2	M6.CEL	CTRL_2
mE3	M9.CEL	TRT_1
mI1	M2.CEL	TRT_1
mI2	M5.CEL	TRT_2
mI3	M8.CEL	TRT_2

Author(s)

rawBoxplotPN 39

rawBoxplotPN

Plotting raw log2 intensities from controls

Description

This function produces a box plot of the log2 raw intensities, extracted directly from CEL files, for positive and negative controls presente in XXXX.control.ps Affymetrix library file. Positive and negative controls are made of housekeeping exon and introns. It gives an idea of signal behaviour before data normalization both in the high and low intensity range

Usage

```
rawBoxplotPN()
```

Author(s)

Raffaele A Calogero

rawpCheck

Raw p-value distribution from limma analysis by a mouse click

Description

This function allow to visualize the histogram of raw p-value distribution generated by limma analysis.

Usage

```
rawpCheck()
```

Details

The histogram of raw p-value distribution will show if raw p-values are uniform in the non significant range and therefore the BH correction can be applied.

Note

BH is the most used method for the correction of type I errors in microarray analysis. However, it has some limitation due to the initial hypotheses: The gene expressions are independent from each other. The raw distribution of p values should be uniform in the non significant range.

Author(s)

Raffaele A Calogero

References

To know more see limma package help

40 reviqrFilter

refseqDownload

Retrieving Reference Sequences from NCBI ftp

Description

This function retieves refeence sequences from NCBI ftp. RefSeq are used for mapping exon-level probe sets to refseq specific isoforms.

Usage

```
refseqDownload()
```

Author(s)

Raffaele A Calogero

retrievePSRseq

This function, given a file containing exon-level probesets retrieves Probe selection Regions, PRS, from RRE db

Description

This function retrieve from RRE the PSR sequences associated to the exon-level probe sets

Usage

```
retrievePSRseq()
```

Author(s)

Raffaele A Calogero

reviqrFilter

Reverse interquantile filtering with a mouse click

Description

This function implements a reverse version of the interquantile filtering proposed by Heydebreck in 2004 to select low variance genelevel probe set. To be used to remove putative differentially expressed genes that will make more difficult the detection of alternative splicing events.

Usage

```
reviqrFilter()
```

runningJetta 41

Details

This function can be used in a analysis focused to the detection of alternative splicing events. The aim of this non specific filtering is to remove the genes that are likely to carry information about the phenotypes under investigation at gene level. This filtering remove genes that show strong changes within the experimental points at the gene level.

Author(s)

Raffaele A Calogero

References

Heydebreck et al. Bioconducotor project Papers 2004

See Also

dabgFilter, crosshybFilter

runningJetta

graphical interface to MADS/jetta R library.

Description

MADS, which stands for Microarray Analysis of Differential Splicing, is a tool to identify differential alternative splicing by exon array. The principle of MADS is to increase the precision of exon-level and gene-level expression estimates by correcting, as much as possible, noise in observed probe intensities due to background and cross-hybridization. MADS incorporates a series of novel algorithms motivated by the probe-rich design of exon-tiling arrays, such as background correction, iterative probe selection and removal of sequence-specific cross-hybridization to off-target transcripts. MADS was published in RNA,2008,14(8): 1470-1479. Junction and Exon array Toolkit for Transcriptome Analysis (JETTA) is compacted version of MADS.

Usage

runningJetta()

Details

Expression indexes are calculated as the order of Background Correction, Normalization and Summarization. In the Summarization step, background corrected and normalized probe intensities of a meta probeset are summarized to expression of the meta probeset. Meta probesets can be defined as gene/transcript clust/exon level.

Background Correction JETTA estimates background signal using background probes and subtracts it from the probe intensity. If the probe intensity is less than the estimated background signal, the background subtracted signal is truncated to 1. Estimation of background signal is based on several models: Median GC: median of background probe signal of the same GC counts MAT: linear model of probe sequence with 80 parameters. see Kapur et al, 2007

Normalization Normalization of JETTA is done for core probes defined in probeset annotation file. If the PSA file is not specified, it considers all probes in the MPS files as core probes. Median scaling: scaling each array so that its median is 100 Quantile: all probes of the same signal quantile have the same signal

42 sample.size.evaluation

Summarization LiWong model: multiplication model of expression and probe effect, see Li and Wong, 2001 Probe selection: select probes based on cross-array correlation of signal. see Xing et al, 2006 Median-polish

Alternative Splicing Detection Detecting alternatively expressed PSR/Exon between two sample groups based on background corrected and normalized probe intensities. It has several criteria to filter out transcript clusts and probes from the analysis. TC expression level: excluding low-expressed transcript clusts TC expression fold change: excluding transcript clusts which have big fold change between two groups Extreme probe signal: excluding probes of which signal is extremely high Cross-hybridized probes: excluding cross hybridized probes, currently pre-calculated results are needed

Author(s)

jseok@stanford.edu

```
sample.size.evaluation1
```

The function executes functions from the sizepower packahe

Description

This function represents a tool for helping users to understand the trade off between sample size and statistical power. To know more about see ssize help.

Usage

```
sample.size.evaluation1()
```

Details

see sizepower help

Author(s)

Raffaele A Calogero

```
sample.size.evaluation
```

The function executes and plots results from ssize and delta fulction from the ssize packahe

Description

This function represents avisual tool for helping users to understand the trade off between sample size and statistical power. To know more about see ssize help.

Usage

```
sample.size.evaluation()
```

showDataset 43

Details

Both ssize and delta outputs are calculated using the BH type I error correction instead of the Bonferroni used as default in the ssize package. Furthermore, instead using the control group variance, this implementation uses the common variance described in Wei et al. BMC Genomics. 2004, 5:87

Main assumptions: A microarray experiment is set up to compare gene expressions between one treatment group and one control group. Microarray data has been normalized and transformed so that the data for each gene is sufficiently close to a normal distribution that a standard 2-sample pooled-variance t-test will reliably detect differentially expressed genes.

Author(s)

Raffaele A Calogero

showDataset

Grabbing info about the available expression set

Description

The size of the normalized expression set can change upon filtering. This function show info about the exact size of the data set.

Usage

```
showDataset()
```

Author(s)

Raffaele A Calogero

showTopTable

 $\label{lem:modification} \textit{Modification of the function immplemented in affylmGUI to generate} \\ \textit{a topTable}$

Description

Modification of the function immplemented in affylmGUI to generate a topTable. To know more about topTable see limma package help

Usage

```
showTopTable(...,export=FALSE)
```

Arguments

export defining the possibility to export data
... Arguments to be passed to methods

Author(s)

44 simFilter

siggenes	The function executes SAM analysis implemented in siggenes biocon-
	ductor library

Description

To know more about SAM in Bioconductor see siggenes help.

Usage

```
siggenes()
```

Author(s)

Raffaele A Calogero

simFilter This function allows filtering on the basis of the average splice index mean or min difference between two groups

Description

Filtering out gene/exon level probe sets associated to average splice index mean or min difference between two groups lower than user defined value

Usage

```
simFilter()
```

Note

This fuction needs the presence of Splice Index.

Author(s)

Raffaele A Calogero

See Also

simFilter

spliceIndex 45

spliceIndex

This function coverts the exon intenisty data in a slice index

Description

Exons intensities are divided for the expression of the corresponding gene, as descrived by Clark et al. Science 2002 May 3;296(5569):907-10.

Usage

```
spliceIndex()
```

Details

The function is not yet optimized, therefore it could take quite a long time to compute spliceIndex if more than 1000 genes are used.

Author(s)

Raffaele A Calogero

See Also

inspecting.splice.index

targetWidget

Widget to create a target file to load .CEL files

Description

Widget to create a target file to load .CEL files to be used with NewLimmaFile function.

Usage

```
targetWidget()
```

Author(s)

46 trainTest

templa $Generating \ a \ template \ A \ to \ be \ uploaded \ in \ Ingenuity \ Pathways \ analysis \ (IPA)$

Description

A template A file to be used in Ingenuity can be generated starting from a topTable containing the full array data.

Usage

templA()

Note

Template A file will contain a column with the gene ID, a column with fold change, a column with true p-value and a column with p-values for discriminating between the set of differentially expressed probe sets and the background. This coulumn is needed to allow IPA to identify the set of enriched functional classes associated to the differentially expressed probe sets.

Author(s)

Raffaele A Calogero

See Also

IPAlistFilter

trainTest

Creating a training set and a test set by a mouse click

Description

This function allows the creation of a training set and a test set to be used for classification purposes.

Usage

```
trainTest()
```

Details

User will be asked to assign names to the available classification parameters. User will be asked to select the number associated to one of the available classification parameters. The training set will be made, using the selected classification parameter and it will be made of 2/3 of the original data set. The test set will be the remaining 1/3.

Author(s)

updateLibs 47

updateLibs	This function allows to update the present installation of Bioconductor
------------	---

Description

The function allows the updating of local installation of Bioconductor. It might be quite long depending on the internet connection speed.

Usage

```
updateLibs()
```

Author(s)

Raffaele A Calogero

variantExons This function is used to generate a table containing exon-level probe set data linked to variant exons

Description

Internal oneChannelGUI annotation data linking exon-level probesets to variant exons, i.e. those exons that ar specific for a subgroup of the isoforms associated to a specific gene, can be upgraded using this function. Annotation RDA files need to be saved in the data subdir of oneChannelGUI dir. Windows users need to exchange the older copies present in Rdata.zip, simply dragging them in the zip file.

Usage

```
variantExons()
```

Author(s)

Raffaele A Calogero

variantSI	This function allows filtering on the basis of variant exons	
-----------	--	--

Description

The function interesects a list of alternative splice exon-level probe sets detected by oneChannelGUI analysis and intersects it to the list of exon-level probe sets associated to variant exons, i.e. these exons that are associated only to a subset of all isoforms associated to a gene. The table of variant Exon is stored in RRE database and is retrieved by the function updating the UCSC tables linking probe set ids with variant exons, located in the General Tools menu. The variantSI function also attaches to the variant exons table the deltaSI associated exon-level probe set.

VennDiagram

Usage

```
variantSI()
```

Note

This fuction needs the presence of Splice Index

Author(s)

Raffaele A Calogero

See Also

variantExons, plotVariantSI

VennDiagram

Venn diagrams using two or three lists

Description

Venn diagrams can be generated using probe sets ids or Entrez gene ids saved in flat files.

Usage

```
VennDiagram()
```

Details

The function asks to the user to select two/three files containing probe set ids or EGs separated by carriage return. Each file should contain only one column and no header.

Author(s)

Index

*Topic IO	plierToZero, 33
oneChannelGUI, 30	plotVariantSI,35
*Topic datasets	rankingConversion, 37
cosieHscore,4	rawBoxplotPN,38
cosieHsfull,5	refseqDownload,39
cosieMmcore, 5	retrievePSRseq, 39
cosieMmfull, 6	simFilter, 43
crosshybhuex.annotation, 8	variantSI,46
${ t crosshyb}{ t moex.annotation,8}$	*Topic methods
${ t crosshybraex.annotation,9}$	buildingLocalAnnotation, 2
geoVSbioc, 15	combineGeoMSF, 3
huex.annotation, 16	createGeoTarget,7
huex.variantexons, 17	crosshybFilter,7
hugene.annotation, 17	dfMAplot, 10
moex.annotation, 27	extractAffyids, 14
moex.variantexons, 27	geneExonLibs, 15
mogene.annotation, 28	GOenrichment, 16
raex.annotation,35	inspecting.one.splice.index,
raex.variantexons, 36	18
ragene.annotation, 36	inspecting.splice.index, 18
*Topic hplot	intensityFilter, 19
dfMAplot, 10	IPAlistFilter, 19
normBoxplot, 29	iqrFilter, 20
plotG0, 34	listFilter, 21
rawpCheck, 38	OpenBeadStudioFiles, 32
VennDiagram,47	PlotOptionsv1,34
*Topic manip	reviqrFilter, 39
AptMidas, 1	showDataset,42
biomartFilter,2	showTopTable, 42
colExtract, 3	spliceIndex,44
consistentFilters,4	templA,45
cosieWrapper,6	updateLibs,46
EG2probeset, 11	variantExons,46
erankProdAltSpl, 12	*Topic package
erankProdAltSplFilter, 11	oneChannelGUI-package, 30
exonsSpecific2as, 13	*Topic utilities
filteringTable, 14	exonContrasts, 12
makeBED15,22	exonTopTableExtract, 13
mapping2ensembl, 22	limmaExons, 21
mapping2exon, 23	masigpro, 24
mapping2RefSeq, 23	masigpro.edesign,24
metaArrayIC, 25	masigpro.view,25
metaArrayMerge, 26	ML.edesign, 26

50 INDEX

myExpresso, 28	erankProdAltSpl, 12
ocPlotHist, 29	erankProdAltSplFilter,11
ocPlotPCA, 29	exonContrasts, 12
OpenLargefiles, 33	exonsSpecific2as, 13
rankProd, 37	exonTopTableExtract, 13
runningJetta,40	ExportfeatureNames
sample.size.evaluation,41	(oneChannelGUI), 30
sample.size.evaluation1,41	ExportNormalizedExpressionValues
siggenes, 43	(oneChannelGUI), 30
targetWidget,44	ExportNormalizedExpressionValues1
trainTest, 45	(oneChannelGUI), 30
.annotation (oneChannelGUI), 30	extractAffyids, 14
.myrk(oneChannelGUI), 30	
	filteringTable, 14
AboutaffylmGUI (oneChannelGUI), 30	
AboutextendedaffylmGUI	geneExonLibs, 15
(oneChannelGUI), 30	genefilterHelp(oneChannelGUI),30
addAnnLib (oneChannelGUI), 30	geoVSbioc, 15
affyPLMHelp(oneChannelGUI), 30	GetNormalizationMethod
affyPlotMA(oneChannelGUI), 30	(oneChannelGUI), 30
aptFolder (oneChannelGUI), 30	GOenrichment, 16
AptMidas, 1	GOstatsHelp(oneChannelGUI), 30
biomartFilter, 2	huex.annotation, 16
buildingLocalAnnotation, 2	huex.variantexons, 17
	hugene.annotation, 17
changeMenu (oneChannelGUI), 30	
changes (oneChannelGUI), 30	initialize.extAffylmGUI
chooseEDir(oneChannelGUI), 30	(oneChannelGUI), 30
colExtract, 3	inspecting.one.splice.index, 18
combineGeoMSF, 3	inspecting.splice.index, 18
ComputeContrasts (oneChannelGUI),	intensityFilter, 19
30	intronicBg(oneChannelGUI), 30
ComputeLinearModelFit	IPAlistFilter, 19
(oneChannelGUI), 30	iqrFilter, 20
consistentFilters,4	
cosieHscore, 4	largedatasetNorm(oneChannelGUI),
cosieHsfull,5	30
cosieMmcore, 5	libraryFilesDir(oneChannelGUI),
cosieMmfull, 6	30
cosieWrapper, 6	limmaExons, 21
createGeoTarget, 7	listFilter, 21
crosshybFilter, 7	
crosshybhuex.annotation, 8	makeBED15, 22
crosshybmoex.annotation, 8	mapping2ensembl, 22
crosshybraex.annotation,9	mapping2exon, 23
	mapping2RefSeq, 23
dabgFilter(oneChannelGUI), 30	masigpro, 24
delete.ML (oneChannelGUI), 30	masigpro.edesign, 24
deleteLocalData (oneChannelGUI),	masigpro.view, 25
30	maSigProHelp (oneChannelGUI), 30
dfMAplot, 10	metaArrayIC, 25
700 1 11	metaArrayMerge, 26
EG2probeset, 11	midasFilter(oneChannelGUI),30

INDEX 51

ML.edesign, 26 moex.annotation, 27 moex.variantexons, 27 mogene.annotation, 28 myExpresso, 28	refseqDownload, 39 retrievePSRseq, 39 reviqrFilter, 39 runningJetta, 40
NewLimmaFile (oneChannelGUI), 30 NormalizeNow (oneChannelGUI), 30 normBoxplot, 29	sample.size.evaluation, 41 sample.size.evaluation1, 41 SaveAsLimmaFile(oneChannelGUI), 30
ocPlotHist, 29 ocPlotPCA, 29 oneChannelGUI, 30 oneChannelGUI—package), 30 oneChannelGUI—package, 30 oneChannelGUI—package, 30 oneChannelGUI.start	SetED (oneChannelGUI), 30 showDataset, 42 showTopTable, 42 siggenes, 43 siggenesHelp (oneChannelGUI), 30 simFilter, 43 sizepowerHelp (oneChannelGUI), 30 spliceIndex, 44 ssizeHelp (oneChannelGUI), 30 targetWidget, 44 templA, 45 trainTest, 45 updateLibs, 46 variantExons, 46 variantSI, 46 VennDiagram, 47 whichKindOfArray (oneChannelGUI), 30
OpenGeoFile (oneChannelGUI), 30 OpenGeoFiles (oneChannelGUI), 30 OpenLargeFile (oneChannelGUI), 30 OpenLargefiles, 33 OpenLimmaFile (oneChannelGUI), 30	
pamrHelp (oneChannelGUI), 30 pdmclassHelp (oneChannelGUI), 30 plierToZero, 33 plotGO, 34 PlotOptionsv1, 34 plotVariantSI, 35	
raex.annotation, 35 raex.variantexons, 36 ragene.annotation, 36 rankingConversion, 37 rankProd, 37 RankProdHelp(oneChannelGUI), 30 rawBoxplotPN, 38 rawpCheck, 38 recoverUnfiltered (oneChannelGUI), 30	