goTools

April 19, 2010

EndNodeList

Reference GO nodes list.

Description

The function <code>EndNodeList</code> builds the default end node list used in <code>ontoCompare</code>. <code>CustomEndNodeList</code> helps you build a list of children nodes starting from a GO id of interest.

Usage

```
EndNodeList()
CustomEndNodeList(id,rank=1)
```

Arguments

 $\label{eq:Valid} \text{Valid GO id: "GO:XXXXXXXX"}.$

rank Number of levels of children of the GO DAG under id you want to add to your

nodes list.

Value

<code>EndNodeList</code> returns a vector of all GO ids 1 level below MF("GO:0003674"), BP("GO:0008150") and CC("GO:0005575"). MF, BP and CC nodes are included. <code>CustomEndNodeList</code> returns a vector of all GO ids children of id, <code>rank</code> levels below it.

Author(s)

Yee Hwa (Jean) Yang, Agnes Paquet

See Also

ontoCompare

2 goTools

Examples

```
## Examples use the probeID dataset. For description type ? probeID.
## library(GO.db)
## EndNodeList()
## MFendnode <- CustomEndNodeList("GO:0003674", rank=2)

## Example (not run)
## data(probeID)
## ontoCompare(affylist, probeType="hgu133a", endnode=MFendnode, goType="MF")</pre>
```

Internal functions Internal goTools functions

Description

Internal goTools functions

Details

These are not to be called by the user.

goTools

Wrapper functions

Description

This functions will allow you to describe and compare sets of oligo ids using Gene Ontology database

Usage

```
ontoCompare(genelist,probeType=c("GO","hgu133a"),
goType="All", endnode, method=c("TGenes", "TIDS", "none"), plot=FALSE,
...)
ontoPlot(objM, names.arg=NULL,beside=TRUE, las=2,legend.text=TRUE, ...)
```

Arguments

genelist list of list of valid probe ids.

 ${\tt method}$

method used to evaluate the percentage of oligos for each end-node. 'TGenes' = for each end node, return the number of direct children found / total number of probe ids. (default). This includes oligos which do not have GO annotations. 'TIDS' = for each end node, return the number of direct children found / total number of GO ids describing the list. 'none' = for each end node, return the number of direct children found.

probeID 3

probeType	type of input given to the function. Valid probe types include GO ids and any probes ids for which a BioC annotation package providing a mapping to GO is available. onto Compare is expecting valid probe ids.
goType	help sort the data by type. If 'All' (default), all oligos are taken into account. 'BP' restricts information to Biological Process, 'CC' to Cellular Component, and 'MF' to Molecular Function.
plot	logical: if 'TRUE', results are output as a graph.
endnode	list of GO ids corresponding to end-nodes of interest.
beside	Logical. If 'TRUE', the bars of the barplot are portrayed as juxtaposed bars. See ?barplot for more details.
las	numeric: if las=2, the axis labels are displayed perpendicular to the axis. See ?par for more details.
legend.text	vector of text used to construct a legend for the plot. See ?barplot for more details.
objM	results from ontoCompare.
names.arg	Labels to use in ontoPlot.
	extra layout parameters to be passed to ontoPlot.

Value

Returns the percentage of probes children of nodes contained in endnode. If 'plot' = TRUE, results are plotted as a pie chart or a bargraph.

Author(s)

Yee Hwa (Jean) Yang, Agnes Paquet

Examples

```
# Examples use the probeID dataset. For description type ?probeID.
# Not run

#library(GO.db)
#data(probeID)
#ontoCompare(affylist, probeType="hgu133a", plot=TRUE)
#res <- ontoCompare(operonlist["L1"], probeType="operon", method="TIDS")
#ontoPlot(res, cex=0.7)</pre>
```

probeID List of probe ids from Affymetrix hgu133a chip and Operon Version 2
Human oligos

Description

The probeID dataset consists of two lists of randomly chosen probe Ids. affylist contains 3 sets of ids from Affymetrix hgu133a. operonlist contains 2 sets of ids from Operon Version 2.

Usage

```
data(probeID)
```

Index

```
*Topic datasets
   probeID, 3
*Topic file
   EndNodeList, 1
   goTools, 2
*Topic hplot
   goTools, 2
*Topic manip
   EndNodeList, 1
*Topic methods
   Internal functions, 2
affylist(probeID), 3
CustomEndNodeList (EndNodeList), 1
EndNodeList, 1
getGOID(Internal functions), 2
getOntology(Internal functions),
goChildren(Internal functions), 2
goParents (Internal functions), 2
goTools, 2
Internal functions, 2
isEndNode(Internal functions), 2
ontoCompare, I
ontoCompare (goTools), 2
ontoCompare.main(Internal
       functions), 2
ontoPlot (goTools), 2
operonlist (probeID), 3
parentsVectWraper(Internal
       functions), 2
probeID, 3
```