Package 'FRASER'

July 27, 2025

Type Package

Title Find RAre Splicing Events in RNA-Seq Data

Version 2.4.3

Date 2025-06-04

Description Detection of rare aberrant splicing events in transcriptome profiles. Read count ratio expectations are modeled by an autoencoder to control for confounding factors in the data. Given these expectations, the ratios are assumed to follow a beta-binomial distribution with a junction specific dispersion. Outlier events are then identified as read-count ratios that deviate significantly from this distribution.

FRASER is able to detect alternative splicing, but also intron retention.

The package aims to support diagnostics in the field of rare diseases where RNA-seq is performed to identify aberrant splicing defects.

biocViews RNASeq, AlternativeSplicing, Sequencing, Software, Genetics, Coverage

License file LICENSE

URL https://github.com/gagneurlab/FRASER

BugReports https://github.com/gagneurlab/FRASER/issues

RoxygenNote 7.3.2 Encoding UTF-8 VignetteBuilder knitr

Depends BiocParallel, data.table, Rsamtools, SummarizedExperiment

Imports AnnotationDbi, BBmisc, Biobase, BiocGenerics, biomaRt, BSgenome, cowplot, DelayedArray (>= 0.5.11), DelayedMatrixStats, extraDistr, generics, GenomeInfoDb, GenomicAlignments, GenomicFeatures, GenomicRanges, IRanges, grDevices, ggplot2, ggrepel, HDF5Array, matrixStats, methods, OUTRIDER, pcaMethods, pheatmap, plotly, PRROC, RColorBrewer, rhdf5, Rsubread, R.utils, S4Vectors, stats, tibble, tools, utils, VGAM, RMTstat, pracma

2 Contents

Suggests magick, BiocStyle, knitr, rmarkdown, testthat, covr, TxDb.Hsapiens.UCSC.hg19.knownGene, org.Hs.eg.db, rtracklayer, SGSeq, ggbio, biovizBase, BSgenome.Hsapiens.UCSC.hg38, BSgenome.Hsapiens.NCBI.GRCh38, BSgenome.Hsapiens.1000genomes.hs37d5, BSgenome.Hsapiens.UCSC.hg19
LinkingTo RcppArmadillo, Rcpp
Collate variables.R getNSetterFuns.R FRASER-package.R FraserDataSet-class.R AllGenerics-definitions.R AllGenerics.R Fraser-pipeline.R annotationOfRanges.R beta-binomial-testing.R calculatePSIValue.R countRNAseqData.R example_functions.R filterExpression.R find_encoding_dimensions.R getURLs.R helper-functions.R mergeExternalData.R saveHDF5Objects.R RcppExports.R autoencoder.R updateD.R updateE.R updateRho.R pvalsNzscore.R makeSimulatedDataset.R fitCorrectionMethods.R plotMethods.R resultAnnotations.R zzz.R
git_url https://git.bioconductor.org/packages/FRASER
git_branch RELEASE_3_21
git_last_commit 9298168
git_last_commit_date 2025-06-04
Repository Bioconductor 3.21
Date/Publication 2025-07-27
Author Christian Mertes [aut, cre] (ORCID:
Maintainer Christian Mertes <mertes@in.tum.de></mertes@in.tum.de>
Contents
annotateRanges assayNames,FraserDataSet-method assays,FraserDataSet-method calculatePSIValues countRNA createTestFraserSettings estimateBestQ filterVariability fit FRASER FraserDataSet FraserDataSet 2 getter_setter_functions 2

annotateRanges 3

	injectOutliers	26
	K	
	length,FraserDataSet-method	28
	loadFraserDataSet	29
	makeSimulatedFraserDataSet	30
	mergeExternalData	31
	plotManhattan	32
	potentialImpactAnnotations	42
	psiTypes	44
	results,FraserDataSet-method	45
	samples	47
	subset.FRASER	51
Index		52
		_

annotate Ranges

Annotates the given FraserDataSet with the HGNC symbol with biomaRt

Description

Annotates the given FraserDataSet with the HGNC symbol with biomaRt

```
annotateRanges(
  fds,
  feature = "hgnc_symbol",
  featureName = feature,
 biotype = list("protein_coding"),
  ensembl = NULL,
 GRCh = 37
)
annotateRangesWithTxDb(
  fds,
  feature = "SYMBOL",
  featureName = "hgnc_symbol",
  keytype = "ENTREZID",
  txdb = NULL,
 orgDb = NULL,
  filter = list()
)
```

4 annotateRanges

Arguments

fds FraserDataSet

feature Defines which feature (default is HGNC symbol) should be annotated. Has to

be the biomaRt feature name or a column name in orgDb.

featureName The column name of the feature in the FraserDataSet mcols.

biotype The biotype for biomaRt.

ensembl The ensembl that should be used. If NULL, the default one is used (hsapi-

ens_gene_ensembl, GRCh37).

GRCh version to connect to. If this is NULL, then the current GRCh38 is used.

Otherwise, this can only be 37 (default) at the moment (see useEnsemb1).

keytype The keytype or column name of gene IDs in the TxDb object (see keytypes for

a list of available ID types).

txdb A TxDb object. If this is NULL, then the default one is used, currently this is

TxDb.Hsapiens.UCSC.hg19.knownGene.

orgDb An orgDb object or a data table to map the feature names. If this is NULL, then

org. Hs. eg. db is used as the default.

filter A named list specifying the filters which should be applied to subset to e.g. only

protein-coding genes for annotation. names (filter) needs to be column names

in the given orgDb object (default: no filtering).

Value

FraserDataSet

```
fds <- createTestFraserDataSet()</pre>
### Two ways to annotage ranges with gene names:
# either using biomart with GRCh38
 fds <- annotateRanges(fds, GRCh=38)</pre>
 rowRanges(fds, type="j")[,c("hgnc_symbol")]
})
# either using biomart with GRCh37
 fds <- annotateRanges(fds, featureName="hgnc_symbol_37", GRCh=37)</pre>
 rowRanges(fds, type="j")[,c("hgnc_symbol_37")]
})
# or with a provided TxDb object
require(TxDb.Hsapiens.UCSC.hg19.knownGene)
txdb <- TxDb.Hsapiens.UCSC.hg19.knownGene
require(org.Hs.eg.db)
orgDb <- org.Hs.eg.db</pre>
fds <- annotateRangesWithTxDb(fds, txdb=txdb, orgDb=orgDb)</pre>
```

```
rowRanges(fds, type="j")[,"hgnc_symbol"]
```

```
assayNames,FraserDataSet-method
```

Returns the assayNames of FRASER

Description

Returns the assayNames of FRASER

Usage

```
## S4 method for signature 'FraserDataSet'
assayNames(x)
```

Arguments

^

FraserDataSet

Value

Character vector

```
assays, FraserDataSet-method
```

Returns the assay for the given name/index of the FraserDataSet

Description

Returns the assay for the given name/index of the FraserDataSet

```
## S4 method for signature 'FraserDataSet'
assays(x, withDimnames = TRUE, ...)

## S4 replacement method for signature 'FraserDataSet,SimpleList'
assays(x, withDimnames = TRUE, HDF5 = TRUE, type = NULL, ...) <- value

## S4 replacement method for signature 'FraserDataSet,list'
assays(x, withDimnames = TRUE, HDF5 = TRUE, type = NULL, ...) <- value

## S4 replacement method for signature 'FraserDataSet,DelayedMatrix'
assays(x, withDimnames = TRUE, HDF5 = TRUE, type = NULL, ...) <- value</pre>
```

6 calculatePSIValues

Arguments

x FraserDataSet

... Parameters passed on to SummarizedExperiment::assays()

HDF5 Logical value indicating whether the assay should be stored as a HDF5 file.

type The psi type.

value The new value to which the assay should be set.

Value

(Delayed) matrix.

calculatePSIValues

PSI value calculation

Description

This function calculates the PSI values for each junction and splice site based on the FraserDataSet object

Usage

```
calculatePSIValues(
  fds,
  types = psiTypes,
  overwriteCts = FALSE,
  BPPARAM = bpparam()
)
```

Arguments

fds A FraserDataSet object

types A vector with the psi types which should be calculated. Default is all of jaccard,

psi5, psi3 and theta.

overwriteCts FALSE or TRUE (the default) the total counts (aka N) will be recalculated based

on the existing junction counts (aka K)

BPPARAM the BiocParallel parameters for the parallelization

Value

FraserDataSet

Examples

```
fds <- createTestFraserDataSet()
fds <- calculatePSIValues(fds, types="jaccard")
### usually one would run this function for all psi types by using:
# fds <- calculatePSIValues(fds)</pre>
```

countRNA

Count RNA-seq data

Description

The FRASER package provides multiple functions to extract and count both split and non-spliced reads from bam files. See Detail and Functions for more information.

```
countRNAData(
  fds,
 NcpuPerSample = 1,
 minAnchor = 5,
  recount = FALSE,
  BPPARAM = bpparam(),
  genome = NULL,
  junctionMap = NULL,
  filter = TRUE,
 minExpressionInOneSample = 20,
  keepNonStandardChromosomes = TRUE,
 countDir = file.path(workingDir(fds), "savedObjects", nameNoSpace(name(fds))),
)
getSplitReadCountsForAllSamples(
  fds,
 NcpuPerSample = 1,
  junctionMap = NULL,
  recount = FALSE,
  BPPARAM = bpparam(),
  genome = NULL,
  countFiles = NULL,
  keepNonStandardChromosomes = TRUE,
  outDir = file.path(workingDir(fds), "savedObjects", nameNoSpace(name(fds)),
    "splitCounts")
)
getNonSplitReadCountsForAllSamples(
  fds,
```

```
splitCountRanges,
  NcpuPerSample = 1,
  minAnchor = 5,
  recount = FALSE,
  BPPARAM = bpparam(),
  longRead = FALSE,
  outDir = file.path(workingDir(fds), "savedObjects", nameNoSpace(name(fds)),
    "nonSplitCounts")
)
addCountsToFraserDataSet(fds, splitCounts, nonSplitCounts)
countSplitReads(
  sampleID,
  fds,
  NcpuPerSample = 1,
  genome = NULL,
  recount = FALSE,
  keepNonStandardChromosomes = TRUE,
  bamfile = bamFile(fds[, sampleID]),
  pairedend = pairedEnd(fds[, sampleID]),
  strandmode = strandSpecific(fds[, sampleID]),
  cacheFile = getSplitCountCacheFile(sampleID, fds),
  scanbamparam = scanBamParam(fds),
  coldata = colData(fds)
)
mergeCounts(
  countList,
  fds,
  junctionMap = NULL,
  assumeEqual = FALSE,
  spliceSiteCoords = NULL,
  BPPARAM = SerialParam()
)
countNonSplicedReads(
  sampleID,
  splitCountRanges,
  fds,
  NcpuPerSample = 1,
 minAnchor = 5,
  recount = FALSE,
  spliceSiteCoords = NULL,
  longRead = FALSE
)
```

Arguments

fds A FraserDataSet object

NcpuPerSample A BiocParallel param object or a positive integer to configure the parallel back-

end of the internal loop per sample

minAnchor Minimum overlap around the Donor/Acceptor for non spliced reads. Default to

5

recount if TRUE the cache is ignored and the bam file is recounted.

BPPARAM the BiocParallel parameters for the parallelization

genome NULL (default) or a character vector specifying the names of the reference

genomes that were used to align the reads for each sample. The names have to be in a way accepted by the getBSgenome function. Available genomes can be listed using the available.genomes function from the BSgenome package. If genome is of length 1, the same reference genome will be used for all samples. If genome is supplied and strandSpecific(fds) == 0L (unstranded), then the strand information will be estimated by checking the dinucleotides found at the intron boundaries (see summarizeJunctions in GenomicAlignments package for details). This can e.g. help to avoid ambiguities when adding gene names

from a gene annotation to the introns in a later step.

junctionMap A object or file containing a map of all junctions of interest across all samples

filter If TRUE, splice sites of introns with low read support in all samples are not

considered when calculating the non-split reads. This helps to speed up the

subsequent steps.

minExpressionInOneSample

The minimal split read count in at least one sample that is required for an intron

to pass the filter.

keepNonStandardChromosomes

Logical value indicating if non standard chromosomes should also be counted.

Defaults to TRUE.

countDir The directory in which the tsv containing the position and counts of the junctions

should be placed.

... Further parameters passed on to Rsubread::featureCounts.

files. Should be a vector of paths to files containing the split read counts for the individual samples. Reading from files is only supported for tsv(.gz) or RDS files containing GRranges objects. The order of the individual sample files

should correspond to the order of the samples in the fds.

outDir The full path to the output folder containing the merged counts. If the given

folder already exists and stores a SummarizedExperiment object, the counts from this folder will be read in and used in the following (i.e. the reads are not recounted), unless the option recount=TRUE is used. If this folder doesn't exist or if recount=TRUE, then it will be created after counting has finished.

splitCountRanges

The merged GRanges object containing the positions of all the introns in the dataset over all samples.

longRead If TRUE, then the isLongRead option of Rsubread::featureCounts is used when counting the non spliced reads overlapping splice sites. The SummarizedExperiment object containing the position and counts of all the splitCounts introns in the dataset for all samples. nonSplitCounts The SummarizedExperiment object containing the position and non split read counts of all splice sites present in the dataset for all samples. sampleID The ID of the sample to be counted. bamfile The BAM file to be used to extract the counts. Defaults to the BAM file defined in the FraserDataSet object. pairedend TRUE or FALSE if the BAM file is paired end. Defaults to the value specified in the FraserDataSet object. strandmode 0 (no, default), 1 (stranded), or 2 (revers) to specify the used protocol for the RNA-seq experiment. cacheFile File path to the cache, where counts are stored. scanbamparam The ScanBamParam object which is used for loading the reads from the BAM file before counting. Defaults to the params stored in the FraserDataSet object. coldata The colData as given by the FraserDataSet object. A list of GRanges objects containing the counts that should be merged into one countList object. assumeEqual Logical indicating whether all objects in countList can be assumed to contain counts for the same ranges. If FALSE, merging of the ranges is performed. spliceSiteCoords A GRanges object containing the positions of the splice sites. If it is NULL, then splice sites coordinates are calculated first based on the positions of the

junctions defined from the split reads.

Details

The functions described in this file extract and count both the split and the non-spliced reads from bam files.

countRNAData is the main function that takes care of all counting steps and returns a FraserDataSet containing the counts for all samples in the fds.

getSplitReadCountsForAllSamples counts split reads for all samples and getNonSplitReadCountsForAllSamples counts non split reads overlapping splice sites for all samples. addCountsToFraserDataSet adds these counts to an existing fds.

countSplitReads calculates the split read counts for a single sample. countNonSplicedReads counts the non split reads overlapping with splice sites for a single sample.

mergeCounts merges the counts from different samples into a single count object, where the counts for junctions that are not present in a sample are set to zero.

Value

countRNAData returns a FraserDataSet.

getSplitReadCountsForAllSamples returns a GRanges object.

```
getNonSplitReadCountsForAllSamples returns a GRanges object.
addCountsToFraserDataSet returns a FraserDataSet.
countSplitReads returns a GRanges object.
mergeCounts returns a SummarizedExperiment object.
countNonSplicedReads returns a GRanges object.
```

Functions

- countRNAData(): This method extracts and counts the split reads and non spliced reads from RNA bam files.
- getSplitReadCountsForAllSamples(): This method creates a GRanges object containing the split read counts from all specified samples.
- getNonSplitReadCountsForAllSamples(): This method creates a GRanges object containing the non split read counts at the exon-intron boundaries inferred from the GRanges object containing the positions of all the introns in this dataset.
- addCountsToFraserDataSet(): This method adds the split read and non split read counts to a existing FraserDataSet containing the settings.
- countSplitReads(): This method counts all split reads in a bam file for a single sample.
- mergeCounts(): This method merges counts for multiple samples into one SummarizedExperiment object.
- countNonSplicedReads(): This method counts non spliced reads based on the given target (acceptor/donor) regions for a single sample.

Examples

```
# On Windows SNOW is the default for the parallele backend, which can be
# very slow for many but small tasks. Therefore, we will use
# for the example the SerialParam() backend.
if(.Platform$OS.type != "unix") {
    register(SerialParam())
}
fds <- countRNAData(createTestFraserSettings())</pre>
```

createTestFraserSettings

Create a test dataset

Description

Create a test case dataset based on the test sample annotation to be used in the vignette and to explore the functionality of the FRASER package. Dependent on the request only the sample annotation or a full fitted model is returned.

12 estimateBestQ

Usage

```
createTestFraserSettings(workingDir = "FRASER_output")
createTestFraserDataSet(
  workingDir = "FRASER_output",
  rerun = FALSE,
  metrics = "jaccard"
)
```

Arguments

workingDir Directory where to store HDF5 and RDS files. Defaults to FRASER_output in

the current working directory.

rerun Defaults to FALSE. If set to TRUE it reruns the full fit of the model.

metrics The splice metrics that should be included in the test fds. One or several of

'jaccard', 'psi5', 'psi3' or 'theta'.

Value

A FraserDataSet object that contains a test case

Examples

```
fds <- createTestFraserSettings()
fds

fds <- createTestFraserDataSet()
fds</pre>
```

estimateBestQ

Find optimal encoding dimension

Description

Finds the optimal encoding dimension by either Optimal Hard Thresholding or injecting artificial splicing outlier ratios while maximizing the precision-recall curve.

```
estimateBestQ(
  fds,
  type = "jaccard",
  useOHT = TRUE,
  implementation = "PCA",
  q_param = getEncDimRange(fds),
  noise_param = 0,
```

estimateBestQ 13

```
minDeltaPsi = 0.1,
iterations = 5,
setSubset = 50000,
injectFreq = 0.01,
BPPARAM = bpparam(),
internalThreads = 1,
plot = TRUE,
delayed = ifelse(ncol(fds) <= 300, FALSE, TRUE),
...
)</pre>
```

Arguments

fds A FraserDataSet object

type The type of PSI (jaccard, psi5, psi3 or theta for theta/splicing efficiency)

useOHT If TRUE (default), Optimal Hard Thresholding (OHT) is used to estimate the op-

timal encoding dimension. OHT is only supported for the Intron Jaccard Index.

implementation The method that should be used to correct for confounders.

q_param Vector specifying which values of q should be tested

Noise_param Vector specifying which noise levels should be tested.

minDeltaPsi Minimal delta psi of an intron to be be considered a variable intron.

iterations The maximal number of iterations. When the autoencoder has not yet converged

after these number of iterations, the fit stops anyway.

setSubset The size of the subset of the most variable introns that should be used for the

hyperparameter optimization.

injectFreq The frequency with which outliers are injected into the data.

BPPARAM the BiocParallel parameters for the parallelization

internalThreads

The number of threads used internally.

plot If TRUE, a plot of the singular values (OHT) or the area under the curve and

the model loss for each evaluated parameter combination (grid-search) will be

displayed after the estimation procedure.

delayed If FALSE, count matrices will be loaded into memory (faster calculations), oth-

erwise the function works on the delayedMatrix representations (more memory efficient). The default value depends on the number of samples in the fds-object.

... Additional parameters passed to injectOutliers.

Value

FraserDataSet

See Also

FRASER

14 filterVariability

Examples

```
# generate data
fds <- makeSimulatedFraserDataSet(m=15, j=20)
fds <- calculatePSIValues(fds)

# run OHT
fds <- estimateBestQ(fds, type="jaccard", useOHT=TRUE)

# run hyperparameter optimization
fds <- estimateBestQ(fds, type="jaccard", useOHT=FALSE, q_param=c(2, 5))

# get estimated optimal dimension of the latent space
bestQ(fds, type="jaccard")
hyperParams(fds, type="jaccard")</pre>
```

filterVariability

Filtering FraserDataSets

Description

This method can be used to filter out introns that are not reliably detected and to remove introns with no variablity between samples.

```
filterVariability(object, ...)
filterExpressionAndVariability(
  object,
 minExpressionInOneSample = 20,
 quantile = 0.75,
  quantileMinExpression = 10,
 minDeltaPsi = 0,
  filter = TRUE,
 delayed = ifelse(ncol(object) <= 300, FALSE, TRUE),</pre>
  filterOnJaccard = TRUE,
 BPPARAM = bpparam()
)
## S4 method for signature 'FraserDataSet'
filterExpression(
  object,
 minExpressionInOneSample = 20,
 quantile = 0.75,
  quantileMinExpression = 10,
  filter = TRUE,
  delayed = ifelse(ncol(object) <= 300, FALSE, TRUE),</pre>
```

filterVariability 15

```
filterOnJaccard = TRUE,
   BPPARAM = bpparam()
)

## S4 method for signature 'FraserDataSet'
filterVariability(
   object,
   minDeltaPsi = 0,
   filter = TRUE,
   delayed = ifelse(ncol(object) <= 300, FALSE, TRUE),
   filterOnJaccard = TRUE,
   BPPARAM = bpparam()
)</pre>
```

Arguments

object A FraserDataSet object

Further parameters passed on to Rsubread::featureCounts.

minExpressionInOneSample

The minimal read count in at least one sample that is required for an intron to

pass the filter.

quantile Defines which quantile should be considered for the filter.

quantileMinExpression

The minimum read count an intron needs to have at the specified quantile to pass

the filter.

minDeltaPsi Only introns for which the maximal difference in the psi value of a sample to

the mean psi of the intron is larger than this value pass the filter.

filter If TRUE, a subsetted fds containing only the introns that passed all filters is

returned. If FALSE, no subsetting is done and the information of whether an

intron passed the filters is only stored in the mcols.

delayed If FALSE, count matrices will be loaded into memory, otherwise the function

works on the delayedMatrix representations. The default value depends on the

number of samples in the fds-object.

filterOnJaccard

If TRUE, the Intron Jaccard Metric is used to define express introns during fit-

lering. Otherwise, the psi5, psi3 and theta metrics are used (default: TRUE).

BPPARAM the BiocParallel parameters for the parallelization

Value

A FraserDataSet with information about which junctions passed the filters. If filter=TRUE, the filtered FraserDataSet is returned.

Functions

• filterExpressionAndVariability(): This functions filters out both introns with low read support and introns that are not variable across samples.

16 fit

• filterExpression(FraserDataSet): This function filters out introns and corresponding splice sites that have low read support in all samples.

• filterVariability(FraserDataSet): This function filters out introns and corresponding splice sites that have low read support in all samples.

Examples

fit

Fitting the denoising autoencoder

Description

This method corrects for confounders in the data and fits a beta-binomial distribution to the introns/splice sites.

For more details please see FRASER.

```
## S3 method for class 'FraserDataSet'
fit(
 object,
 implementation = c("PCA", "PCA-BB-Decoder", "AE", "AE-weighted", "PCA-BB-full",
  "fullAE", "PCA-regression", "PCA-reg-full", "PCA-BB-Decoder-no-weights", "BB"),
  type = psiTypes,
  rhoRange = c(-30, 30),
 weighted = FALSE,
  noiseAlpha = 1,
  convergence = 1e-05,
  iterations = 15,
  initialize = TRUE,
  control = list(),
 BPPARAM = bpparam(),
 nSubset = 15000,
 minDeltaPsi = 0.1,
)
```

Arguments

object A FraserDataSet object

implementation The method that should be used to correct for confounders.

The encoding dimensions to be used during the fitting proceadure. Should be

fitted using estimateBestQ if unknown. If a named vector is provided it is used

for the different splicing types.

type The type of PSI (jaccard, psi5, psi3 or theta for theta/splicing efficiency)

rhoRange Defines the range of values that rho parameter from the beta-binomial distribu-

tion is allowed to take. For very small values of rho, the loss can be instable, so

it is not recommended to allow rho < 1e-8.

weighted If TRUE, the weighted implementation of the autoencoder is used

noiseAlpha Controls the amount of noise that is added for the denoising autoencoder.

convergence The fit is considered to have converged if the difference between the previous

and the current loss is smaller than this threshold.

iterations The maximal number of iterations. When the autoencoder has not yet converged

after these number of iterations, the fit stops anyway.

initialize If FALSE and a fit has been previoulsy run, the values from the previous fit will

be used as initial values. If TRUE, (re-)initialization will be done.

control List of control parameters passed on to optim().

BPPARAM the BiocParallel parameters for the parallelization

nSubset The size of the subset to be used in fitting if subsetting is used.

minDeltaPsi Minimal delta psi of an intron to be be considered a variable intron.

... Currently not used

Value

FraserDataSet

See Also

FRASER

FRASER: Find RAre Splicing Events in RNA-seq data

Description

This help page describes the FRASER function which can be used run the default FRASER pipeline. This pipeline combines the beta-binomial fit, the computation of Z scores and p values as well as the computation of delta-PSI values.

Usage

```
FRASER(
  fds,
  type = fitMetrics(fds),
  implementation = c("PCA", "PCA-BB-Decoder", "AE-weighted", "AE", "BB"),
  iterations = 15,
 BPPARAM = bpparam(),
  correction,
  subsets = NULL,
)
calculateZscore(fds, type = currentType(fds), logit = TRUE)
calculatePvalues(
  fds,
  type = currentType(fds),
  implementation = "PCA",
  BPPARAM = bpparam(),
  distributions = c("betabinomial"),
  capN = 5 * 1e+05
)
calculatePadjValues(
  fds,
  type = currentType(fds),
  method = "BY",
  rhoCutoff = NA,
  geneLevel = TRUE,
  geneColumn = "hgnc_symbol",
  subsets = NULL,
 BPPARAM = bpparam()
)
calculatePadjValuesOnSubset(
  genesToTest,
  subsetName,
  type = currentType(fds),
  method = "BY",
  geneColumn = "hgnc_symbol",
  BPPARAM = bpparam()
)
```

Arguments

fds A FraserDataSet object

q	The encoding dimensions to be used during the fitting proceadure. Should be fitted using <code>estimateBestQ</code> if unknown. If a named vector is provided it is used for the different splicing types.
type	The type of PSI (jaccard, psi5, psi3 or theta for theta/splicing efficiency)
implementation	The method that should be used to correct for confounders.
iterations	The maximal number of iterations. When the autoencoder has not yet converged after these number of iterations, the fit stops anyway.
BPPARAM	A BiocParallel object to run the computation in parallel
correction	Deprecated. The name changed to implementation.
subsets	A named list of named lists specifying any number of gene subsets (can differ per sample). For each subset, FDR correction will be limited to genes in the subset, and the FDR corrected pvalues stored as an assay in the fds object in addition to the transcriptome-wide FDR corrected pvalues. See the examples for how to use this argument.
• • •	Additional parameters passed on to the internal fit function
logit	Indicates if z scores are computed on the logit scale (default) or in the natural (psi) scale.
distributions	The distribution based on which the p-values are calculated. Possible are beta-binomial, binomial and normal.
capN	Counts are capped at this value to speed up the p-value calculation
method	The p.adjust method that should be used for genome-wide multiple testing correction.
rhoCutoff	The cutoff value on the fitted rho value (overdispersion parameter of the betabinomial) above which junctions are masked with NA during p value adjustment (default: NA, no masking).
geneLevel	Logical value indiciating whether gene-level p values should be calculated. Defaults to TRUE.
geneColumn	The column name of the column that has the gene annotation that will be used for gene-level pvalue computation.
genesToTest	A named list with the subset of genes to test per sample. The names must correspond to the sampleIDs in the given fds object.
subsetName	The name under which the resulting FDR corrected pvalues will be stored in metadata(fds).

Details

All computed values are returned as an FraserDataSet object. To have more control over each analysis step, one can call each function separately.

- fit to control for confounding effects and fit the beta binomial model parameters
- calculatePvalues to calculate the nominal p values
- calculatePadjValues to calculate adjusted p values (per sample)
- calculateZscore to calculate the Z scores

Available methods to correct for the confounders are currently: a denoising autoencoder with a BB loss ("AE" and "AE-weighted"), PCA ("PCA"), a hybrid approach where PCA is used to fit the latent space and then the decoder of the autoencoder is fit using the BB loss ("PCA-BB-Decoder"). Although not recommended, it is also possible to directly fit the BB distribution to the raw counts ("BB").

Value

FraserDataSet

Functions

- FRASER(): This function runs the default FRASER pipeline combining the beta-binomial fit, the computation of Z scores and p values as well as the computation of delta-PSI values.
- calculateZscore(): This function calculates z-scores based on the observed and expected logit psi.
- calculatePvalues(): This function calculates two-sided p-values based on the beta-binomial distribution (or binomial or normal if desired). The returned p values are not yet adjusted with Holm's method per donor or acceptor site, respectively.
- calculatePadjValues(): This function adjusts the previously calculated p-values per sample for multiple testing. First, the previously calculated junction-level p values are adjusted with Holm's method per donor or acceptor site, respectively. Then, if gene symbols have been annotated to junctions (and not otherwise requested), gene-level p values are computed.
- calculatePadjValuesOnSubset(): This function does FDR correction only for all junctions in a certain subset of genes which can differ per sample. Requires gene symbols to have been annotated to junctions. As with the full FDR correction across all junctions, first the previously calculated junction-level p values are adjusted with Holm's method per donor or acceptor site, respectively. Then, gene-level p values are computed.

Author(s)

Christian Mertes <mertes@in.tum.de>
Ines Scheller <scheller@in.tum.de>

See Also

fit

```
# set default parallel backend
register(SerialParam())

# preprocessing
fds <- createTestFraserDataSet()

# filtering not expressed introns
fds <- calculatePSIValues(fds)
fds <- filterExpressionAndVariability(fds)</pre>
```

FraserDataSet 21

```
# Run the full analysis pipeline: fits distribution and computes p values
fds <- FRASER(fds, q=2, implementation="PCA")</pre>
# afterwards, the fitted fds-object can be saved and results can
# be extracted and visualized, see ?saveFraserDataSet, ?results and
# ?plotVolcano
# The functions run inside the FRASER function can also be directly
# run themselves.
# To directly run the fit function:
fds <- fit(fds, implementation="PCA", q=2, type="jaccard")</pre>
# To directly run the nomial and adjusted p value and z score
# calculation, the following functions can be used:
fds <- calculatePvalues(fds, type="jaccard")</pre>
head(pVals(fds, type="jaccard"))
fds <- calculatePadjValues(fds, type="jaccard", method="BY")</pre>
head(padjVals(fds, type="jaccard"))
fds <- calculateZscore(fds, type="jaccard")</pre>
head(zScores(fds, type="jaccard"))
# example of restricting FDR correction to subsets of genes of interest
genesOfInterest <- list("sample1"=c("TIMMDC1"), "sample2"=c("MCOLN1"))</pre>
fds <- calculatePadjValues(fds, type="jaccard",</pre>
                  subsets=list("exampleSubset"=genesOfInterest))
padjVals(fds, type="jaccard", subsetName="exampleSubset")
padjVals(fds, type="jaccard", level="gene", subsetName="exampleSubset")
fds <- calculatePadjValues(fds, type="jaccard",</pre>
                  subsets=list("anotherExampleSubset"=c("TIMMDC1")))
padjVals(fds, type="jaccard", subsetName="anotherExampleSubset")
# only adding FDR corrected pvalues on a subset without calculating
# transcriptome-wide FDR again:
fds <- calculatePadjValuesOnSubset(fds, genesToTest=genesOfInterest,</pre>
         subsetName="setOfInterest", type="jaccard")
padjVals(fds, type="jaccard", subsetName="setOfInterest")
```

FraserDataSet

The FRASER dataset object

Description

Constructs an FRASER object based on the given input. It can take only the annotation (colData) or count tables (junctions/spliceSites).

```
FraserDataSet(colData = NULL, junctions = NULL, spliceSites = NULL, ...)
```

22 FraserDataSet-class

Arguments

colData A DataFrame containing the annotation of the samples junctions, spliceSites

A data.frame like object containing the raw counts for each junction or splice site. It requires the columns startID and endID for the junctions and spliceSiteID and type for the splice sites. Those columns identifies the corresponding splice site for the given junction and map to the splice site. For each sample the counts are saved in a corresponding column with the same name. It can also be a GRange object.

Any parameters corresponding to the slots and their possible values. See Fraser-DataSet

Value

A FraserDataSet object.

Author(s)

Christian Mertes <mertes@in.tum.de>

Examples

FraserDataSet-class FraserDataSet

Description

This class is designed to store the whole FRASER data set needed for an analysis of a disease cohort

Author(s)

Christian Mertes <mertes@in.tum.de>

getter_setter_functions 23

```
getter_setter_functions
```

Getter/Setter functions

Description

This is a collection of small accessor/setter functions for easy access to the values within the FRASER model.

```
featureExclusionMask(fds, type = currentType(fds))
featureExclusionMask(fds, type = currentType(fds)) <- value</pre>
rho(fds, type = currentType(fds))
zScores(fds, type = currentType(fds), byGroup = FALSE, ...)
pVals(
  fds,
  type = currentType(fds),
  level = "site",
  filters = list(),
  dist = "BetaBinomial",
padjVals(
  fds,
  type = currentType(fds),
  dist = c("BetaBinomial"),
  level = "site",
  subsetName = NULL,
  filters = list(),
)
availableFDRsubsets(fds)
predictedMeans(fds, type = currentType(fds))
deltaPsiValue(fds, type = currentType(fds))
currentType(fds)
currentType(fds) <- value</pre>
```

```
fitMetrics(fds)
fitMetrics(fds) <- value

pseudocount(value = NULL)
hyperParams(fds, type = currentType(fds), all = FALSE)
bestQ(fds, type = currentType(fds))
dontWriteHDF5(fds)
dontWriteHDF5(fds) <- value
verbose(fds)</pre>
```

Arguments

fds	An FraserDataSet object.
type	The type of psi (psi5, psi3

type The type of psi (psi5, psi3 or theta)
value The new value to be assigned.

byGroup If TRUE, aggregation by donor/acceptor site will be done.

... Internally used parameters.

level Indicates if the retrieved p values should be adjusted on the donor/acceptor site-

level (default) or if unadjusted junction-level p values should be returned.

filters A named list giving the filters that were applied for masking during p value

correction. Used for storing and retrieving the correct set of requested p values.

dist Distribution for which the p-values should be extracted.

subsetName The name of a subset of genes of interest for which FDR corrected pvalues

were previously computed. Default is NULL (using transcriptome-wide FDR

corrected pvalues).

all Logical value indicating whether hyperParams(fds) should return the results

of all evaluated parameter combinations or only for the optimal parameter com-

bination.

Value

A (delayed) matrix or vector dependent on the type of data retrieved.

Functions

• featureExclusionMask(): Retrieves a logical vector indicating for each junction whether it is included or excluded during the fitting procedure.

getter_setter_functions 25

• featureExclusionMask(fds, type = currentType(fds)) <- value: To remove certain junctions from being used in the train step of the encoding dimension we can set the featureExclusion vector to FALSE. This can be helpfull if we have local linkage between features which we do not want to model by the autoencoder.

- rho(): Returns the fitted rho values for the beta-binomial distribution
- zScores(): This returns the calculated z-scores.
- pVals(): This returns the calculated p-values.
- padjVals(): This returns the adjusted p-values.
- availableFDRsubsets(): This returns the names of FDR subsets for which adjusted p values have been calculated.
- predictedMeans(): This returns the fitted mu (i.e. psi) values.
- deltaPsiValue(): Returns the difference between the observed and the fitted psi values.
- currentType(): Returns the psi type that is used within several methods in the FRASER package (defaults to jaccard).
- currentType(fds) <- value: Sets the psi type that is to be used within several methods in the FRASER package.
- fitMetrics(): Returns the splice metrics that will be fitted (defaults to jaccard, used within several methods in the FRASER package).
- fitMetrics(fds) <- value: Sets the splice metrics that will be fitted (used within several methods in the FRASER package).
- pseudocount(): Sets and returns the pseudo count used within the FRASER fitting procedure.
- hyperParams(): This returns the results of OHT or the hyperparameter optimization. It returns NULL if estimateBestQ was not run yet.
- bestQ(): This returns the optimal latent space dimension estimated by OHT or a hyperparameter optimization. If estimateBestQ has not been run before, a simple estimate of about a tenth of the number of samples is returned.
- dontWriteHDF5(): Gets the current value of whether the assays should be stored as hdf5 files.
- dontWriteHDF5(fds) <- value: Sets whether the assays should be stored as hdf5 files.
- verbose(): Dependent on the level of verbosity the algorithm reports more or less to the user. 0 means being quiet and 10 means everything.
- verbose(fds) <- value: Sets the verbosity level to a value between 0 and 10. 0 means being quiet and 10 means reporting everything.

```
fds <- createTestFraserDataSet()

# should assays be saved as hdf5?
dontWriteHDF5(fds)
dontWriteHDF5 <- TRUE

# get/set the splice metric for which results should be retrieved currentType(fds) <- "jaccard"</pre>
```

26 injectOutliers

```
currentType(fds)
# get fitted parameters
bestQ(fds)
predictedMeans(fds)
rho(fds)
# get statistics
pVals(fds)
padjVals(fds)
# zscore not calculated by default
fds <- calculateZscore(fds, type="jaccard")</pre>
zScores(fds)
# set and get pseudocount
pseudocount(4L)
pseudocount()
# retrieve or set a mask to exclude certain junctions in the fitting step
featureExclusionMask(fds, type="jaccard") <- sample(</pre>
        c(FALSE, TRUE), nrow(mcols(fds, type="jaccard")), replace=TRUE)
featureExclusionMask(fds, type="jaccard")
# controlling the verbosity level of the output of some algorithms
verbose(fds) <- 2</pre>
verbose(fds)
```

injectOutliers

Inject artificial outliers in an existing fds

Description

Inject artificial outliers in an existing fds

```
injectOutliers(
  fds,
  type = psiTypes,
  freq = 0.001,
  minDpsi = 0.2,
  minCoverage = 2,
  deltaDistr = "uniformDistr",
  verbose = FALSE,
  method = c("samplePSI", "meanPSI", "simulatedPSI"),
  BPPARAM = bpparam()
)
```

K 27

Arguments

fds FraserDataSet type The psi type

freq The injection frequency.

minDpsi The minimal delta psi with which outliers will be injected.

minCoverage The minimal total coverage (i.e. N) required for a junction to be considered for

injection of an outlier.

deltaDistr The distribution from which the delta psi value of the injections is drawn (de-

fault: uniform distribution).

verbose Should additional information be printed during computation?

method Defines by which method the new psi of injections is computed, i.e. to which

value the delta psi of the injection is added: "meanPSI" for adding to the mean psi of the junction over all samples or "samplePSI" to add to the psi value of the junction in the specific sample. "simulatedPSI" is only possible if a simulated

dataset is used.

BPPARAM A BiocParallel object to run the computation in parallel

Value

FraserDataSet

Examples

```
# A generic dataset
fds <- makeSimulatedFraserDataSet()
fds <- calculatePSIValues(fds)
fds <- injectOutliers(fds, minDpsi=0.2, freq=1E-3)</pre>
```

Κ

Getter/setter for count data

Description

Getter/setter for count data setter for count data

```
K(fds, type = currentType(fds))
N(fds, type = currentType(fds))
## S4 method for signature 'FraserDataSet'
counts(object, type = currentType(object), side = c("ofInterest", "otherSide"))
```

```
## S4 replacement method for signature 'FraserDataSet,ANY'
counts(
  object,
  type = currentType(object),
  side = c("ofInterest", "otherSide"),
  ...
) <- value</pre>
```

Arguments

fds, object	FraserDataSet
type	The psi type.
side	"ofInterest" for junction counts, "other" for sum of counts of all other junctions at the same donor site (psi5) or acceptor site (psi3), respectively.
	Further parameters that are passed to assays(object,)
value	An integer matrix containing the counts.

Value

FraserDataSet

Examples

```
fds <- createTestFraserDataSet()

counts(fds, side="ofInterest")
counts(fds, type="jaccard", side="other")
head(K(fds))
head(K(fds, type="psi5"))
head(K(fds, type="psi3"))
head(N(fds, type="theta"))</pre>
```

```
length,FraserDataSet-method
```

retrieve the length of the object (aka number of junctions)

Description

retrieve the length of the object (aka number of junctions)

```
## S4 method for signature 'FraserDataSet'
length(x)
```

loadFraserDataSet 29

Arguments

x FraserDataSet

Value

Length of the object.

loadFraserDataSet

Loading/Saving FraserDataSets

Description

This is a convenient function to load and save a FraserDataSet object. It looks and saves the FraserDataSet objects and HDF5 files on disk under the given working dir. Internally it uses HDF5 files for all assays.

Usage

```
loadFraserDataSet(dir, name = NULL, file = NULL, upgrade = FALSE)
saveFraserDataSet(fds, dir = NULL, name = NULL, rewrite = FALSE)
```

Arguments

dir A path where to save the objects (replaces the working directory)

name The analysis name of the project (saved within the 'dir')

file The file path to the fds-object.RDS file that should be loaded.

upgrade Should the version of the loaded object be updated?

fds A FraserDataSet object ot be saved

rewrite logical if the object should be rewritten. This makes sense if you have filtered

or subsetted the object and want to save only the subsetted version

Value

FraserDataSet

```
fds <- createTestFraserSettings()
name(fds) <- "saveing_test"

# make sure the object is saved to disc
dontWriteHDF5(fds) <- FALSE
fdsSaved <- saveFraserDataSet(fds)
fdsSaved</pre>
```

```
# load object from disc
fdsLoaded <- loadFraserDataSet(dir=workingDir(fds), name=name(fds))
fdsLoaded
all.equal(fdsSaved, fdsLoaded)</pre>
```

makeSimulatedFraserDataSet

Create an simulated example data set for FRASER

Description

Simulates a data set based on random counts following a beta binomial (or Dirichlet-Multinomial) distribution.

Usage

```
makeSimulatedFraserDataSet(
   m = 100,
   j = 500,
   q = 10,
   distribution = c("BB", "DM"),
   ...
)
```

Arguments

m	Number of simulated samples
j	Number of simulated junctions
q	number of simulated latent variables.
distribution	Either "BB" for a beta-binomial simulation or "DM" for a dirichlet-multinomial simulation.
	Further arguments used to construct the FraserDataSet.

Value

An FraserDataSet containing an example dataset based on simulated data

```
# A generic dataset
fds1 <- makeSimulatedFraserDataSet()
fds1

# A generic dataset with specificed sample size and injection method
fds2 <- makeSimulatedFraserDataSet(m=10, j=100, q=3)
fds2</pre>
```

mergeExternalData 31

Description

To boost its own sequencing data, one can download existing and precounted data. This function merges the existing FraserDataSet with external count data.

Usage

```
mergeExternalData(fds, countFiles, sampleIDs, annotation = NULL)
```

Arguments

fds	A FraserDataSet
countFiles	A character vector of file names pointing to the external count data. The vector has to be names or the files have to start with k_j , k_t
sampleIDs	The samples to be merged from the external data.
annotation	A sample annotation of the external data (optional).

Details

For more details on existing datasets have a look at: https://github.com/gagneurlab/drop#datasets Since FRASER can not hand NA values, the merge will return only the intersecting regions and will drop any non overlapping features. This has to be kept in mind when analysing rare disease samples.

Value

Merged FraserDataSet object.

plotManhattan

Visualization functions for FRASER

Description

The FRASER package provides multiple functions to visualize the data and the results of a full data set analysis.

Plots the p values over the delta psi values, known as volcano plot. Visualizes per sample the outliers. By type and aggregate by gene if requested.

Plot the number of aberrant events per samples

Plots the observed split reads of the junction of interest over all reads coming from the given donor/acceptor.

Plots the observed values of the splice metric across samples for a junction of interest.

Plots the expected psi value over the observed psi value of the given junction.

Plots the quantile-quantile plot

Histogram of the geometric mean per junction based on the filter status

Histogram of minimal delta psi per junction

Count correlation heatmap function

```
plotManhattan(object, ...)
## S4 method for signature 'FraserDataSet'
plotVolcano(
 object,
  sampleID,
  type = fitMetrics(object),
  basePlot = TRUE,
  aggregate = FALSE,
 main = NULL,
 label = NULL,
  deltaPsiCutoff = 0.1,
  padjCutoff = 0.1,
  subsetName = NULL,
)
## S4 method for signature 'FraserDataSet'
plotAberrantPerSample(
  object,
 main,
  type = fitMetrics(object),
  padjCutoff = 0.1,
```

```
deltaPsiCutoff = 0.1,
  aggregate = TRUE,
  subsetName = NULL,
 BPPARAM = bpparam(),
)
plotExpression(
  fds,
  type = fitMetrics(fds),
  idx = NULL,
  result = NULL,
  colGroup = NULL,
  basePlot = TRUE,
 main = NULL,
  label = "aberrant",
  subsetName = NULL,
)
plotSpliceMetricRank(
  type = fitMetrics(fds),
  idx = NULL,
  result = NULL,
  colGroup = NULL,
  basePlot = TRUE,
 main = NULL,
 label = "aberrant",
  subsetName = NULL,
)
plotExpectedVsObservedPsi(
  fds,
  type = fitMetrics(fds),
  idx = NULL,
  result = NULL,
  colGroup = NULL,
 main = NULL,
 basePlot = TRUE,
  label = "aberrant",
  subsetName = NULL,
)
## S4 method for signature 'FraserDataSet'
plotQQ(
```

```
object,
  type = NULL,
  idx = NULL,
  result = NULL,
  aggregate = FALSE,
  global = FALSE,
 main = NULL,
  conf.alpha = 0.05,
  samplingPrecision = 3,
 basePlot = TRUE,
 label = "aberrant",
 Ncpus = min(3, getDTthreads()),
  subsetName = NULL,
)
## S4 method for signature 'FraserDataSet'
plotEncDimSearch(object, type = psiTypes, plotType = c("sv", "auc", "loss"))
plotFilterExpression(
 fds,
 bins = 200,
 legend.position = c(0.8, 0.8),
  onlyVariableIntrons = FALSE
)
plotFilterVariability(
  fds,
 bins = 200,
 legend.position = c(0.8, 0.8),
 onlyExpressedIntrons = FALSE
)
## S4 method for signature 'FraserDataSet'
plotCountCorHeatmap(
 object,
  type = psiTypes,
  logit = FALSE,
  topN = 50000,
  topJ = 5000,
 minMedian = 1,
 minCount = 10,
 main = NULL,
 normalized = FALSE,
  show_rownames = FALSE,
  show_colnames = FALSE,
 minDeltaPsi = 0.1,
  annotation_col = NA,
```

```
annotation_row = NA,
  border_color = NA,
  nClust = 5,
  plotType = c("sampleCorrelation", "junctionSample"),
  sampleClustering = NULL,
 plotMeanPsi = TRUE,
 plotCov = TRUE,
)
plotBamCoverage(
  fds,
 gr,
  sampleID,
 control_samples = sample(samples(fds[, which(samples(fds) != sampleID)]), min(3,
    ncol(fds) - length(sampleID))),
  txdb = NULL,
 min_junction_count = 20,
 highlight_range = NULL,
  highlight_range_color = "firebrick",
  color_annotated = "gray",
  color_novel = "goldenrod3",
  color_sample_interest = "firebrick",
  color_control_samples = "dodgerblue4",
  toscale = c("exon", "gene", "none"),
 mar = c(2, 10, 0.1, 5),
  curvature_splicegraph = 1,
  curvature_coverage = 1,
  splicegraph_labels = c("genomic_range", "id", "name", "none"),
  splicegraph_position = c("top", "bottom"),
)
plotBamCoverageFromResultTable(
  fds,
  result,
  show_full_gene = FALSE,
  txdb = NULL,
 orgDb = NULL,
  res_gene_col = "hgncSymbol",
  res_geneid_type = "SYMBOL",
  txdb_geneid_type = "ENTREZID",
  left_extension = 1000,
 right_extension = 1000,
)
```

```
## S4 method for signature 'FraserDataSet'
plotManhattan(
  object,
  sampleID,
  value = "pvalue",
  type = fitMetrics(object),
  chr = NULL,
  main = paste0("sample: ", sampleID),
  chrColor = c("black", "darkgrey"),
  subsetName = NULL,
  ...
)
```

Arguments

object, fds An FraserDataSet object.

... Additional parameters passed to plot() or plot_ly() if not stated otherwise in the

details for each plot function

sampleID A sample ID which should be plotted. Can also be a vector. Integers are treated

as indices.

type The psi type: either psi5, psi3 or theta (for SE).

basePlot if TRUE (default), use the R base plot version, else use the plotly framework.

aggregate If TRUE, the pvalues are aggregated by gene (default), otherwise junction level

pvalues are used (default for Q-Q plot).

main Title for the plot, if missing a default title will be used.

label Indicates the genes or samples that will be labelled in the plot (only for basePlot=TRUE).

Setting label="aberrant" will label all aberrant genes or samples. Labelling can be turned off by setting label=NULL. The user can also provide a custom

list of gene symbols or sampleIDs.

padjCutoff, deltaPsiCutoff

Significance or delta psi cutoff to mark outliers

subsetName The name of a subset of genes of interest for which FDR corrected pvalues

were previously computed. Those FDR values on the subset will then be used to determine aberrant status. Default is NULL (using transcriptome-wide FDR

corrected pvalues).

BPPARAM BiocParallel parameter to use.

idx A junction site ID or gene ID or one of both, which should be plotted. Can also

be a vector. Integers are treated as indices.

result The result table to be used by the method.

colGroup Group of samples that should be colored.

global Flag to plot a global Q-Q plot, default FALSE

conf.alpha If set, a confidence interval is plotted, defaults to 0.05

samplingPrecision

Plot only non overlapping points in Q-Q plot to reduce number of points to plot.

Defines the digits to round to.

Ncpus Number of cores to use.

plotType The type of plot that should be shown as character string. For plotEncDim-

Search, it has to be either "sv" for a plot of the singular values against their rank, "auc" for a plot of the area under the curve (AUC) or "loss" for the model loss. For the correlation heatmap, it can be either "sampleCorrelation" for a sample-sample correlation heatmap or "junctionSample" for a junction-

sample correlation heatmap.

bins Set the number of bins to be used in the histogram.

legend.position

Set legend position (x and y coordinate), defaults to the top right corner.

onlyVariableIntrons

Logical value indicating whether to show only introns that also pass the vari-

ability filter. Defaults to FALSE.

onlyExpressedIntrons

Logical value indicating whether to show only introns that also pass the expres-

sion filter. Defaults to FALSE.

logit If TRUE, the default, psi values are plotted in logit space.

topN Top x most variable junctions that should be used for the calculation of sample

x sample correlations.

topJ Top x most variable junctions that should be displayed in the junction-sample

correlation heatmap. Only applies if plotType is "junctionSample".

minMedian, minCount, minDeltaPsi

Minimal median $(m \ge 1)$, delta psi $(|\Delta \psi| > 0.1)$, read count $(n \ge 10)$ value of

a junction to be considered for the correlation heatmap.

normalized If TRUE, the normalized psi values are used, the default, otherwise the raw psi

values

show_rownames, show_colnames

Logical value indicating whether to show row or column names on the heatmap

axes.

annotation_col, annotation_row

Row or column annotations that should be plotted on the heatmap.

border_color Sets the border color of the heatmap

nClust Number of clusters to show in the row and column dendrograms.

sampleClustering

A clustering of the samples that should be used as an annotation of the heatmap.

plotMeanPsi, plotCov

If TRUE, then the heatmap is annotated with the mean psi values or the junction

coverage.

gr A GRanges object indicating the genomic range that should be shown in plotBamCoverage.

control_samples

The sampleIDs of the samples used as control in plotBamCoverage.

txdb A TxDb object giving the gene/transcript annotation to use.

min_junction_count

The minimal junction count across samples required for a junction to appear in the splicegraph and coverage tracks of plotBamCoverage.

highlight_range

A GenomicRanges or GenomicRangesList object of ranges to be highlighted in the splicegraph of plotBamCoverage.

highlight_range_color

The color of highlighted ranges in the splicegraph of plotBamCoverage.

color_annotated

The color for exons and junctions present in the given annotation (in the splicegraph of plotBamCoverage).

color_novel The color for novel exons and junctions not present in the given annotation (in the splicegraph of plotBamCoverage).

color_sample_interest

The color in plotBamCoverage for the sample of interest.

color_control_samples

The color in plotBamCoverage for the samples used as controls.

toscale

In plotBamCoverage, indicates which part of the plotted region should be drawn to scale. Possible values are 'exon' (exonic regions are drawn to scale), 'gene' (both exonic and intronic regions are drawn to scale) or 'none' (exonic and intronic regions have constant length) (see SGSeq package).

The margin of the plot area for plotBamCoverage (b,l,t,r). mar

curvature_splicegraph

The curvature of the junction arcs in the splicegraph in plotBamCoverage. Decrease this value for flatter arcs and increase it for steeper arcs.

curvature_coverage

The curvature of the junction arcs in the coverage tracks of plotBamCoverage. Decrease this value for flatter arcs and increase it for steeper arcs.

cex

For controlling the size of text and numbers in plotBamCoverage. splicegraph_labels

> Indicated the format of exon/splice junction labels in the splicegraph of plotBamCoverage. Possible values are 'genomic_range' (gives the start position of the first exon and the end position of the last exon that are shown), 'id' (format E1,... J1,...), 'name' (format type:chromosome:start-end:strand for each feature), 'none' for

no labels (see SGSeq package).

splicegraph_position

The position of the splicegraph relative to the coverage tracks in plotBamCoverage. Possible values are 'top' (default) and 'bottom'.

show_full_gene Should the full genomic range of the gene be shown in plotBamCoverageFromResultTable (default: FALSE)? If FALSE, only a certain region (see parameters left extension and right_extension) around the outlier junction is shown.

orgDb A OrgDb object giving the mapping of gene ids and symbols.

res_gene_col The column name in the given results table that contains the gene annotation.

res_geneid_type

The type of gene annotation in the results table in res_gene_col (e.g. SYM-BOL or ENTREZID etc.). This information is needed for mapping between the results table and the provided annotation in the txdb object.

txdb_geneid_type

The type of gene_id present in genes(txdb) (e.g. ENTREZID). This information is needed for mapping between the results table and the provided annotation in the txdb object.

left_extension Indicating how far the plotted range around the outlier junction should be extended to the left in plotBamCoverageFromResultTable.

right_extension

Indicating how far the plotted range around the outlier junction should be extended to the right in plotBamCoverageFromResultTable.

value Indicates which assay is shown in the manhattan plot. Defaults to 'pvalue'.

Other options are 'deltaPsi' and 'zScore'.

chr Vector of chromosome names to show in plotManhattan. The default is to

show all chromosomes.

chrColor Interchanging colors by chromosome for plotManhattan.

Details

This is the list of all plotting function provided by FRASER:

plotAberrantPerSample()

- plotVolcano()
- plotExpression()
- plotQQ()
- plotExpectedVsObservedPsi()
- plotCountCorHeatmap()
- plotFilterExpression()
- plotFilterVariability()
- plotEncDimSearch()
- plotBamCoverage()
- plotBamCoverageFromResultTable()
- plotManhattan()
- plotSpliceMetricRank()

For a detailed description of each plot function please see the details. Most of the functions share the same parameters.

plotAberrantPerSample: The number of aberrant events per sample are plotted sorted by rank. The ... parameters are passed on to the aberrant function.

plotVolcano: the volcano plot is sample-centric. It plots for a given sample and psi type the negative log10 nominal P-values against the delta psi values for all splice sites or aggregates by gene if requested.

plotExpression: This function plots for a given site the read count at this site (i.e. K) against the total coverage (i.e. N) for the given psi type $(\psi_5, \psi_3, or\theta \text{ (SE)})$ for all samples.

plotQQ: the quantile-quantile plot for a given gene or if global is set to TRUE over the full data set. Here the observed P-values are plotted against the expected ones in the negative log10 space.

plotExpectedVsObservedPsi: A scatter plot of the observed psi against the predicted psi for a given site.

plotSpliceMetricRank: This function plots for a given intron the observed values of the selected splice metrix against the sample rank.

plotCountCorHeatmap: The correlation heatmap of the count data either of the full data set (i.e. sample-sample correlations) or of the top x most variable junctions (i.e. junction-sample correlations). By default the values are log transformed and row centered. The ... arguments are passed to the pheatmap function.

plotFilterExpression: The distribution of FPKM values. If the FraserDataSet object contains the passedFilter column, it will plot both FPKM distributions for the expressed introns and for the filtered introns.

plotFilterVariability: The distribution of maximal delta Psi values. If the FraserDataSet object contains the passedFilter column, it will plot both maximal delta Psi distributions for the variable introns and for the filtered (i.e. non-variable) introns.

plotEncDimSearch: Visualization of the hyperparameter optimization. It plots the encoding dimension against the achieved loss (area under the precision-recall curve). From this plot the optimum should be choosen for the q in fitting process.

plotManhattan: A Manhattan plot showing the junction pvalues by genomic position. Useful to identify if outliers cluster by genomic position.

plotBamCoverage: A sashimi plot showing the read coverage from the underlying bam files for a given genomic range and sampleIDs.

plotBamCoverageFromResultTable: A sashimi plot showing the read coverage from the underlying bam files for a row in the results table. Can either show the full range of the gene with the outlier junction or only a certain region around the outlier.

Value

If base R graphics are used nothing is returned else the plotly or the gplot object is returned.

Examples

```
# create full FRASER object
fds <- makeSimulatedFraserDataSet(m=40, j=200)
fds <- calculatePSIValues(fds)
fds <- filterExpressionAndVariability(fds, filter=FALSE)
# this step should be done for more dimensions in practice
fds <- estimateBestQ(fds, type="jaccard", useOHT=FALSE, q_param=c(2,5,10,25))
# assign gene names to show functionality on test dataset
# use fds <- annotateRanges(fds) on real data
mcols(fds, type="j")$hgnc_symbol <-
    paste0("gene", sample(1:25, nrow(fds), replace=TRUE))
# fit and calculate pvalues
genesOfInterest <- rep(list(paste0("gene", sample(1:25, 10))), 4)
names(genesOfInterest) <- c("sample1", "sample6", "sample15", "sample23")
fds <- FRASER(fds, subsets=list("testSet"=genesOfInterest))</pre>
```

```
# QC plotting
plotFilterExpression(fds)
plotFilterVariability(fds)
plotCountCorHeatmap(fds, "jaccard")
plotCountCorHeatmap(fds, "jaccard", normalized=TRUE)
plotEncDimSearch(fds, type="jaccard")
# extract results
plotAberrantPerSample(fds, aggregate=FALSE)
plotAberrantPerSample(fds, aggregate=TRUE, subsetName="testSet")
plotVolcano(fds, "sample2", "jaccard", label="aberrant")
plotVolcano(fds, "sample1", "jaccard", aggregate=TRUE, subsetName="testSet")
# dive into gene/sample level results
res <- as.data.table(results(fds))</pre>
plotExpression(fds, result=res[1])
plotQQ(fds, result=res[1])
plotExpectedVsObservedPsi(fds, res=res[1])
plotSpliceMetricRank(fds, res=res[1])
# other ways to call these plotting functions
plotExpression(fds, idx=10, sampleID="sample1", type="jaccard")
plotExpression(fds, result=res[1], subsetName="testSet")
plotQQ(fds, idx=10, sampleID="sample1", type="jaccard")
plotQQ(fds, result=res[1], subsetName="testSet")
plotExpectedVsObservedPsi(fds, idx=10, sampleID="sample1", type="jaccard")
plotExpectedVsObservedPsi(fds, result=res[1], subsetName="testSet")
plotSpliceMetricRank(fds, idx=10, sampleID="sample1", type="jaccard")
plotSpliceMetricRank(fds, result=res[1], subsetName="testSet")
# create manhattan plot of pvalues by genomic position
if(require(ggbio)){
    plotManhattan(fds, type="jaccard", sampleID="sample10")
# plot splice graph and coverage from bam files in a given region
if(require(SGSeq)){
    fds <- createTestFraserSettings()</pre>
    gr <- GRanges(seqnames="chr19",</pre>
        IRanges(start=7587496, end=7598895),
        strand="+")
    plotBamCoverage(fds, gr=gr, sampleID="sample3",
        control_samples="sample2", min_junction_count=5,
        curvature_splicegraph=1, curvature_coverage=1,
        mar=c(1, 7, 0.1, 3))
    # plot coverage from bam file for a row in the result table
    fds <- createTestFraserDataSet()</pre>
    require(TxDb.Hsapiens.UCSC.hg19.knownGene)
    txdb <- TxDb.Hsapiens.UCSC.hg19.knownGene
    require(org.Hs.eg.db)
    orgDb <- org.Hs.eg.db
```

potentialImpactAnnotations

Additional result annotations

Description

These functions work on the result table and add additional annotations to the reported introns: the type of potential impact on splicing (e.g. exon skipping, exon truncation, ...), potential occurence of frameshift, overlap with UTR regions as well as a flag for introns that are located in blacklist regions of the genome.

annotateIntronReferenceOverlap adds basic annotations to the fds for each intron based on the overlap of the intron's location with the reference annotation. Has to be run before the result table is created so that the new column can be included in it (see examples).

annotatePotentialImpact annotates each intron in the results table with the type of potential impact on splicing and potential occurence of frameshift (likely, unlikely, inconclusive). Can also calculate overlap with annotated UTR regions. Potential impact can be: annotatedIntron_increasedUsage, annotatedIntron_reducedUsage, exonTruncation, exonElongation, exonTruncation&Elongation, exonSkipping, splicingBeyondGene, multigenicSplicing, downstreamOfNearestGene, upstreamOfNearestGene, complex (everything else). Splice sites (theta metric) annotations indicate how the splice site is located with respect to the reference annotation. The annotated types are: annotated-SpliceSite, exonicRegion, intronicRegion.

flagBlacklistRegions flags introns in the results table on whether or not they are located in a blacklist region of the genome. By default, the blacklist regions as reported in *Amemiya*, *Kundaje* & Boyle (2019) and downloaded from here are used.

Usage

```
annotateIntronReferenceOverlap(fds, txdb, BPPARAM = bpparam())
annotatePotentialImpact(
```

```
result,
  txdb,
  fds,
  addPotentialImpact = TRUE,
  addUTRoverlap = TRUE,
  minoverlap = 5,
  BPPARAM = bpparam()
)

flagBlacklistRegions(
  result,
  blacklist_regions = NULL,
  assemblyVersion = c("hg19", "hg38"),
  minoverlap = 5
)
```

Arguments

fds A FraserDataSet

txdb A txdb object providing the reference annotation.

BPPARAM For controlling parallelization behavior. Defaults to bpparam().

result A result table as generated by FRASER, including the column annotatedJunction

as generated by the function annotateIntronReferenceOverlap.

addPotentialImpact

Logical, indicating if the type of the potential impact should be added to the

results table. Defaults to TRUE.

addUTRoverlap Logical, indicating if the overlap with UTR regions should checked and added

to the results table. Defaults to TRUE.

minoverlap Integer value defining the number of base pairs around the splice site that need

to overlap with UTR or blacklist region, respectivly, to be considered matching.

Defaults to 5 bp.

blacklist_regions

A BED file that contains the blacklist regions. If NULL (default), the BED files that are packaged with FRASER are used (see Details for more information).

assemblyVersion

Indicates the genome assembly version of the intron coordinates. Only used if blacklist_regions is NULL. For other versions, please provide the BED file containing the blacklist regions directly.

Value

An annotated FraserDataSet or results table, respectively

Functions

• annotateIntronReferenceOverlap(): This method calculates basic annotations based on overlap with the reference annotation (start, end, none, both) for the full fds. The overlap type is added as a new column annotatedJunction in mcols(fds).

44 psiTypes

• annotatePotentialImpact(): This method annotates the splice event type to junctions in the given results table.

• flagBlacklistRegions(): This method flags all introns and splice sites in the given results table for which at least one splice site (donor or acceptor) is located in a blacklist region. Blacklist regions of the genome are determined from the provided BED file.

Examples

```
# get data, fit and compute p-values and z-scores
fds <- createTestFraserDataSet()</pre>
# load reference annotation
library(TxDb.Hsapiens.UCSC.hg19.knownGene)
txdb <- TxDb.Hsapiens.UCSC.hg19.knownGene
# add basic annotations for overlap with the reference annotation
# run this function before creating the results table
fds <- annotateIntronReferenceOverlap(fds, txdb)</pre>
# extract results: for this small example dataset, no cutoffs used
# to get some results
res <- results(fds, padjCutoff=NA, deltaPsiCutoff=NA)
# annotate the type of potential impact on splicing and UTR overlap
res <- annotatePotentialImpact(result=res, txdb=txdb, fds=fds)</pre>
# annotate overlap with blacklist regions
res <- flagBlacklistRegions(result=res, assemblyVersion="hg19")</pre>
# show results table containing additional annotations
res
```

psiTypes

Available splice metrics

Description

Available splice metrics

Usage

psiTypes

Format

An object of class character of length 4.

Examples

```
# to show all available splice metrics:
psiTypes
```

results, FraserDataSet-method

Extracting results and aberrant splicing events

Description

The result function extracts the results from the given analysis object based on the given options and cutoffs. The aberrant function extracts aberrant splicing events based on the given cutoffs.

Usage

```
## S4 method for signature 'FraserDataSet'
results(
  object,
  sampleIDs = samples(object),
  padjCutoff = 0.1,
  deltaPsiCutoff = 0.1,
  rhoCutoff = NA,
  aggregate = FALSE,
  collapse = FALSE,
 minCount = 5,
  psiType = psiTypes,
  geneColumn = "hgnc_symbol",
  all = FALSE,
  returnTranscriptomewideResults = TRUE,
  additionalColumns = NULL,
 BPPARAM = bpparam()
)
## S4 method for signature 'FraserDataSet'
aberrant(
  object,
  type = fitMetrics(object),
  padjCutoff = 0.1,
  deltaPsiCutoff = 0.1,
 minCount = 5,
  rhoCutoff = NA,
  by = c("none", "sample", "feature"),
  aggregate = FALSE,
  geneColumn = "hgnc_symbol",
  subsetName = NULL,
  all = FALSE,
```

)

Arguments

object A FraserDataSet object

sampleIDs A vector of sample IDs for which results should be retrieved

padjCutoff The FDR cutoff to be applied or NA if not requested.

deltaPsiCutoff The cutoff on delta psi or NA if not requested.

rhoCutoff The cutoff value on the fitted rho value (overdispersion parameter of the betabi-

nomial) above which junctions are filtered

aggregate If TRUE the returned object is aggregated to the feature level (i.e. gene level).

collapse Only takes effect if aggregate=TRUE. If TRUE, collapses results across the dif-

ferent psi types to return only one row per feature (gene) and sample.

minCount The minimum count value of the total coverage of an intron to be considered as

significant. result

psiType The psi types for which the results should be retrieved.

geneColumn The column name of the column that has the gene annotation that will be used

for gene-level pvalue computation.

all By default FALSE, only significant introns (or genes) are listed in the results.

If TRUE, results are assembled for all samples and introns/genes regardless of

significance.

 ${\tt returnTranscriptomewideResults}$

If FDR corrected pvalues for subsets of genes of interest have been calculated, this parameter indicates whether additionally the transcriptome-wide results should be returned as well (default), or whether only results for those subsets should be

retrieved.

additionalColumns

Character vector containing the names of additional columns from mcols(fds) that should appear in the result table (e.g. ensembl_gene_id). Default is NULL,

so no additional columns are included.

BPPARAM The BiocParallel parameter.

type Splicing type (psi5, psi3 or theta)

by By default none which means no grouping. But if sample or feature is speci-

fied the sum by sample or feature is returned

subsetName The name of a subset of genes of interest for which FDR corrected pvalues

were previously computed. Those FDR values on the subset will then be used to determine aberrant status. Default is NULL (using transcriptome-wide FDR

corrected pvalues).

... Further arguments can be passed to the method. If "n", "padjVals", "dPsi" or

"rhoVals" are given, the values of those arguments are used to define the aberrant

events.

Value

For results: GRanges object containing significant results. For aberrant: Either a of logical values of size introns/genes x samples if "by" is NA or a vector with the number of aberrant events per sample or feature depending on the vaule of "by"

Examples

```
# get data, fit and compute p-values and z-scores
fds <- createTestFraserDataSet()</pre>
# extract results: for this example dataset, no cutoffs are used to
# show the output of the results function
res <- results(fds, all=TRUE)</pre>
res
# aggregate the results by genes (gene symbols need to be annotated first
# using annotateRanges() function)
results(fds, padjCutoff=NA, deltaPsiCutoff=0.1, aggregate=TRUE)
# aggregate the results by genes and collapse over all psi types to obtain
# only one row per gene in the results table
results(fds, padjCutoff=NA, deltaPsiCutoff=0.1, aggregate=TRUE,
        collapse=TRUE)
# get aberrant events per sample: on the example data, nothing is aberrant
# based on the adjusted p-value
aberrant(fds, type="jaccard", by="sample")
# get aberrant events per gene (first annotate gene symbols)
fds <- annotateRangesWithTxDb(fds)</pre>
aberrant(fds, type="jaccard", by="feature", padjCutoff=NA, aggregate=TRUE)
# find aberrant junctions/splice sites
aberrant(fds, type="jaccard")
# retrieve results limiting FDR correction to only a subset of genes
# first, we need to create a list of genes per sample that will be tested
geneList <- list('sample1'=c("TIMMDC1"), 'sample2'=c("MCOLN1"))</pre>
fds <- calculatePadjValues(fds, type="jaccard",</pre>
                 subsets=list("exampleSubset"=geneList))
results(fds, all=TRUE, returnTranscriptomewideResults=FALSE)
```

samples

Getter/Setter methods for the FraserDataSet

Description

The following methods are getter and setter methods to extract or set certain values of a Fraser-DataSet object.

samples sets or gets the sample IDs; condition; nonSplicedReads return a RangedSummarizedExperiment object containing the counts for the non spliced reads overlapping splice sites in the fds.

Usage

```
samples(object)
samples(object) <- value</pre>
condition(object)
condition(object) <- value</pre>
bamFile(object)
bamFile(object) <- value</pre>
name(object)
name(object) <- value</pre>
strandSpecific(object)
strandSpecific(object) <- value</pre>
pairedEnd(object)
pairedEnd(object) <- value</pre>
workingDir(object)
workingDir(object) <- value</pre>
scanBamParam(object)
scanBamParam(object) <- value</pre>
nonSplicedReads(object)
nonSplicedReads(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
samples(object)
## S4 replacement method for signature 'FraserDataSet'
samples(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
```

```
condition(object)
## S4 replacement method for signature 'FraserDataSet'
condition(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
bamFile(object)
## S4 replacement method for signature 'FraserDataSet'
bamFile(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
name(object)
## S4 replacement method for signature 'FraserDataSet'
name(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
workingDir(object)
## S4 replacement method for signature 'FraserDataSet'
workingDir(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
strandSpecific(object)
## S4 replacement method for signature 'FraserDataSet'
strandSpecific(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
pairedEnd(object)
## S4 replacement method for signature 'FraserDataSet'
pairedEnd(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
scanBamParam(object)
## S4 replacement method for signature 'FraserDataSet'
scanBamParam(object) <- value</pre>
## S4 method for signature 'FraserDataSet'
nonSplicedReads(object)
## S4 replacement method for signature 'FraserDataSet'
nonSplicedReads(object) <- value</pre>
FRASER.mcols.get(x, type = NULL, ...)
```

```
FRASER.rowRanges.get(x, type = NULL, ...)
mapSeqlevels(fds, style = "UCSC", ...)
```

Arguments

object A FraserDataSet object.

value The new value that should replace the current one.

x A FraserDataSet object.

type The psi type (psi3, psi5 or theta)

... Further parameters. For mapSeqLevels: further parameters passed to Genome-InfoDb::mapSeqlevels().

fds FraserDataSet

Value

style

Getter method return the respective current value.

Author(s)

```
Christian Mertes <mertes@in.tum.de>
Ines Scheller <scheller@in.tum.de>
```

Examples

```
fds <- createTestFraserDataSet()</pre>
samples(fds)
samples(fds) <- 1:dim(fds)[2]</pre>
condition(fds)
condition(fds) <- 1:dim(fds)[2]</pre>
bamFile(fds) # file.paths or objects of class BamFile
bamFile(fds) <- file.path("bamfiles", samples(fds), "rna-seq.bam")</pre>
name(fds)
name(fds) <- "My Analysis"</pre>
workingDir(fds)
workingDir(fds) <- tempdir()</pre>
strandSpecific(fds)
strandSpecific(fds) <- TRUE</pre>
strandSpecific(fds) <- "reverse"</pre>
strandSpecific(fds)
scanBamParam(fds)
scanBamParam(fds) <- ScanBamParam(mapqFilter=30)</pre>
nonSplicedReads(fds)
rowRanges(fds)
rowRanges(fds, type="theta")
mcols(fds, type="psi5")
mcols(fds, type="theta")
```

The style of the chromosome names.

subset.FRASER 51

```
seqlevels(fds)
seqlevels(mapSeqlevels(fds, style="UCSC"))
seqlevels(mapSeqlevels(fds, style="Ensembl"))
seqlevels(mapSeqlevels(fds, style="dbSNP"))
```

subset.FRASER

Subsetting by indices for junctions

Description

Providing subsetting by indices through the single-bracket operator

Usage

```
## S3 method for class 'FRASER'
subset(x, i, j, by = c("j", "ss"), ..., drop = FALSE)
## S4 method for signature 'FraserDataSet, ANY, ANY, ANY'
x[i, j, by = c("j", "ss"), ..., drop = FALSE]
```

Arguments

x	A FraserDataSet object
i	A integer vector to subset the rows/ranges
j	A integer vector to subset the columns/samples
by	a character (j or ss) defining if we subset by junctions or splice sites
	Parameters currently not used or passed on
drop	No dimension reduction is done. And the drop parameter is currently not used at all.

Value

A subsetted FraserDataSet object

Examples

```
fds <- createTestFraserDataSet()
fds[1:10,2:3]
fds[,samples(fds) %in% c("sample1", "sample2")]
fds[1:10,by="ss"]</pre>
```

Index

* datasets	<pre>bamFile,FraserDataSet-method(samples),</pre>
psiTypes, 44	47
[,FraserDataSet,ANY,ANY,ANY-method	bamFile<-(samples), 47
(subset.FRASER), 51	bamFile<-,FraserDataSet-method
<pre>'featureExclusionMask<-'</pre>	(samples), 47
(getter_setter_functions), 23	bestQ(getter_setter_functions), 23
aberrant, 39	calculatePadjValues (FRASER), 17
aberrant,FraserDataSet-method	calculatePadjValuesOnSubset (FRASER), 17
(results, Fraser Data Set-method),	calculatePSIValues, 6
45	calculatePvalues (FRASER), 17
addCountsToFraserDataSet, 10, 11	calculateZscore (FRASER), 17
addCountsToFraserDataSet (countRNA), 7	condition (samples), 47
annotateIntronReferenceOverlap, 42	condition,FraserDataSet-method
annotateIntronReferenceOverlap	(samples), 47
<pre>(potentialImpactAnnotations),</pre>	condition<- (samples), 47
42	condition<-,FraserDataSet-method
annotatePotentialImpact, 42	(samples), 47
annotatePotentialImpact	countNonSplicedReads, 10, 11
(potentialImpactAnnotations),	countNonSplicedReads (countRNA), 7
42	countRNA, 7
annotateRanges, 3	countRNAData, 10
annotateRangesWithTxDb	countRNAData (countRNA), 7
(annotateRanges), 3	counts, FraserDataSet-method (K), 27
assayNames, FraserDataSet-method, 5	<pre>counts<-,FraserDataSet,ANY-method(K),</pre>
assays, FraserDataSet-method, 5	27
assays<-,FraserDataSet,DelayedMatrix-method	countSplitReads, 10, 11
(assays,FraserDataSet-method),	<pre>countSplitReads (countRNA), 7 createTestFraserDataSet</pre>
J	
<pre>assays<-,FraserDataSet,list-method</pre>	(createTestFraserSettings), 11 createTestFraserSettings, 11
(assays, i raser bacaset-illection),	<pre>createrestrrasersettings, fr currentType (getter_setter_functions),</pre>
assays<-,FraserDataSet,SimpleList-method	23
(assays, FraserDataSet-method),	currentType<-
5	(getter_setter_functions), 23
available.genomes, 9	(getter_setter_runetions), 25
availableFDRsubsets	deltaPsiValue
(getter_setter_functions), 23	(getter_setter_functions), 23
(5	dontWriteHDF5
bamFile (samples), 47	(getter_setter_functions), 23

INDEX 53

dontWriteHDF5<-	K, 27
$(\texttt{getter_setter_functions}), 23$	keytypes, 4
estimateBestQ, 12, 17, 19, 25	length,FraserDataSet-method,28 loadFraserDataSet,29
fds-methods (samples), 47	Toddi i dsei bataset, 29
featureExclusionMask	makeSimulatedFraserDataSet, 30
(getter_setter_functions), 23	mapSeqlevels (samples), 47
featureExclusionMask,	mergeCounts, 10, 11
(getter_setter_functions), 23	mergeCounts (countRNA), 7
featureExclusionMask<-	mergeExternalData, 31
(getter_setter_functions), 23	
filterExpression,FraserDataSet-method	N(K), 27
(filterVariability), 14	name (samples), 47
filterExpressionAndVariability	name, FraserDataSet-method(samples), 47
(filterVariability), 14	name<- (samples), 47
filtering (filterVariability), 14	<pre>name<-,FraserDataSet-method(samples),</pre>
filterVariability, 14	47
filterVariability,FraserDataSet-method	nonSplicedReads(samples),47
(filterVariability), 14	nonSplicedReads,FraserDataSet-method
fit, 16, 20	(samples), 47
<pre>fitMetrics (getter_setter_functions), 23</pre>	<pre>nonSplicedReads<- (samples), 47</pre>
<pre>fitMetrics<- (getter_setter_functions),</pre>	nonSplicedReads<-,FraserDataSet-method
23	(samples), 47
flagBlacklistRegions, 42	
flagBlacklistRegions	<pre>padjVals (getter_setter_functions), 23</pre>
<pre>(potentialImpactAnnotations),</pre>	<pre>padjVals, (getter_setter_functions), 23</pre>
42	pairedEnd(samples),47
FRASER, 13, 16, 17, 17	pairedEnd,FraserDataSet-method
FRASER.mcols.get (samples), 47	(samples), 47
FRASER.rowRanges.get (samples), 47	pairedEnd<-(samples),47
FraserDataSet, 6, 9, 10, 13, 15, 17, 18, 21,	pairedEnd<-,FraserDataSet-method
22, 36, 46	(samples), 47
FraserDataSet-class, 22	pheatmap, 40
	<pre>plotAberrantPerSample (plotManhattan),</pre>
getBSgenome, 9	32
<pre>getNonSplitReadCountsForAllSamples, 10,</pre>	$\verb plotAberrantPerSample,FraserDataSet-method \\$
11	(plotManhattan), 32
getNonSplitReadCountsForAllSamples	plotBamCoverage(plotManhattan), 32
(countRNA), 7	plotBamCoverageFromResultTable
<pre>getSplitReadCountsForAllSamples, 10</pre>	(plotManhattan), 32
getSplitReadCountsForAllSamples	plotCountCorHeatmap(plotManhattan), 32
(countRNA), 7	$\verb plotCountCorHeatmap,FraserDataSet-method \\$
<pre>getter_setter_functions, 23</pre>	(plotManhattan), 32
	plotEncDimSearch (plotManhattan), 32
<pre>hyperParams (getter_setter_functions),</pre>	plotEncDimSearch,FraserDataSet-method
23	(plotManhattan), 32
	plotExpectedVsObservedPsi
injectOutliers, 26	(plotManhattan), 32

54 INDEX

plotExpression (plotManhattan), 32	summarizeJunctions, 9
plotFilterExpression (plotManhattan), 32	
plotFilterVariability(plotManhattan),	useEnsembl, 4
32	
plotFunctions (plotManhattan), 32	verbose (getter_setter_functions), 23
plotManhattan, 32	verbose<- (getter_setter_functions), 23
plotManhattan,FraserDataSet-method	1: 5: (1) 47
(plotManhattan), 32	workingDir (samples), 47
plotQQ (plotManhattan), 32	workingDir,FraserDataSet-method
plotQQ,FraserDataSet-method	(samples), 47
(plotManhattan), 32	workingDir<- (samples), 47
plotSpliceMetricRank (plotManhattan), 32	workingDir<-,FraserDataSet-method
plotVolcano (plotManhattan), 32	(samples), 47
plotVolcano,FraserDataSet-method	
(plotManhattan), 32	zScore, (getter_setter_functions), 23
potentialImpactAnnotations, 42	zScores (getter_setter_functions), 23
predictedMeans	
(getter_setter_functions), 23	
<pre>pseudocount (getter_setter_functions),</pre>	
	
psiTypes, 44	
pVals (getter_setter_functions), 23	
pVals, (getter_setter_functions), 23	
nearly Francisco mathed 45	
results, FraserDataSet-method, 45	
rho (getter_setter_functions), 23	
rho, (getter_setter_functions), 23	
samples, 47	
samples, FraserDataSet-method (samples),	
47	
**	
samples<- (samples), 47	
<pre>samples<-,FraserDataSet-method</pre>	
(samples), 47	
saveFraserDataSet (loadFraserDataSet),	
29	
scanBamParam (samples), 47	
scanBamParam, FraserDataSet-method	
(samples), 47	
scanBamParam<- (samples), 47	
scanBamParam<-,FraserDataSet-method	
(samples), 47	
strandSpecific (samples), 47	
strandSpecific,FraserDataSet-method	
(samples), 47	
strandSpecific<- (samples), 47	
strandSpecific<-,FraserDataSet-method	
(samples), 47	
subset.FRASER, 51	