Applying the Pathway Activity Profiling - PAPi

Raphael Aggio

May 3, 2016

Introduction

This document describes how to use the Pathway Activity Profiling - *PAPi*. PAPi is an R package for predicting the activity of metabolic pathways based solely on metabolomics data. It relates metabolites' abundances with the activity of metabolic pathways. See Aggio, R.B.M; Ruggiero, K. and Villas-Boas, S.G. (2010) - Pathway Activity Profiling (PAPi): from metabolite profile to metabolic pathway activity. Bioinformatics.

1 Step 1 - Building a local KEGG database: build-Database

PAPi makes use of the information available at the Kyoto Encyclopedia of Genes and Genomes - KEGG database (http://www.genome.jp/kegg/). PAPi can be applied in two ways: online and offline. When applying PAPi online, an internet connection is required to collect online biochemical data from KEGG database. On the other hand, when applied offline, PAPi makes use of a local database and no internet connection is required. As KEGG database is constantly updated, PAPi may produce different results when applied online on different days. Offline, PAPi produces always the same results when using the same local database. PAPi is significantly faster when applied offline. As default, PAPi brings a local database and the function buildDatabase can be used to build new local databases. buildDatabase uses the internet connection to create and install new databases inside of the folder R.home("library/PAPi/data/databases/"). If save = TRUE, the new database is also installed in a folder defined by the user. The local database consists on two CSV files containing the required data about metabolites and metabolic pathways. buildDatabase may take up to ten hours depending on the speed of the internet connection. We highly recommend to build a local database once and do all the required analysis using the same database. Doing so, the results obtained from PAPi can be reproduced at anytime and the time required for each analysis is enormously faster.

2 Step 2 - Adding KEGG codes to metabolomics data: addKeggCodes

We consider a typical metabolomics data set: a data frame containing the name of identified metabolites in the first column and their abundances in the different samples in the following columns. Bellow you have an example of a metabolomics data set.

> library(PAPi)

> library(svDialogs)

> data(metabolomicsData)

> print(metabolomicsData)

	Names	Sample1	Sample2	Sample3	Sample4	Sample5	Sample6
1	Replicates	cond1	cond1	cond1	cond2	cond2	cond2
2	Glucose	0.8	0.6	0.7	1.2	1.1	1.5
3	Fructose 6-phosphate	0.3	0.2	0.4	0.6	0.7	0.7
4	Glyceraldehyde 3-phosphate	1.1	1.5	1.2	<na></na>	<na></na>	<na></na>
5	Glycerone phosphate	1.2	1.5	1.3	0.2	0.3	0.2
6	3-Phospho-D-glycerate	0.2	0.3	0.5	1.1	1	0.9

For applying *PAPi*, the name of compounds in the metabolomics data set must be substituted by their respective KEGG codes. This process can be performed manually, however, it may be considerable time consuming according to the size of the input data. The function addKeggCodes automates this process using a KEGG codes list, which is a data frame containing all potentially identifiable compound and their respective KEGG codes. See bellow:

```
> data(keggLibrary)
> print(keggLibrary)
kegg Name
1 C00031 Glucose
2 C05345 Fructose 6-phosphate
3 C00118 Glyceraldehyde 3-phosphate
4 C00111 Glycerone phosphate
5 C00197 3-Phospho-D-glycerate
6 absent Citrate
```

The KEGG codes list can be built as a CSV file containing a list of KEGG codes in the first column and their respective compound names in the second column, following the same format as showed above. Alternatively, the data frame data(keggLibrary) can be used. Ideally, the KEGG codes list should contain every compound potentially identifiable by the protocol you use to analyze and identify metabolites. For example, I have in my KEGG codes list the name of every compound present in my GC-MS library. The KEGG code of a compound can be searched at KEGG website: http://www.genome.jp/kegg/, however, if the argument addCodes = TRUE, the function addKeggCodes will automatically recognize if you have missing compounds - compounds in your input data that are not present in your KEGG codes list. Then, it will ask if you would like to add those compounds to the KEGG codes list in use. If you click on "yes", addKeggCodes will automatically search the KEGG database (using the INTERNET connection) for missing compounds. For each missing compound, addKeggCodes will present to you a list with all the potential matches from KEGG database. If the right compound is in the list, you just have to click on the compound and click on "Ok". The KEGG code and the name of this respective compound will be automatically aded to the bottom of your KEGG codes list. On the other hand, if the desired compound is not in the list presented to you, you have two options: go to the end of the list and click on "SKIP"; or go to the end of the list and click on "OTHER", which will open a new dialog box allowing you to manually add the respective KEGG code. If the compounds' names in the input data are similar to the names used by KEGG database, addKeggCodes will find the KEGG code of most compounds. After looking for KEGG codes, addKeggCodes gives you the option to save the new KEGG codes list as CSV file and generates a data frame as the input data, however, now with KEGG codes instead of compounds' names. Here is an example of how to used addKeggCodes and the result it generates:

> print(metabolomicsData)

	Names	Sample1	Sample2	Sample3	Sample4	Sample5	Sample6	
1	Replicates	cond1	cond1	cond1	cond2	cond2	cond2	
2	Glucose	0.8	0.6	0.7	1.2	1.1	1.5	
3	Fructose 6-phosphate	0.3	0.2	0.4	0.6	0.7	0.7	
4	Glyceraldehyde 3-phosphate	1.1	1.5	1.2	<na></na>	<na></na>	<na></na>	
5	Glycerone phosphate	1.2	1.5	1.3	0.2	0.3	0.2	
6	3-Phospho-D-glycerate	0.2	0.3	0.5	1.1	1	0.9	
>	<pre>print(keggLibrary)</pre>							
	kegg	Name						
1	C00031 (Glucose						
2	C05345 Fructose 6-pho	osphate						
3	C00118 Glyceraldehyde 3-pho	osphate						
4	4 C00111 Glycerone phosphate							
5	C00197 3-Phospho-D-gl	ycerate						
6	absent	Citrate						
>	AddedKegg <- addKeggCodes(
+	metabolomicsData,							
+	keggLibrary,							
+	save = FALSE,							
+	addCodes = TRUE							
+)							
-	[1] "Data frame loaded" [1] "keggCodes - Data frame loaded"							
	 "Every compound in the in "The final data frame was 	nputData	was four				-	FALSE"
>	<pre>print(AddedKegg)</pre>							

	Name	Sample1	Sample2	Sample3	Sample4	Sample5	Sample6
1	Replicates	cond1	cond1	cond1	cond2	cond2	cond2
2	C00197	0.2	0.3	0.5	1.1	1	0.9
3	C05345	0.3	0.2	0.4	0.6	0.7	0.7
4	C00031	0.8	0.6	0.7	1.2	1.1	1.5
5	C00118	1.1	1.5	1.2	<na></na>	<na></na>	<na></na>
6	C00111	1.2	1.5	1.3	0.2	0.3	0.2

3 Step 3 - Applying PAPi: papi

Once the names of metabolites were substituted by their respective KEGG codes, the input data is ready to be analyzed by the function **papi**. Here is an example of how to use **papi** and the results it generates:

```
> data(papiData)
```

```
> print(papiData)
```

	Name	Sample1	Sample2	Sample3	Sample4	Sample5	Sample6	
1	Replicates	cond1	cond1	cond1	cond2	cond2	cond2	
2	C00197	0.2	0.3	0.5	1.1	1	0.9	
3	C05345	0.3	0.2	0.4	0.6	0.7	0.7	
4	C00031	0.8	0.6	0.7	1.2	1.1	1.5	
5	C00118	1.1	1.5	1.2	<na></na>	<na></na>	<na></na>	
6	C00111	1.2	1.5	1.3	0.2	0.3	0.2	

```
> #papiResults <- papi(papiData, save = FALSE, offline = TRUE, localDatabase = "default")
> data(papiResults)
```

```
> head(papiResults)
```

			pathwayname	Sample1	Sample2	Sample3	Sample4
1			Replicates	cond1	cond1	cond1	cond2
22			Bile secretion	140	105	122.5	210
10			Methane metabolism	92	120	100	16
20			ABC transporters	79.2	59.4	69.3	118.8
13	Pentose	and gluc	uronate interconversions	63.25	82.5	68.75	11
12	Nicoti	inate and	nicotinamide metabolism	56.4	70.5	61.1	9.4
	Sample5	Sample6					
1	cond2	cond2					
22	192.5	262.5					
10	24	16					
20	108.9	148.5					
13	16.5	11					
12	14.1	9.4					

When offline = TRUE, PAPi is performed using a local database. When local-Database = "choose", a list of all installed local databases is presented to the user to choose which one to be used.

4 Step 4 - Applying t-test or ANOVA: papiHtest

papiHtest performs a t-test or ANOVA on the input data, which is a data frame produced by the function papi. For identifying to which experimental condition each sample belongs, the first row of the input data may receive the value "Replicates" in the first cell of the first column and a character string indicating to which experimental condition each sample belongs in the following columns, such as described below:

```
> data(papiResults)
```

```
> head(papiResults)
```

			pathwayname	e Sample1	Sample2	Sample3	Sample4
1			Replicates	s cond1	cond1	cond1	cond2
22			Bile secretion	n 140	105	122.5	210
10			Methane metabolism	n 92	120	100	16
20			ABC transporters	5 79.2	59.4	69.3	118.8
13	Pentose	and glucu	ronate interconversions	63.25	82.5	68.75	11
12	Nicoti	nate and	nicotinamide metabolism	n 56.4	70.5	61.1	9.4
	Sample5	Sample6					
1	cond2	cond2					
22	192.5	262.5					
10	24	16					
20	108.9	148.5					
13	16.5	11					
12	14.1	9.4					

Alternatively, papiHtest will ask the user to indicate the samples belonging to each experimental condition. It is all interactive. The argument StatTest of papiHtest is used to determine if a t-test or an ANOVA should be performed. If StatTest = "T-TEST", "T-test", "t-test", "t-TEST", "t" or "T", a t-test is performed. If StatTest = "ANOVA", "Anova", "A" or "a", a ANOVA is performed. If StatTest is missing, t-test will be performed for input data sets showing two experimental conditions and ANOVA will be performed for input data sets showing more than two experimental conditions. papiHtest can be applied as follows:

> head(papiResults)

				pathwayname	Sample1	Sample2	Sample3	Sample4
1				Replicates	cond1	cond1	cond1	cond2
22				Bile secretion	140	105	122.5	210
10			Me	ethane metabolism	92	120	100	16
20				ABC transporters	79.2	59.4	69.3	118.8
13	Pentose	and gluci	ironate	interconversions	63.25	82.5	68.75	11
12	Nicoti	inate and	nicotir	namide metabolism	56.4	70.5	61.1	9.4
:	Sample5	Sample6						
1	cond2	cond2						
22	192.5	262.5						
10	24	16						
20	108.9	148.5						
13	16.5	11						
12	14.1	9.4						

			pat	hwayname	Sample1	Sample2	Sample3	Sample4
1			Re	cond1	cond1	cond1	cond2	
12	Carbon f	fixation	in photosynthetic c	organisms	26.45	34.5	28.75	4.6
7		(Clavulanic acid bios	synthesis	11	15	12	<na></na>
23 Fructose and mannose metabolism 41.6 51.2 46.4 1							19.2	
33	33 Inositol phosphate metabolism 48.3						52.5	8.4
13	13 Methane metabolism 92 120 100						16	
	Sample5	Sample6	pvalues					
1	cond2	cond2	bonferroni					
12	6.9	4.6	0.0391647314284331					
7	<na></na>	<na> <na> O</na></na>						
23	24	21.6	0.0335649923259184					
33	12.6	8.4	0.0391647314284331					
13	24	16	0.0391647314284331					

5 Step 5 - Generating PAPi graph: papiLine

The results produced by PAPi can be automatically plotted in a line graph using the function papiLine. The input data is a data frame such as produced by papi and papiHtest. papiLine was developed to automatically detect graphical parameters, however, some particular data sets may end up in graphs with the legend in the wrong position or axis labels too small. For this reason, papiLine has 20 arguments that allow the user to customize most of the graphical parameters. If save = TRUE, a png file will be saved in a folder defined by the user. If save = FALSE, a new window will pop up with the graph and the user can use the R options for saving the graph to a PDF file. Here is an example of how to use papiLine:

> head(papiResults)

22 Bile secretion 140 105 122.5 210				pathwayname	Sample1	Sample2	Sample3	Sample4
	1			Replicates	cond1	cond1	cond1	cond2
10 Methane metabolism 92 120 100 16	22			Bile secretion	140	105	122.5	210
	10			Methane metabolism	92	120	100	16
20 ABC transporters 79.2 59.4 69.3 118.8	20			ABC transporters	79.2	59.4	69.3	118.8
13 Pentose and glucuronate interconversions 63.25 82.5 68.75 11	13	Pentose	and gluc	uronate interconversions	63.25	82.5	68.75	11
12 Nicotinate and nicotinamide metabolism 56.4 70.5 61.1 9.4	12	Nicot	inate and	nicotinamide metabolism	56.4	70.5	61.1	9.4
Sample5 Sample6		Sample5	Sample6					
1 cond2 cond2	1	cond2	cond2					
22 192.5 262.5	22	192.5	262.5					
10 24 16	10	24	16					
20 108.9 148.5	20	108.9	148.5					

13	16.5	11
12	14.1	9.4
> pa	apiLine	(
+	I	papiResults,
+	1	relative = TRUE,
+	5	<pre>setRef.cond = TRUE,</pre>
+	I	Ref.cond = "cond1",
+	5	save = FALSE
+)		
[1]	"Data 1	frame loaded"